



Upload image

Aug 10, 2020

Preparation of staphylococcal protein-A conjugated to horseradish peroxidase by the periodate method.

Angel A Justiz-Vaillant¹¹University of the West Indies St. Augustine**1** Works for me dx.doi.org/10.17504/protocols.io.bjg8kzw[University of the West Indies](#) angel.vaillant@sta.uwi.eduAngel Justiz-Vaillant
University of the West Indies St. Augustine

ABSTRACT

This reagent can be used in ELISA, Western blotting and Dot blot to detect antigens and antibodies. It is important in the immunodiagnosis of infectious diseases and other problems. I find this useful in the detection of anti-HIV antibodies by ELISA.

DOI

dx.doi.org/10.17504/protocols.io.bjg8kzw

PROTOCOL CITATION

Angel A Justiz-Vaillant 2020. Preparation of staphylococcal protein-A conjugated to horseradish peroxidase by the periodate method.. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bjg8kzw>

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Aug 09, 2020

LAST MODIFIED

Aug 10, 2020

PROTOCOL INTEGER ID

40192

MATERIALS

NAME	CATALOG #	VENDOR
Ammonium Sulfate		P212121
Sodium periodate	SB0875.SIZE.100g	Bio Basic Inc.
sodium borohydride	452882	Sigma Aldrich
Horseradish Peroxidase (HRP) type IV	P8375-25KU	Sigma Aldrich
Staphylococcal Protein-A		Sigma Aldrich

MATERIALS TEXT

Pipettes
20ml to 1000 ml glass
Scale
Incubator
Refrigerator

Freezer
Centrifuges

SAFETY WARNINGS

Pay attention to all details as the times of reactions among the proteins involved in this preparation. It will prevent over-oxidation. The average time of preparation is 18 hours.

BEFORE STARTING

All reagents but specially the enzyme and more importantly the sodium periodate solution has to be prepared freshly before mixing it with the enzyme.

- 1 Horseradish peroxidase (500 µg in 50 µl NaCO₃, pH 9.6) is mixed with freshly made sodium periodate solution (1.71 mg/ml) followed by incubation in the dark for 2 h.
- 2 Mix 500 µg of staphylococcal protein-A (SpA) with an equal amount (500 micrograms) of a mix of horseradish peroxidase-sodium periodate.
- 3 The mixture is incubated for 3 hours at 4°C with gentle agitation.
- 4 Forty µl of freshly prepared NaBH₄ solution (5 mg NaBH₄ /ml 0.1 mM NaOH) is then added to the preparation.
- 5 The preparation is incubated for 90 min at 4°C in the dark with gentle agitation.
- 6 Cold 50% saturated ammonium sulphate solution (pH 7.4) is added drop by drop in the ratio 1:1 (v/v).
- 7 The mixture is then centrifuged for 25 min at 4°C and recover the pellet at the bottom of the tube.
- 8 The pellets is re-suspended in 200 µl of PBS pH=7.4 and dialysed against 1L of PBS for 24 h with 3 buffer changes.
- 9 An equal volume of glycerol is added to the dialysate followed by 100 µl of bovine serum albumin, BSA (20 mg/ ml).
- 10 The conjugate is then stored at -20°C until further used.

11

