



2 ▼

Feb 15, 2022

Immunoprecipitation using Protein A/G Magnetic Beads V.2

New England Biolabs¹¹New England Biolabs

1

dx.doi.org/10.17504/protocols.io.bddai22e**New England Biolabs (NEB)**Tech. support phone: **+1(800)632-7799** email: **info@neb.com****New England Biolabs**
New England Biolabs

This protocol explains immunoprecipitation using Protein A/G Magnetic Beads.

DOI

dx.doi.org/10.17504/protocols.io.bddai22e<https://www.neb.com/protocols/0001/01/01/immunoprecipitation-using-protein-ag-magnetic-beads>

New England Biolabs 2022. Immunoprecipitation using Protein A/G Magnetic Beads. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.bddai22e>

New England Biolabs



Protein A magnetic beads, Protein G magnetic beads, immunoprecipitating, non-specific binding to beads, pre-clearing crude cell extract of proteins, Immunoprecipitation, IP

protocol ,

Mar 07, 2020

Feb 15, 2022

33922

MATERIALS

[☒ Sodium Orthovanadate \(Vanadate\) - 1 ml New England](#)

Biolabs Catalog #P0758S

[☒ Protein A Magnetic Beads - 1 ml New England](#)

Biolabs Catalog #S1425S

[☒ Protein G Magnetic Beads - 1 ml New England](#)

Biolabs Catalog #S1430S

[☒ PMSF Sigma](#)

Aldrich Catalog #P7626

[☒ Bromophenol blue Bio Basic](#)

Inc. Catalog #BB2230.SIZE.25g

[☒ Glycerol Bio Basic](#)

Inc. Catalog #GB0232.SIZE.500ml

[☒ SDS Bio Basic](#)

Inc. Catalog #SB0485.SIZE.100g

[☒ DTT \(Dithiothreitol\) \(> 99% pure\) Protease free Gold](#)

Biotechnology Catalog #DTT

[☒ EGTA Gold](#)

Biotechnology Catalog #E-217

[☒ Triton X-100](#)

Sigma Catalog #93426

[☒ Tris-HCl Life](#)

Technologies Catalog #AM9855

[☒ 2-Mercaptoethanol Sigma](#)

Aldrich Catalog #M3148

[☒ EDTA Fisher](#)

Scientific Catalog #16 004Y

Immunoprecipitation Buffer:

A	B
NaCl	150 mM
Tris-HCl (pH 7.4)	10 mM
EDTA	1 mM
EGTA (pH 8.0)	1 mM
Sodium ortho-vanadate	0.2 mM
PMSF	0.2 mM
Triton X-100	1%
NP-40	0.50%

3X SDS Sample Loading Buffer:

A	B
Tris-HCl (pH 6.8)	187.5 mM
SDS	6%(w/v)
Glycerol	30%
DTT	150 mM
Bromophenol blue	0.03% (w/v)
β-mercaptoethanol)	2%

Please refer to the Safety Data Sheets (SDS) for health and environmental hazards.

Use **25 µl of Protein A/G Magnetic Beads per 200 µl of crude cell lysate** containing 200-500 µg of total protein in a standard immunoprecipitation protocol. It is important to increase the volume of beads proportionately for larger cell lysate volumes.

Prepare *Immunoprecipitation buffer* with the following reagents:

A	B
NaCl	150 mM
Tris-HCl (pH 7.4)	10 mM
EDTA	1 mM
EGTA (pH 8.0)	1 mM
Sodium ortho-vanadate	0.2 mM
PMSF	0.2 mM
Triton X-100	1%
NP-40	0.50%

Prepare *3X SDS Sample Loading Buffer* using the following reagents:

A	B
Tris-HCl (pH 6.8)	187.5 mM
SDS	6%(w/v)
Glycerol	30%
DTT	150 mM
Bromophenol blue	0.03% (w/v)
β-mercaptoethanol)	2%





Cell Lysis

- 1 Rinse a 60 mm culture dish of confluent cells with PBS.

- 2 Lyse the cells with  **0.5 mL cold Immunoprecipitation Buffer** .

Immunoprecipitation buffer is prepared with the following reagents:

A	B
NaCl	150 mM
Tris-HCl (pH 7.4)	10 mM
EDTA	1 mM
EGTA (pH 8.0)	1 mM
Sodium ortho-vanadate	0.2 mM
PMSF	0.2 mM
Triton X-100	1%
NP-40	0.50%

- 3 Maintain constant agitation for  **00:30:00** at  **4 °C** .
- 4 Scrape the cells from the dish.
- 5 Sonicate  **On ice** for  **00:00:05** ; repeat 4 more times:

5.1 Sonicate  **On ice** for  **00:00:05** (1/4).

5.2 Sonicate  **On ice** for  **00:00:05** (2/4).

5.3 Sonicate  **On ice** for  **00:00:05** (3/4).

5.4 Sonicate  On ice for  00:00:05 (4/4).


6



Centrifuge for  00:05:00 at  4 °C .

7



Assay for total protein then adjust concentration to approximately  1 mg/ml with Immunoprecipitation Buffer.

The supernatant is the crude cell lysate.

Immunoprecipitation

8



In a 1.5 ml microcentrifuge tube, add  25 µL protein A/G Magnetic Beads to  200 µL crude cell extract .

Steps 8-12 pre-clear crude cell extract of proteins which can bind non-specifically to the beads.

9

Gently vortex.

10



Incubate at  4 °C for  01:00:00 .

11 Apply magnetic field for ⌚ 00:00:30 to pull beads to the side of the tube.

12 

Pipette supernatant to a clean 1.5 ml microcentrifuge tube and discard the beads.

13 Add 🧴 1 µg - 🧴 5 µg of desired antibody to crude cell lysate.

14 Vortex.

15 

Incubate at 🌡 4 °C for ⌚ 01:00:00 .

If monoclonal antibodies are used, add 🧴 5 µg rabbit anti-mouse IgG antibody .
Vortex and incubate an additional ⌚ 00:30:00 at 🌡 4 °C . Alternatively, Protein G Magnetic Beads ([NEB #S1430S](#)) can be used for immunoprecipitations with monoclonal antibodies.

16 

Add 🧴 25 µL Protein A/G Magnetic Beads suspension .

17 Gently vortex.

18 

Incubate with agitation for ⌚ 01:00:00 at 🌡 4 °C .

19 Apply magnetic field to pull beads to the side of the tube.

20 

Carefully pipette to remove supernatant.

21 

Wash with  **500 µL Immunoprecipitation Buffer** by gentle vortex.

22 

Apply magnetic field, then remove supernatant and discard.

23 

Repeat wash steps two more times:

23.1 

Wash with  **500 µL Immunoprecipitation Buffer** by gentle vortex. (1/2)

23.2 

Apply magnetic field, then remove supernatant and discard. (1/2)

23.3 

Wash with  **500 µL Immunoprecipitation Buffer** by gentle vortex. (2/2)

23.4 

Apply magnetic field, then remove supernatant and discard. (2/2)

24 

Resuspend bead pellet in  **30 µL 3X SDS Sample Loading Buffer** .

3X SDS Sample Loading Buffer is prepared using the following reagents:

A	B
Tris-HCl (pH 6.8)	187.5 mM
SDS	6%(w/v)
Glycerol	30%
DTT	150 mM
Bromophenol blue	0.03% (w/v)
β-mercaptoethanol)	2%

25



Incubate sample at  **70 °C** for  **00:05:00** .

26

Apply magnetic field to sample, then load supernatant on SDS-PAGE gel and electrophorese.