



# Preparation of electrocompetent *Escherichia coli* V.2

lewis.bingle<sup>1</sup>

<sup>1</sup>University of Sunderland

Version 2

Oct 09, 2020

1

Works for me

This protocol is published without a DOI.

SunMedMicro

lewis.bingle

## ABSTRACT

Preparation of competent *Escherichia coli* cells for electroporation.

## PROTOCOL CITATION

lewis.bingle 2020. Preparation of electrocompetent *Escherichia coli*. protocols.io  
<https://protocols.io/view/preparation-of-electrocompetent-escherichia-coli-bm74k9qw>  
Version created by lewis.bingle

## LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Oct 09, 2020

## LAST MODIFIED

Oct 09, 2020

## PROTOCOL INTEGER ID

42972

## GUIDELINES

1. After harvesting and chilling cells (steps 1-2) they should be kept as cold as possible during the rest of the procedure e.g. keep all tubes and buffers on ice whenever possible; do not hold in the hand for long periods; handle tubes by top (lid) rather than bottom.
2. As cells are washed, the pellet sticks increasingly poorly to the wall of the tube after centrifugation – it may be necessary to increase the speed or duration of the centrifugation step.

## MATERIALS TEXT

10 % glycerol in distilled / purified water.

## ABSTRACT

Preparation of competent *Escherichia coli* cells for electroporation.

## BEFORE STARTING

1. Prechill 10 % glycerol solution (or SDW) on ice or in fridge and keep on ice during use.
2. Prechill centrifuge tubes on ice. Prechill centrifuge (if refrigerated) or put rotor / buckets in fridge.

1 Inoculate starter culture of the desired strain from single colony (e.g. 5 mL LB broth in universal bottle). Incubate (37°C with shaking at 200-220 rpm) Overnight .

2 Dilute overnight culture 1 / 100 into 25 mL fresh LB broth (supplemented with appropriate antibiotic if required)

and incubate (37°C with shaking at 200-220 rpm) to mid-log phase growth (OD<sub>600nm</sub> of between 0.5 and 0.7).

- 3 Decant bacteria into 50 ml Falcon (centrifuge) tube and chill on ice (approx. ⌚00:20:00).
- 4 Centrifuge 🌀4000 x g, 4°C for ⌚00:10:00 (use refrigerated centrifuge if possible, or pre-chill centrifuge rotor / bucket). \* Discard supernatant.
- 5 Wash the cells 3x with ice-cold, sterile aqueous glycerol solution [M]10 % (w/v) (can use sterile distilled water) in decreasing volumes of 📏25 mL, 📏12 mL and 📏6 mL (centrifuge as in step 4 between each wash).
- 6 Resuspend final cell pellet in 📏300 µl of sterile 10% (w/v) glycerol (or water) 🧊 On ice.
- 7 Store electrocompetent cells on ice for immediate use or snap-freeze in liquid nitrogen for long-term storage at 🧊 -80 °C.