

May 03, 2024

RAFT Based Synthesis of In-house Polymers

DOI

dx.doi.org/10.17504/protocols.io.ewov19noylr2/v1

Caroline Brown^{1,2,3}, Snehasish Ghosh^{1,2,3}, Kallol Gupta^{1,2,3}

¹Nanobiology Institute, Yale University, West Haven, CT, USA;

²Department of Cell Biology, Yale University School of Medicine, New Haven, CT, USA;

³Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, 20815



Caroline Brown

Yale University

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.ewov19noylr2/v1

Protocol Citation: Caroline Brown, Snehasish Ghosh, Kallol Gupta 2024. RAFT Based Synthesis of In-house Polymers. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.ewov19noylr2/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 03, 2024

Last Modified: May 03, 2024

Protocol Integer ID: 99193

Abstract





The protocol gives details for synthesizing the ChloroSMA series of membrane active polymers.



Chemicals and reagents


- 1 All reagents were purchased from Sigma-Aldrich and used without further purification unless specified otherwise.

Synthesis of ChloroSMA20

- 2 To a  50 mL round bottom flask add  1.52 millimolar (mM) chloro-styrene ,  1.52 millimolar (mM) maleic anhydride , and  0.076 millimolar (mM) 2-(dodecylthiocarbonothioylthio) propionic acid .

- 3 Dissolve in dimethylformamide.

- 4 Add  0.004 millimolar (mM) 1,1-Azobis(cyclohexanecarbonitrile) .

- 5 Cap the flask with a rubber stopper and bubble with nitrogen for  00:15:00 minutes.

15m





- 6 Heat to  90 °C for  16:00:00 hours.

16h

- 7 Precipitate in isopropanol or water, filter, and dry in vacuo.

Synthesis of ChloroSMA 40

16h 15m

- 8 To a  50 mL round bottom flask add  3.1 millimolar (mM) chloro-styrene ,  3.1 millimolar (mM) maleic anhydride , and  0.076 millimolar (mM) 2-(dodecylthiocarbonothioylthio) propionic acid .

- 9 Dissolve in dimethylformamide.



10 Add [M] 0.004 millimolar (mM) 1,1-Azobis(cyclohexanecarbonitrile) .

11 Cap the flask with a rubber stopper and bubble with nitrogen for ⌚ 00:15:00 minutes.

15m

12 Heat to 🌡️ 90 °C for ⌚ 16:00:00 hours.

16h

13 Precipitate in isopropanol or water, filter, and dry in vacuo.

Synthesis of ChloroSMA 60

16h 15m

14 To a 🧪 50 mL round bottom flask add [M] 4.56 millimolar (mM) chloro-styrene ,
[M] 4.56 millimolar (mM) maleic anhydride , and
[M] 0.076 millimolar (mM) 2-(dodecylthiocarbonothioylthio) propionic acid .

15 Dissolve in dimethylformamide.

16 Add [M] 0.004 millimolar (mM) 1,1-Azobis(cyclohexanecarbonitrile) .

17 Cap the flask with a rubber stopper and bubble with nitrogen for ⌚ 00:15:00 minutes.

15m

18 Heat to 🌡️ 90 °C for ⌚ 16:00:00 hours.

16h

19 Precipitate in isopropanol or water, filter, and dry in vacuo.

Synthesis of ChloroSMA 80

16h 15m

20 To a 🧪 50 mL round bottom flask add [M] 6.1 millimolar (mM) chloro-styrene ,
[M] 6.1 millimolar (mM) maleic anhydride , and



[M] 0.076 millimolar (mM) 2-(dodecylthiocarbonothioylthio) propionic acid .

21 Dissolve in dimethylformamide.

22 Add [M] 0.004 millimolar (mM) 1,1-Azobis(cyclohexanecarbonitrile) .

23 Cap the flask with a rubber stopper and bubble with nitrogen for ⌚ 00:15:00 minutes.

15m

24 Heat to 🔥 90 °C for ⌚ 16:00:00 hours.

16h

25 Precipitate in isopropanol or water, filter, and dry in vacuo.

End-group removal of poly(4 chlorostyrene-alt-maleic anhydride)

6h 15m

26 Dissolve precipitated polymer in dimethylformamide and add [M] 9.9 millimolar (mM) benzoyl peroxide to a round bottom flask.

27 Seal with a rubber stopper and bubble with nitrogen for ⌚ 00:15:00 minutes.

15m

28 Leave the escape needle in the flask and heat to 🔥 90 °C for ⌚ 06:00:00 hours.

6h

29 Precipitate polymer with isopropanol or water and dry in vacuo.

Hydrolysis to produce poly(4 chlorostyrene-alt-maleic acid)

1d 16h

30 In a round bottom flask dissolve 🧪 800 mg anhydride polymer in 🧪 2 g KOH and 🧪 20 mL water .

31 Reflux for ⌚ 04:00:00 hours. The mixture should become clear.

4h



32 Dialyze the hydrolyzed polymer using a 3.5KDa mwco membrane for  36:00:00 hours.

1d 12h

33 Add water to precipitate the polymer.

34 Wash polymer using 0.1N HCl, 3 times, and dry the polymer under nitrogen.