



JAN 06, 2023

OPEN ACCESS

Protocol Citation: Andreas Sagen 2023. Phosphate Buffered Saline (PBS). [protocols.io](https://protocols.io/view/phosphate-buffered-saline-pbs-cmdqu25w) <https://protocols.io/view/phosphate-buffered-saline-pbs-cmdqu25w>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Jan 06, 2023

Last Modified: Jan 06, 2023


PROTOCOL integer ID:
74896

Keywords: Buffer, PBS, phosphate buffered saline

🌐 Phosphate Buffered Saline (PBS)

Andreas Sagen¹

¹University of Oslo



Andreas Sagen

University of Oslo, The National Institute of Occupational H...

ABSTRACT

Phosphate-buffered saline (PBS) is a buffer solution (pH ~7.4) commonly used in biological research. It is a water-based salt solution containing disodium hydrogen phosphate, sodium chloride and optionally, potassium chloride and potassium dihydrogen phosphate. The buffer helps to maintain a constant pH. PBS has many uses because it is isotonic and non-toxic to most cells^[1].

GUIDELINES

Follow step by step, unless stated otherwise. Equipment needed should be standard to a microbiology lab.

MATERIALS

Analytical scale, autoclave, bottle, weight vessel, LAF bench

SAFETY WARNINGS



When removing autoclaved components, be sure to take care as this can be very hot.





BEFORE START INSTRUCTIONS

Prepare glassware by cleaning it, and ensure that scale is sufficiently calibrated

- 1 All compounds are measured using a high precision analytical scale from powdered compounds. Each compound is measured to within 1% of the target weight. All compounds are mixed in a Duran bottle

100 mL 10x PBS w/o Ca₂++, Mg₂++

- 1.1 Fill the bottle with 50 mL double-distilled water

1.2 Measure  8000 mg Sodium chloride,  200 mg Potassium chloride,  144 mg Disodium phosphate and  245 mg Monopotassium phosphate


Powdered compounds:

 Sodium chloride **Sigma-aldrich Catalog #S9625**

 Potassium chloride **Sigma-aldrich Catalog #P3911**


 Disodium phosphate **Sigma-aldrich Catalog #S9763**

 Monopotassium phosphate **Sigma-aldrich Catalog #P5379**

1.3 Add powdered solids into bottle, and use a magnetic mixer with a stir bar to mix for  00:05:00

5m

1.4 Adjust pH while mixing to  7.4 using concentrated sodium hydroxide

1.5 Add distilled water to a total of  100 mL

1.6 Autoclave liquid at  121 °C for  00:15:00

15m