



Aug 05, 2020

Protein Extraction

avinash.kale¹¹domnic colvin

1

Works for me

dx.doi.org/10.17504/protocols.io.bjdnki5e

UM-DAE Centre for Excellence in Basic Sciences



avinash.kale

DOI

dx.doi.org/10.17504/protocols.io.bjdnki5e

DOCUMENT CITATION

avinash.kale 2020. Protein Extraction. **protocols.io**<https://dx.doi.org/10.17504/protocols.io.bjdnki5e>

LICENSE

This is an open access document distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Aug 05, 2020

LAST MODIFIED

Aug 05, 2020

DOCUMENT INTEGER ID

40078

Protein Extraction

The four *Bacillus* bacteria were inoculated in 200 ml LB medium (Serva) and incubated at 37°C overnight. The cells were centrifuged at 13000 rpm for 15 minutes at 4°C to pellet the cells. The cell pellet was sonicated at amplitude 30%, pulse 10 seconds, gap 10 seconds, for a total time of 2 minutes: 30 seconds or 15 Cycles to obtain the cell lysate. The TRIzol extraction protocol was used with minor modifications(24). In brief, to the cell lysate, 1 mL of TRIzol reagent and 200 µL of chloroform were added. The solution was mixed vigorously and was centrifuged 14,000g for 15 minutes at 4°C. The Upper clear phase was carefully decanted and to the remaining phase, 300 µL of ethanol was added and incubated at room temperature for 3 minutes followed by centrifugation at 2000g for 15 minutes at 4 °C. To the supernatant, a 4-fold volume of ice-cold acetone was added and kept for overnight incubation at -20°C. The protein pellet was washed thrice first with 0.3 M Guanidium-HCl in 95 % ethanol followed by ice-cold acetone two times. All the protein pellets were air dried and wer subsequently dissolved in rehydration buffer (7 M urea, 2 M thiourea, 2 % CHAPS, Distilled water). Protein concentrations were estimated using Bradfor assay.