


JAN 20, 2023

CUMC-TMC-10X-Gene-Expression

Shaunice Grier¹

¹Columbia University

Cellular Senescence Network (SenNet) Method Development Community



Ivy Kosater

OPEN ACCESS

Protocol Citation: Shaunice Grier 2023. CUMC-TMC-10X-Gene-Expression.

protocols.io

<https://protocols.io/view/cumc-tmc-10x-gene-expression-cm4zu8x6>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Jan 20, 2023

Last Modified: Jan 20, 2023

PROTOCOL integer ID:
75641

ABSTRACT

Post Mortem Human Fresh Frozen Dorsolateral Prefrontal Cortex tissues are cryo-sectioned onto 10X Visium slides.

Tissue Optimization is performed and determines optimal Gene Expression testing conditions.

cDNA samples are processed and generated into libraries. Libraries are qc'd for size and concentration with the Agilent Bioanalyzer 2100, Qubit Fluorometer and qPCR with KAPA Library Quantification Kit. QC'd libraries are submitted for sequencing.

Tissue QC

1

Section

2

H&E Stain

3 Image

4 Assess

Tissue Optimization

5 Permeabilization & cDNA Synthesis

6 Tissue Removal

Gene Expression

7 cDNA Synthesis

8 Second Strand Synthesis & Denaturation

9 cDNA Amplification & QC

Sequencing

- 11** Submitted libraries are evaluated for size distribution on the Fragment Analyzer and quantified using Picogreen and qPCR with the Universal KAPA Library Quantification Kit. Library pools are loaded between 1.2 – 1.5 nM using the standard NovaSeq workflow as outlined in Illumina's NovaSeq 6000 system guide.