



Sep 02, 2020

Modified 1D Native Barcoding genomic DNA protocol from the Temperton Lab (University of Exeter)

Michael Henson¹, Ben Temperton², Cameron Thrash³¹University of Chicago; ²University of Exeter; ³University of Southern California**1** Works for me dx.doi.org/10.17504/protocols.io.bixskfne

Thrash Lab

Michael Henson

ABSTRACT

Modified 1D Native Barcoding genomic DNA protocol (v. NBE_9006_v103_revO_21Dec2016 for SQK-LSK108) from the Temperton Lab (University of Exeter)

DOI

dx.doi.org/10.17504/protocols.io.bixskfne

PROTOCOL CITATION

Michael Henson, Ben Temperton, Cameron Thrash 2020. Modified 1D Native Barcoding genomic DNA protocol from the Temperton Lab (University of Exeter). **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bixskfne>

KEYWORDS

MinION sequencing, 1D Native Barcoding genomic DNA protocol, NBE_9006_v103_revO_21Dec2016, Genome sequencing

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Jul 23, 2020

LAST MODIFIED

Sep 02, 2020

PROTOCOL INTEGER ID

39634

MATERIALS

NAME	CATALOG #	VENDOR
Blunt/TA Ligase Master Mix - 50 rxns	M0367S	New England Biolabs
NEBNext Ultra II End Repair/dA-Tailing Module - 24 rxns	E7546S	New England Biolabs
Ampure XP beads	A63881	Beckman Coulter
Native Barcoding Kit 1D (EXP-NBD103)		
Ligation Sequencing Kit 1D	SQK-LSK108	

NAME	CATALOG #	VENDOR
Qubit® dsDNA HS assay kit, 100 reactions	Q32851	Life Technologies
g-TUBE	520079	Covaris
DNA LoBind Tubes, 1.5 mL	0030108051	Eppendorf
Qubit® dsDNA BR Assay Kit	Q32853	Thermo Fisher

DNA Fragmentation

- 1 Load 46 µL of genomic DNA into a Covaris g-TUBE.
- 2 Spin the g-TUBE for 1 minute in an Eppendorf 5424 centrifuge.

2.1 To achieve ~ XXX bp, spin the tube at 6000 x g.

- 3 Invert the tube and spin again for 1 minute.
- 4 Transfer volume to sterile 1.5 mL Eppendorf DNA LoBind tubes.
*Prior to the next step be sure to analyze your DNA for size, quantity, and quality.

Library Preparation

- 5 Combine: 45 µL genomic DNA, 7 µL Ultra II End-prep reaction buffer, 3 µL Ultra II End-prep enzyme mix, 5 µL of Nuclease-free water.
- 6 Mix tube gently by flicking, then spin down.
- 7 Transfer volume to a new 0.2 mL PCR tube.
- 8 Incubate for 30 minutes at 20°C and 30 minutes at 65°C.
- 9 Add 60 µL of AmPure XP beads and mix by pipetting.
- 10 Incubate for 5 minutes at room temperature (RT), while flicking occasionally.

- 11 Spin down the sample and pellet on the magnet.
- 12 Pipette off supernatant while on the magnet, avoiding pellet.
- 13 Wash beads with 200 μ L of 80% ethanol (EtOH). Do not disturb the pellet. Remove EtOH.
- 14 Repeat step 13
- 15 Spin down the tube and place the tube back in the magnetic rack.
- 16 Remove any residual ethanol.
- 17 Allow beads to briefly dry.
- 18 Remove the tube from the magnetic rack, and resuspend in 25 μ L of prewarmed 55°C nuclease-free water.
- 19 Incubate for 2 minutes at 55°C.
- 20 Pellet beads on magnet until eluate is clear.
- 21 Pipette off 25 μ L of eluate into 1.5 mL Eppendorf LoBind DNA tube.
- 22 Quantify 1 μ L of eluate using Qubit HS dsDNA kit.
- 23 Thaw barcodes at RT.

- 24 Add 22.5 μ L of end-prepped DNA, 2.5 μ L Native barcode, and 25 μ L of Blunt/TA ligase Master Mix.
- 25 Mix gently by flicking tube and spin down.
- 26 Incubate reaction mixture for 10 minutes at RT.
- 27 Add 50 μ L of AMPure XP beads to the reaction mixture and mix by pipetting.
- 28 Incubate for 5 minutes at 55°C while occasionally flicking to mix.
- 29 Spin down the sample and pellet on the magnetic rack.
- 30 Wash beads with 200 μ L of 80% EtOH. Do not disturb the pellet. Remove EtOH.
- 31 Pipette off supernatant, avoiding pellet.
- 32 Repeat step 31
- 33 Spin down the tube and place the tube back in the magnetic rack.
- 34 Remove any residual EtOH.
- 35 Allow beads to briefly dry.
- 36 Remove the tube from the magnetic rack, and resuspend in 25 μ L of prewarmed 55°C nuclease-free water.

- 37 Incubate for 2 minutes at 55°C.
- 38 Pipette off 25 µL of eluate into 1.5 mL Eppendorf LoBind DNA tube.
- 39 Quantify 1 µL of eluate using the Qubit HS dsDNA kit.
- 40 Mix together 22.5 µL of end-prepped DNA, 2.5 µL of Native barcode, and 25 µL of Blunt/TA ligase Master Mix (in that order).
 - 40.1 Mix by gently flicking.
- 41 Spin the tube down.
- 42 Incubate the reaction for 10 minutes at RT.
- 43 Add 50 µL of resuspended AMPure XP beads to the reaction and mix by pipetting.
- 44 Incubate at 55°C for 5 minutes while occasionally flicking to mix.
- 45 Spin down sample(s) and pellet on the magnetic rack.
- 46 Pipette off the supernatant.
- 47 Wash beads with 200 µL of 80% EtOH. Do not disturb the pellet. Remove EtOH.
- 48 Repeat step 47

- 49 Spin down the tube and place the tube back in the magnetic rack.
- 50 Remove any residual EtOH.
- 51 Allow the beads to briefly dry.
- 52 Remove the tube from the magnetic rack, and resuspend in 26 μ L of nuclease-free water.
- 53 Incubate for 2 minutes at RT.
- 54 Pellet beads on the magnet until clear and colorless.
- 55 Transfer 26 μ L into clear, sterile 1.5 mL Eppendorf DNA LoBind tube.
- 56 Quantify 1 μ L of eluate using Qubit HS dsDNA kit.
- 57 Pool equimolar amounts of each barcoded sample into a 1.5 mL Eppendorf DNA LoBind tube (< 2 μ g).
- 58 Quantify 1 μ L of eluate using Qubit HS dsDNA kit.
- 59 Dilute pooled sample to 50 μ L in Nuclease-free water.

Pooled Library Preparation

- 60 Mix together 50 μ L pooled barcoded samples sample, 20 μ L barcode Adapter Mix (BAM), 20 μ L Blunt/TA Master Mix, 10 μ L Quick T4 DNA ligase.

60.1 Mix by flicking between each sequential addition.

- 61 Mix by flicking the tube.
- 62 Spin down the sample.
- 63 Incubate the reaction for 10 minutes at RT.
- 64 Add 40 μ L of resuspended AMPure XP beads to the adapter ligation reaction and mix by pipetting.
- 65 Incubate for 5 minutes at RT, occasionally mixing by flicking.
- 66 Place on a magnetic rack and allow beads to pellet.
- 67 Pipette off supernatant.
- 68 Wash beads with 140 μ L of ABB buffer.
- 69 Resuspend beads in ABB buffer by flicking the tube.
- 70 Pellet beads on the magnetic rack.
- 71 Remove residual ABB buffer.
- 72 Repeat step 68 - 71
- 73 Remove the tube from the magnetic tube and resuspend pellet in 15 μ L of Elution buffer.

- 74 Incubate for 10 minutes at RT.
 - 75 Pellet beads on the magnetic tube rack until the eluate is clear and colorless.
 - 76 Pipet up and retain the 15 μ L of eluate into a clean 1.5 mL Eppendorf DNA LoBind tube
 - 77 Place the tube of eluate (library) on ice until required for library loading
-
- 77.1 Use 1 μ L of DNA library to Quantify using Qubit HS dsDNA assay kit.