

Making E8 medium
For ES/iPS culture

Aug 12, 2021

Making E8 medium for ES/iPS culture

Jiuchun Zhang¹, Harper JW¹¹Harvard Medical School

1 Works for me

Share

dx.doi.org/10.17504/protocols.io.bsacnaaw

Felix Kraus

ABSTRACT

This protocol is about making E8 medium for ES/iPS culture.

ATTACHMENTS

[MakingE8fromDMEM-F12.pdf](#)

DOI

dx.doi.org/10.17504/protocols.io.bsacnaaw

PROTOCOL CITATION

Jiuchun Zhang, Harper JW 2021. Making E8 medium for ES/iPS culture. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bsacnaaw>

KEYWORDS

E8, E8 medium, ES, iPS

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Feb 10, 2021

LAST MODIFIED

Aug 12, 2021

OWNERSHIP HISTORY

Feb 10, 2021 Julia Rossmannith protocols.io

Feb 15, 2021 Felix Kraus

PROTOCOL INTEGER ID

47140

MATERIALS TEXT

Preparation of TGFb1:

- Recombinant Human TGF-beta 1 (R&D 240-B-001MG/CF)
[Albumin human Sigma](#)
- **Aldrich Catalog #A7223** Step 1
- 11.6 M HCl
- molecular grade water
- Millipore filters

Preparation of Transferrin:

- [holo-Transferrin human Sigma](#)
- **Aldrich Catalog #T0665-1G** Step 13
- Millipore Steriflip filter
- PBS

Additional reagents for making E8 medium:

- sodium bicarbonate
- sodium selenite
- [L-Ascorbic acid 2-phosphate sesquimagnesium salt hydrate Sigma](#)
- **Aldrich Catalog #A8960** Step 26
- DMEM/F12.
- Insulin solution (10 mg/ml)
- FGF2 (100 µg/ml)

SAFETY WARNINGS

For hazard information and safety warnings, please refer to the SDS (Safety Data Sheet).

Preparation of Recombinant Human TGF-beta 1 30m

1 [Albumin human Sigma](#)

Aldrich Catalog #A7223

**Most temperature sensitive, work fast

***This is for large number of aliquots. Change volume according to your own scale.

Order 1 mg but have them provide it in 2 x 500 µg aliquots plus a sample for testing if they will do this. 1000X = 1.745 µg/mL

- For 500 µg protein: bring to a total volume of 286.5 mL in buffer
- Use **[M]4 Milimolar (mM) HCl** containing **[M]1 mg/mL recombinant human albumin**
- Chilled Buffer (make about 300-500 mL of buffer)

2

Prepare buffer:

For 100 mL:

1. **34.4 µl HCl (11.6M)** (in fume hood cabinet) for final concentration of **[M]4 Milimolar (mM)**

2. Add **2 mL HSA (50 mg/mL)** for final **1 mg/mL** concentration
3. Bring vol up to 100 mL with molecular grade water

For 500 mL: **172 µl HCl** + **10 mL HSA** + **489.83 mL H₂O**

This can be made at any point and stored at **4 °C**.

30m

- 3 Bring Solution to **4 °C** **On ice** for **00:30:00** – do not need to do this if buffer already cold.

- 4 Set up tubes in hood – will need 2 hoods (over 500 aliquots).

- 5 To a 250 mL Millipore filter, add **100 mL chilled buffer**.



45m

Then add your **500 µg bottle of TGFβ from R&D Systems** (note the total volume – use pipet to measure) – thaw **On ice**, will take around **00:30:00** - **00:45:00**.

- 7 For TGFβ coming in the form of lyophilized powder, dissolve in the chilled buffer at **100 µg/ml**, then add to appropriate amount of buffer.

- 8 Bring to a total volume of 286.5 mL.

- 9 Filter.

- 10 Aliquot **500 µl** /tube (get some help because that's over 500 aliquots and you want to work quickly to get them frozen).

- 11 Label.

- 12 Freeze @ **-20 °C** for immediate use or **-80 °C** for long term storage.

Preparation of Holo-Transferrin

13 [holo-Transferrin human Sigma](#)
Aldrich Catalog #T0665-1G

Change scale accordingly if you are making a different amount.

- 1000X = 10.67 mg/mL à 426.8 mg (+/- 5mg) / 40 ml
- Use chilled PBS to dissolve this

14 Static zap your Falcon Tube and weigh paper before transferring the crystals.

15 On weigh paper, measure out **426.8 mg Holo-Transferrin** .

16 Add your **40 mL chilled PBS** .

17 Vortex gently until in solution, avoid foaming.

18 Filter with Steriflip.

19 Aliquot **500 µl** / tube.

20 Label.

21 Parafilm bottle of holo-transferrin when done.

The solution will be red in color.

- If you freeze it at **-80 °C** it will be orange in color
- If you freeze it at **-20 °C** it will turn clear
- If you freeze it at **-80 °C** and move it to **-20 °C** it will slowly turn from orange to clear
- The color is an indicator of the oxidation state, either way it's all right to use

Making E8 medium for ES/iPS culture

22 Dissolve **1.358 g sodium bicarbonate** in **50 mL water**. Warm the water at **37 °C** to completely dissolve the salt.

23 Dissolve **10 mg sodium selenite** in **1 mL water**. Then dilute this solution 1:100 with water to make **100 µg/ml** solution.

24 

Add **350 µl 100 µg/ml sodium selenite solution** to the **50 mL sodium bicarbonate solution**. Then filter sterilize.

25 Aliquot into 10 ml aliquots. Each aliquot is enough for 500 ml of DMEM/F12.

26 Dissolve **160 mg L-Ascorbic acid 2-phosphate sesquimagnesium salt hydrate** (**L-Ascorbic acid 2-phosphate sesquimagnesium salt hydrate Sigma Aldrich Catalog #A8960**) in **5 mL water**. Filter sterilize. Then aliquot into 1 ml aliquot. Each aliquot is enough for **500 mL DMEM/F12**.

27 Freeze all aliquots at **-80 °C**.

28 

To make E4 media. Add **10 mL sodium bicarbonate and sodium selenite solution** and **1 mL VitC solution** from above to **500 mL DMEM/F12**.

29 E4 media made this way should have osmolarity around 310 and pH around **pH7.4**.


30 

To make E6 media, add to the above E4 medium the following: **0.25 mL insulin solution (10 mg/ml)** and **0.5 mL transferrin (10 mg/ml)**. Mix well.

(Note: The original E8 formula contains **20 µg/ml insulin**. This recipe contains **5 µg/ml insulin**. Both should work well for ES/iPS culture).

31 

To make complete media (E8), add to the above E4 medium the following:

▣ **25 mL insulin solution (10 mg/ml)** , ▣ **0.5 mL FGF2 (100 µg/ml)** , ▣ **0.5 mL TGFb1 (1.7 µg/ml)** ,
▣ **0.5 mL transferrin (10 mg/ml)** . Mix well. Thaw all additives on cold-rack  **On ice** . Add cold base-media and wait until it is thawed, then add.

32 The complete medium is good for 2 weeks at  **4 °C** . **DO NOT** warm the media before feeding. Warming up the media will destabilize the growth factors in the media.

33 E4 and E6 media is good for at least 3 months at  **4 °C** , possibly up to a year.