




May 18, 2022

Immunophenotyping for NHPs, containment protocol

Jonathan Audet¹, Courtney Meilleur¹¹Public Health Agency of Canada

1

dx.doi.org/10.17504/protocols.io.x54v9j2b1g3e/v1 **Jonathan Audet**
Public Health Agency of Canada

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to protocols.io is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with protocols.io, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

This is a protocol used to perform immunophenotyping of whole blood (collected in EDTA tubes) or cells from bronchoalveolar lavage (BAL). We have successfully used this protocol for rhesus macaques, cynomolgus macaques, and African green monkeys (although antibody mix provided was titrated on AGM only). It allows the characterization of T cells (CD69/CD25 for activation; CCR7/CD45RA for memory phenotype; α/β vs γ/δ), B cells, Monocytes, Neutrophils, Basophils, Eosinophils, and potentially NK cells (CD56 does not work for AGM).

This protocol was designed to deal with samples coming from containment labs (> CL2) at our facility. If you are using this protocol at a different facility please ensure that proper testing and approvals are in place. The protocol can be used for experiments completed entirely in CL2, simply go from step 17 directly to 31.

DOI

dx.doi.org/10.17504/protocols.io.x54v9j2b1g3e/v1

Jonathan Audet, Courtney Meilleur 2022. Immunophenotyping for NHPs, containment protocol. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.x54v9j2b1g3e/v1>



flow cytometry, nonhuman primate, immunophenotyping, whole blood, broncho-alveolar lavage

_____ protocol ,

Dec 30, 2020

May 18, 2022

45922

For BAL fluid, make sure to set FSC voltage much lower than for whole blood.

 [Human TruStain](#)

FcX BioLegend Catalog #422302 Step 4

 [BD FACS Lysing Solution \(10X\) BD](#)

Biosciences Catalog #349202 Step 2

 [BD Cytofix/Cytoperm BD](#)

Biosciences Catalog #554722 In 2 steps

 [Ghost Dye Red 780 Tonbo](#)

Biosciences Catalog #13-0865-T500 In 2 steps

FACS Tube
Tube
Falcon 14-959-2A

2 ml screw-cap tubes
Microtubes
Sarstedt 72.694.006

 [PBS Contributed by users](#) In 3 steps

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to protocols.io is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company

nor any of the authors, contributors, administrators, or anyone else associated with protocols.io, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

The antibody mix described in the methods was tested on African green monkey whole blood and BAL fluid. CD56 does not stain AGM NK cells. All antibodies should cross-react with cynomolgus and rhesus macaques but the panel might need to be re-titrated.

Preparations 1h 10m

1 Prepare the staining mix. (for 1 sample:)

| A | B | C | D | E |
|-----------------|--------------------|------------------|--------------|-----------------|
| Supplier | Antibody | Clone | Channel | Volume per test |
| BD Biosciences | CD45 | D058-1283 | BUV395 | 1.25 |
| BD Biosciences | CD3 | SP34-2 | BUV496 | 5 |
| BD Biosciences | CD8 | RPA-T8 | BUV563 | 1.25 |
| BD Biosciences | CD16 | 3G8 | BUV737 | 5 |
| BD Biosciences | CD45RA | 5H9 | BV421 | 0.625 |
| BD Biosciences | CD49d | 9F10 | BV480 | 5 |
| BD Biosciences | CD4 | L200 | BV605 | 0.625 |
| BD Biosciences | CD14 | M5E2 | BV650 | 5 |
| BD Biosciences | CD123 | 7G3 | BV786 | 1.25 |
| BD Biosciences | CD25 | M-A251 | BB515 | 5 |
| Miltenyi Biotec | CD66abce | TET2 | PerCP-Vio700 | 1 |
| BD Biosciences | CD56 | MY31 or NCAM16.2 | BV711 | 5 |
| BD Biosciences | CD163 | GHI/61 | PE | 20 |
| BD Biosciences | CCR7 | 2-L1-A | PE-CF594 | 0.625 |
| BioLegend | HLA-DR | L243 | PE/Fire640 | 1.25 |
| BD Biosciences | CD69 | FN50 | PE-Cy7 | 5 |
| BD Biosciences | TCRgd | B1 | APC | 5 |
| BD Biosciences | CD20 | 2H7 | Alexa 700 | 0.625 |
| | | | | |
| | Total | | | 68.5 |
| | | | | |
| | BD Brilliant Stain | | | 31.5 |

- 2 Prepare the 1X FACS Lysing solution by diluting the 10X stock with Milli-Q water. You will need 2 ml of 1X solution for each sample.

[🔗BD FACS Lysing Solution \(10X\) BD](#)

Biosciences Catalog #349202

- 3 Prepare a 1:100 dilution of viability dye combining  **495 µL PBS** and  **5 µL**

[🔗Ghost Dye Red 780 Tonbo](#)

Biosciences Catalog #13-0865-T500

For whole blood, use the Ghost Dye undiluted.

Surface Staining 30m

- 4 Put 5 µl of TruStain FcX in FACS tubes (1 tube per sample).

[🔗Human TruStain](#)

FcX BioLegend Catalog #422302

FACS Tube

Tube

Falcon 14-959-2A

- 5 Add 100 µl of whole blood (EDTA blood) or BAL cells to the correct tube. Mix.

 **100 µL Sample**

- 6 Incubate 10 min at Room Temperature (RT). ⌚ **00:10:00** 🌡 **Room temperature**

10m

- 7 Add 5 µl of the diluted

[🔗Ghost Dye Red 780 Tonbo](#)

Biosciences Catalog #13-0865-T500



If staining whole blood: Add 5 µl of undiluted Ghost Dye.

- 8 Incubate 20 min at RT in the dark. 🌡 **Room temperature**

20m

⌚ **00:20:00**


9 Add stain mix. Mix.  100 µL


10 Incubate 20 min at RT in the dark.  Room temperature
 00:20:00

20m

RBC Lysis

11 Add 2 ml of 1X FACS Lysing solution. Vortex immediately, but gently.  2 mL


12 Incubate no more than 12 min at RT in the dark.  00:10:00 No more than 12 minutes^{10m}

13 Spin at 300 x g for 5 min.  300 x g, 20°C, 00:05:00

14 Decant supernatant.

15 Add 2 ml of PBS. Vortex.  2 mL

 PBS Contributed by users

16 Spin at 300 x g for 5 min.  300 x g, 20°C, 00:05:00

17 Decant supernatant.

Sample Inactivation 30m

18 Resuspend in Cytotfix/Cytoperm.

 [BD Cytofix/Cytoperm](#) **BD**

Biosciences Catalog #554722

18.1 100 µl per 5×10^5 cells.

18.2 Use at least 400 µl for easy decanting.

19 Incubate at least 30 min at RT in the dark. ⌚ 00:30:00 🌡 **Room temperature In the dark.**^{30m}

20 Spin at 500 x g for 8-10 min. ⌚ **500 x g, 20°C, 00:08:00**

21 On a clean bench, decant supernatant.

22 Use same volume of Cytofix/Cytoperm as before to resuspend the cells.

 [BD Cytofix/Cytoperm](#) **BD**

Biosciences Catalog #554722

23 Transfer in a 2 ml screwcap tube. Shake the tube to cover all surfaces with Cytofix/Cytoperm.

2 ml screw-cap tubes

Microtubes

Sarstedt 72.694.006

24 Transfer tubes from containment space to CL2 space according to the facility's approved protocols/SOPs.

- 25 Tubes can be opened in CL2 (in a BSC) no less than 30 min after the resuspension (step 17). (Samples are generally processed the next day; keep at 4 C overnight, in the dark)

Final Wash & Run 30m

- 26 Give tubes a quick spin in a tabletop centrifuge.
- 27 Ensure all tubes have at least a few hundred microliters of Cytotfix/Cytoperm.
- 28 Transfer the samples into 1 ml of PBS in FACS tubes. [🔗PBS Contributed by users](#)
- 29 Spin at 500 x g for 8-10 min. 🌀 **500 x g, 20°C, 00:08:00**
- 30 Decant supernatant.
- 31 Resuspend in 200 µl of PBS. [🔗PBS Contributed by users](#)
- 32 Run on FACSymphony A5.
- 33 Gating strategy for whole blood:

 [Gating_blood.pdf](#)