



Upload image

Jul 10, 2020

# Protocol for RPA-PCR Couple

Mustafa Munawar<sup>1</sup><sup>1</sup>University of Eastern Finland

1 Works for me dx.doi.org/10.17504/protocols.io.bifrkbm6

Mustafa Munawar

DOI

[dx.doi.org/10.17504/protocols.io.bifrkbm6](https://dx.doi.org/10.17504/protocols.io.bifrkbm6)

DOCUMENT CITATION

Mustafa Munawar 2020. Protocol for RPA-PCR Couple. **protocols.io**  
[dx.doi.org/10.17504/protocols.io.bifrkbm6](https://dx.doi.org/10.17504/protocols.io.bifrkbm6)

LICENSE

This is an open access document distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Jul 10, 2020

LAST MODIFIED

Jul 10, 2020

DOCUMENT INTEGER ID

39121

## Reagents

TwistAmp Liquid Basic (TwistDx Inc., United Kingdom)

40 mM dNTPs mix (total) (VWR Life Science)

RPA primers

Routine PCR tubes strips with cap strips

Routine PCR Reagents and primers/ probes

## Procedure

### Design RPA Assay

- 1) Pick several forward and reverse RPA primers of 30 to 35 nucleotides length from the region flanking your nucleic acid target.
- 2) Screen RPA primers to find optimal primer pair using TwistAmp Liquid Basic (TwistDx Inc.) and following related instructions from TwistDx manuals. Keep the nucleic acid template concentration high in the screening process.

### Design PCR Assay

- 3) Design and optimize a PCR assay from your nucleic acid target.

### Prepare RPA Reaction Mix (for slower reaction kinetics)

- 4) Prepare pre-master mix in routine 1.5 mL tube. The recipe for one reaction or 1X is following

2x Reaction Buffer = 25  $\mu$ L10x Basic E-mix = 2  $\mu$ L40 mM dNTPs (total) = 1  $\mu$ LOptimum forward RPA primer (10  $\mu$ M) = 1.5  $\mu$ LOptimum reverse RPA primer (10  $\mu$ M) = 1.5  $\mu$ L

Water = 15.21  $\mu$ L

- 5) Mix and spin the pre-master mix briefly.
- 6) Bring 20x Core Reaction Mix to room temperature and pipette mix it.
- 7) Add 20x Core Reaction Mix to the lid of pre-master mix vial and invert the vial 10 times to make the RPA master mix ready. Keep the volume of 20x Core Reaction Mix at only 1  $\mu$ L per 1X recipe.

#### Prepare PCR Reaction Mix

- 8) Prepare PCR reactions as per manufacturer recommendations.

#### Crude Maceration

- 9) Takes plant tissue weighing from 25-100 mg in round bottom 2 ml Eppendorf tube and initially macerate in 100  $\mu$ L water autoclaved distilled water with a plastic pestle.
- 10) Add more water to prepare suitable final dilution of macerate such as 1:10 or 1:20 (w/v) tissue weight by volume of water. Vortex and leave for particulate settlement for around one minute.

#### RPA-PCR couple

- 11) Pipette 47.21  $\mu$ L of the RPA master mix to tube strips and add 1  $\mu$ L macerate supernatant/template and 1.79  $\mu$ L of 280mM MgOAc in lid strips.
- 12) Start RPA reaction by inverting tubes (six inverts at least) and incubate tubes at 39  $^{\circ}$ C for 20 min, with an agitation step after four minutes. Incubation can be accomplished in PTC-0200 DNA Engine Cycler with no lid heating.
- 13) After incubation, heat reaction at 85  $^{\circ}$ C for 1 min to denature enzymes and then cool at 4  $^{\circ}$ C.
- 14) Transfer 1  $\mu$ L reaction directly to the PCR and store the rest of the RPA reaction at -20  $^{\circ}$ C for future use.
- 15) Run PCR in any PCR machine, keeping the number of cycles around 30.

#### **Equipment**

PTC-0200 DNA Engine Cycler or any other suitable small tubes incubator  
Mx3000P QPCR System (Agilent, Germany) or any PCR machine  
Routine centrifuges and vortexer