



Jan 12, 2021

Islet Embedding for Histology

IIDP-HIPP ¹¹Integrated Islet Distribution Program and Human Islet Phenotyping Program**1** Works for me dx.doi.org/10.17504/protocols.io.bqunmwve

Integrated Islet Distribution Program and Human Islet Phenotyping Program
Tech. support email: heather.durai@vumc.org



ABSTRACT

This Standard Operating Procedure (SOP) is based on the Vanderbilt Human Islet Phenotyping Program (HIPP) Immobilization of Islets in 3-D Gel for Histology. This SOP provides HIPP procedure for embedding islets in Collagen I gel for cryosectioning and histological staining.

This SOP defines the assay method used by the Human Islet Phenotyping Program (HIPP) for quantitative and qualitative determination of the Purified Human Pancreatic Islet product, post-shipment, manufactured for use in the National Institute of Diabetes and Digestive and Kidney Diseases (NIDDK)-sponsored research in the Integrated Islet Distribution Program (IIDP).

This Standard Operating Procedure (SOP) #: HIPP-06-v02

DOI

dx.doi.org/10.17504/protocols.io.bqunmwve

PROTOCOL CITATION

IIDP-HIPP 2021. Islet Embedding for Histology. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bqunmwve>

KEYWORDS

HIPP, Islet Embedding, cryosection, histological stain

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Dec 17, 2020

LAST MODIFIED

Jan 12, 2021

PROTOCOL INTEGER ID

45678

GUIDELINES


- **Integrated Islet Distribution Program (IIDP) (RRID:SCR_014387):** The IIDP is a grant funded program commissioned and funded by the NIDDK to provide quality human islets to the diabetes research community to advance scientific discoveries and translational medicine. The IIDP consists of the NIDDK Project Scientist and Program Official, the External Scientific Panel and the CC at City of Hope (COH). The IIDP CC integrates an interactive group of academic laboratories including the subcontracted IIDP centers.
- **IIDP Coordinating Center (CC):** Joyce Niland, Ph.D. and Carmella Evans-Molina, M.D., Ph.D. serve as Co-Principal Investigators (Co-PIs) for the IIDP Program located within the Department of Diabetes and Cancer Discovery Science at COH to coordinate the activities of the IIDP and Human Islet Phenotyping Program (HIPP). Dr. Niland, contact PI, oversees the daily activity of the IIDP staff, provides informatics/ biostatistical input, and subcontracts with the Islet Isolation Centers (IICs) to ensure the delivery of the highest quality human islets to IIDP approved investigators. Dr. Evans-Molina serves as the liaison to the HIPP, interacting closely to ensure that extensive, high quality phenotypic data are collected on islets distributed by the IICs. She also facilitates the delivery of this information to both the IICs and the IIDP-approved investigators, while responding to questions, issues, or suggestions for further HIPP enhancements.
- **Human Islet Phenotyping Program (HIPP):** The HIPP is a subcontracted entity of the IIDP through the COH and Vanderbilt University. The HIPP is directed by Marcela Brissova, Ph.D. and is responsible for performing specific standardized quality control assays agreed upon by both the IIDP and the HIPP, in order to provide enhanced, quality data on the human islets post-shipment, to the IIDP. The results of these assays will be approved by the CC and posted on the IIDP website for both the centers and the approved investigators.


References:


Dai C, Brissova M, Hang Y, Thompson C, Poffenberger G, Shostak A, et al. Islet-enriched gene expression and glucose-induced insulin secretion in human and mouse islets. *Diabetologia*. 2012 Mar;55(3):707–18. PMID: PMC3268985.
<https://pubmed.ncbi.nlm.nih.gov/22167125/>

Guo S, Dai, Guo M, Taylor B, Harmon JS, Sander M, et al. Inactivation of specific β cell transcription factors in type 2 diabetes. *J Clin Invest*. 2013 Aug;123(8):3305–16. PMID: PMC3726150.
<https://pubmed.ncbi.nlm.nih.gov/23863625/>

MATERIALS TEXT

1. PBS (phosphate buffered saline) with no Ca/Mg, 1X (Invitrogen 14190-144)
2. Collagen I stock (BD 354249)
3. DMEM, 5X (Dulbecco modified Eagle's minimal essential medium, Gibco 31600-034)
4. HEPES (N-2-hydroxyethylpiperazine-N-2-ethane sulfonic acid, Sigma H0887)
5. NaHCO₃ (sodium bicarbonate, Sigma S8761)
6. Sucrose (Fisher Scientific BP220-1)
7. 16% Paraformaldehyde (Electron Microscopy Science 15710)
8.  **1.5 mL** Centrifuge Tube (Fisher Scientific 05-408-129)
9. P-200 Pipet Tips, Sterile (Fisher Scientific P-2069)
10. Wide-bore pipet tips for transferring islets (Fisher Scientific 212362C)
11. 96-well Plate (Falcon 353072)

12. Needle, 28G, insulin syringe (Exel International 26027)
13. 12-well plate (Fisher Scientific 351143)
14. Cryomolds (Tissue-Tek 4557)
15. OCT (Optimum Cutting Temperature compound, Tissue-Tek 4583)
16. Plus Gold Slides (Fisher Scientific 15-188-48)
17. Cryostat (Leica)
18. Adjustable tilt rocker (LabNet)
19. Stereomicroscope (Olympus)
20. Digital Camera
21.  **-80 °C** Freezer (Thermo Scientific)

 [PBS \(phosphate buffered saline\) with no Ca/Mg](#) **Thermo Fisher**

Scientific Catalog #Invitrogen 14190-144

 [Collagen I](#)

stock Corning Catalog #BD 354249

 [DMEM 5X](#) **Thermo Fisher**

Scientific Catalog #Gibco 31600-034

 [HEPES](#) **Sigma**


Aldrich Catalog #Sigma H0887

 [NaHCO3](#) **Sigma**

Aldrich Catalog #Sigma S8761

 [Sucrose](#) **Fisher**

Scientific Catalog #BP220-1

 [16% Paraformaldehyde](#) **Electron Microscopy**

Sciences Catalog #15710

Centrifuge Tube

1.5mL

Fisher Scientific

05-408-129



Pipet Tips

P-200

Sterile

Fisher Scientific P-2069



Wide-bore pipet tips

Fisher Scientific

212362C



96-well Plate

Falcon

353072



Needle, insulin syringe

28G

Exel International

26027



12-well plate

Falcon


351143



Cryomolds

Tissue-Tek

4557



OCT (Optimum Cutting Temperature compound)

Tissue-Tek

4583



Plus Gold Slides

Fisher Scientific

15-188-48



Cryostat

Leica

Cryostat



Adjustable tilt rocker

LabNet

35 Deluxe Rocking Platform



Stereomicroscope

Olympus

Stereomicroscope



Freezer

-80°C

Thermo
Scientific

ultra-low freezers
(ULT)



SAFETY WARNINGS

Sucrose (*Fisher Scientific BP220-1*)



May cause skin, eye, and respiratory tract irritation. Wear personal protective equipment. Ensure adequate ventilation. Avoid contact with skin, eyes and clothing. Avoid ingestion and inhalation. Use personal protective equipment. Wash off immediately with plenty of water for at least 15 minutes. Get medical attention immediately if symptoms occur. Rinse immediately with plenty of water, also under the eyelids, for at least 15 minutes. Obtain medical attention. Move to fresh air. If breathing is difficult, give oxygen. If breathing is difficult, give oxygen. Get medical attention immediately if symptoms occur.

16% Paraformaldehyde (*Electron Microscopy Science 15710*)



16% Paraformaldehyde 15710.pdf

ABSTRACT

This Standard Operating Procedure (SOP) is based on the Vanderbilt Human Islet Phenotyping Program (HIPP) Immobilization of Islets in 3-D Gel for Histology. This SOP provides HIPP procedure for embedding islets in Collagen I gel for cryosectioning and histological staining.

This SOP defines the assay method used by the Human Islet Phenotyping Program (HIPP) for quantitative and qualitative determination of the Purified Human Pancreatic Islet product, post-shipment, manufactured for use in the National Institute of Diabetes and Digestive and Kidney Diseases (NIDDK)-sponsored research in the Integrated Islet Distribution Program (IIDP).

This Standard Operating Procedure (SOP) #: HIPP-06-v02


Procedures

1 Islet Immobilization and Fixation

- 1.1 Prepare Collagen I working solution (1 mL) by combining 375 µl Collagen I stock, 355 µl sterile water, 20 µl HEPES, 50 µl NaHCO₃, 200 µl 5X DMEM.
- 1.2 Using a P-1000 pipette and wide-bore tips, transfer an aliquot of islet suspension containing approximately 500 IEQs into a 1.5 mL centrifuge tube and centrifuge at 1000 rpm for 1 min. Aspirate the supernatant. If islet medium is rich in serum, wash the islets once with 1X PBS.
- 1.3 Add 150 µl of Collagen I working solution to the islet pellet and transfer the mixture into a 96-well plate. Place the plate in a tissue culture incubator set to 37 °C . Incubate for 90 min and Collagen I will form into a gel. While gel is solidifying, islets will settle by gravity at the bottom of the gel and form a monolayer.
- 1.4 Place 96-well plate on ice and add ice-cold 4% paraformaldehyde/1X PBS on top of the gel to fill up the well. Fix for 10 min on ice.
- 1.5 Using a 28 G needle (0.3 mL insulin syringe) loosen up the gel containing islets from the sides of the well. Then, using a fine spatula transfer the gel into a 12-well plate containing 3 mL of ice-cold 4% paraformaldehyde/1X PBS per well. ***Be very careful while you do this because the gel is relatively soft.*** Fix the islets on ice for an additional 15 min under very mild agitation using an adjustable tilt rocker (low setting).
- 1.6 Gently aspirate the fixative and wash the gel with 3 mL of 1X PBS on ice for 20 min under very mild agitation using an adjustable tilt rocker (low setting).
- 1.7 Repeat step 1.6 two times.

2 Islet Embedding and Sectioning

- 2.1 After the last wash, add 3 mL of 30% sucrose/1X PBS. Allow the gel to equilibrate for 2-3 hours or overnight at 4 °C .
- 2.2 Lift the gel from the dish using a fine spatula and place it into a standard cryomold that is half-way filled with OCT compound. While doing this, orient the gel such that **the bottom of the gel (where the islets are) is facing the bottom of the cryomold**. Add more OCT to fill up the cryomold. Take an image of the gel on an Olympus SZX12 microscope at 16x. Example of islets embedded in Collagen I is shown in *Figure 1*.



- 2.3 Place the embedded gel in a -80°C freezer. When OCT compound is frozen, wrap the block in pre-labeled aluminum foil, place it a zip-lock bag, seal and store at -80°C .
- 2.4 Cut 8- μm cryosections and plate two sections per slide on *Plus Gold* slides. Transfer slides and blocks on dry ice back to the HIPP and store in a -80°C freezer.
- 2.5 Enter slides and blocks into the Vanderbilt HIPP electronic inventory system and upload images of islet gels to the IIDP-HIPP database.

Deviations and Resolution

- 3 Document any deviations that occurred during this protocol that affect the final results and report with the analysis of the assay.