

Jun 27, 2024

# Mutant generation in Streptococcus mitis strain B6

DOI

#### dx.doi.org/10.17504/protocols.io.261ge5ezyg47/v1

Samantha King<sup>1</sup>

<sup>1</sup>Nationwide Children's Hospital

King lab



### Samantha King

Nationwide Children's Hospital

# OPEN ACCESS



DOI: dx.doi.org/10.17504/protocols.io.261ge5ezyg47/v1

Protocol Citation: Samantha King 2024. Mutant generation in Streptococcus mitis strain B6. protocols.io https://dx.doi.org/10.17504/protocols.io.261ge5ezyg47/v1

License: This is an open access protocol distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working We use this protocol and it's

working

**Created:** June 26, 2024

Last Modified: June 27, 2024

Protocol Integer ID: 102478

Keywords: Transformation, mutagenesis, Streptococcus mitis B6

**Funders Acknowledgement: American Heart Association** Grant ID: 19TPA34760750



#### Abstract

This protocol is the methodology that we have successfully employed to generate and confirm insertion deletion mutants in Streptococcus mitis strain B6. Attached to the protocol is a file that includes primers for the MonX mutation

## **Attachments**



**B6** mutagenesis of Mo...

14KB

#### **Materials**

#### **Table 1: Primers Used in This Study**

A	В	С	D
Target	Name	Sequence 5' to 3'	Location (accession no.)
Spec	S1	CGATTTTCGTTCGTGAATAC	5418-5399 (KM009065)
	S2	TATGCAAGGGTTTATTGTTTTC	4265-4286 (KM009065)
pJet	P1	CGACTCACTATAGGGAGAGCGGC	310-332 (EF694056.1)
	P2	AAGAACATCGATTTTCCATGGCAG	405-428 (EF694056.1)

## Safety warnings



Appropriate biosafety procedures need to be followed



## Streptococcus mitis B6 mutagenesis protocol

1

#### Creation of a plasmid construct

- 1.1 The regions upstream and downstream of the region to be deleted were amplified using primers 1 and 2, and 3 and 4, respectively. These primers were designed to contain appropriate overhangs to allow In-fusion with pJET 1.2/blunt and the antibiotic resistance cassette.
- 1.2 The spectinomycin resistance cassette (*aad*9) was amplified using primers S1 and S2
- 1.3 PCR products were then purified with a Qiagen PCR Purification Kit
- 1.4 The three fragments for each mutant construct were cloned into pJET 1.2/Blunt (Thermo Fisher Scientific) using In-Fusion Snap Assembly (Takara) and transformed into Escherichia coli stellar competent cells.
- 1.5 Transformants were selected on LB agar plates supplemented with ampicillin (100  $\mu$ g/ml) and incubated at 37°C overnight
- 1.6 The resulting colonies were confirmed as ampicillin resistant by streaking on a new LB agar plate supplemented with ampicillin (100 μg/ml).
- 1.7 Transformants were screened by colony PCR using pJET 1.2/blunt Fwd and Rev primers.
- 1.8 For transformants giving an appropriate PCR product, a 5 mL LB culture supplemented with ampicillin (100  $\mu$ g/ml) was grown overnight at 37 °C with shaking at 200 rpm.
- 1.9 The plasmid was then purified using Qiagen Miniprep Kit and confirmed by sequencing
- 2 Transformation of *S. mitis* B6
- 2.1 Strains were growth at 37°C in C+Y pH8 [1] starting at a low inoculum i.e. from a plate or diluting from a culture 1:100 (starting optical density at 600nm  $[OD_{600}] = 0.03$  to 0.05).



- 2.2 When the culture was close to  $OD_{600}$  =0.1, 950  $\mu$ l of C+Y pH 8.0 medium was added to a 1.5 ml tube with 10 µl of 100 mM CaCl2, 2 µl of competence stimulating peptide (CSP) (EMRKPDGALFNLFRRR - 1 mg/ml), and 150 ng of DNA. A no DNA control tube is included to account for potential contamination. These tubes were prewarmed in a waterbath to 37°C.
- 2.3 When the culture reached an  $OD_{600}$  of 0.10 to 0.12, 50  $\mu$ l of culture was added to the prewarmed tubes.
- 2.4 Tubes were incubated in a waterbath at 37 °C for 2 hr.
- 2.5 Reactions were pelleted by centrifugation and resuspended in approximately 100 µl of media. This was plated on selective Tryptic Soy Agar (TSA) plates spread with 5000 U catalase (Worthington Biochemical Corporation).
- 2.6 Plates were incubated at 37oC in 5% CO2 overnight and then patched onto selective plates.
- 3 **Confirmation of putative transformants**
- 3.1 Putative transformants were grown in tryptic soy broth and DNA prepared as previously described [2].
- 3.2 The mutations were confirmed by PCR and sequencing (using primers S1 and S2) or genome sequencing.

#### Protocol references

- LACKS S, HOTCHKISS RD. A study of the 1. genetic material determining an enzyme in Pneumococcus. Biochim Biophys Acta. 1960;39:508-18. doi: 10.1016/0006-3002(60)90205-5. PubMed PMID: 14413322.
- Gaytán MO, Singh AK, Woodiga SA, Patel 2. SA, An SS, Vera-Ponce de León A, et al. A novel sialic acid-binding adhesin present in multiple species contributes to the pathogenesis of Infective endocarditis. PLoS Pathog. 2021;17(1):e1009222. Epub 20210119. doi: 10.1371/journal.ppat.1009222. PubMed PMID: 33465168; PubMed Central PMCID: PMCPMC7846122.