



Sep 28, 2022

# Molecular Cloning

Lizhen Zhu<sup>1</sup>, Ruixin Lin<sup>2</sup>, Xiaoqi Ye<sup>2</sup>, Ruiyan zeng<sup>2</sup>, 18814306949 Minkun Huang<sup>2</sup>, 624787660<sup>2</sup>

<sup>1</sup>Shenzhen Polytechnic; <sup>2</sup>ShenZhen Polytechnic

1 Works for me

Share

[dx.doi.org/10.17504/protocols.io.e6nvw59o9vmk/v1](https://dx.doi.org/10.17504/protocols.io.e6nvw59o9vmk/v1)

Ruixin Lin  
ShenZhen Polytechnic

## ABSTRACT

SZPT-CHINA Experimental methods of molecular cloning

## DOI

[dx.doi.org/10.17504/protocols.io.e6nvw59o9vmk/v1](https://dx.doi.org/10.17504/protocols.io.e6nvw59o9vmk/v1)

## PROTOCOL CITATION

Lizhen Zhu, Ruixin Lin, Xiaoqi Ye, Ruiyan zeng, 18814306949 Minkun Huang, 624787660 2022. Molecular Cloning. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.e6nvw59o9vmk/v1>

## LICENSE

———— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Oct 08, 2020

## LAST MODIFIED

Sep 28, 2022

## PROTOCOL INTEGER ID

42914

## GUIDELINES

advisers WeiJia Wu  
LiJun Zhang (PhD)  
Yongjun Tang(PhD)

## MATERIALS TEXT

### MATERIALS

☒ Q5 High-Fidelity 2X Master Mix - 500 rxns **New England**

**Biolabs Catalog #M0492L**

☒ Nuclease-free Water **New England**

**Biolabs Catalog #E6317**

☒ NEBNext Second Strand Synthesis Enzyme Mix **New England**

**Biolabs Catalog #E7425**

☒ Nuclease-free Water **New England**







**Biolabs Catalog #E7764**

☒ T4 DNA Ligase Buffer (10X) **Thermo**

**Fisher Catalog #B69**

☒ 2×PCR Master Mix **Solarbio**

### Vector linearization and target fragment preparation

- 1 Q5 DNA polymerase PCR system(  50  $\mu\text{L}$  )
  -  25  $\mu\text{L}$  Q5 High-Fidelity 2X Master Mix
  -  2  $\mu\text{L}$  10 $\mu\text{L}$  Forward Primer
  -  2  $\mu\text{L}$  10  $\mu\text{M}$  Reverse Primer
  -  1  $\mu\text{L}$  Template DNA
  -  20  $\mu\text{L}$  Nuclease-Free Water







After mixing gently, centrifuge to collect all the liquid to the bottom.

### PCR Program

STEP	TEMP	TIME
Initial Denaturation	98°C	30 seconds
Denaturation	98°C	5 seconds
Annealing	*50–72°C	10 seconds
Extension	72°C	15 seconds/kb
25–35 Cycles		
Final Extension	72°C	2 minutes
Hold	4–10°C	

### Enzyme digestion Connection

- 2 Enzyme digestion Connection

- 2.1 Enzyme digestion and connection system(  10  $\mu\text{L}$  )
  -  3  $\mu\text{L}$  Vector
  -  1  $\mu\text{L}$  Insert
  -  1  $\mu\text{L}$  T4 DNA Ligase Buffer (10X)
  -  0.5  $\mu\text{L}$  NEB Golden Gate Enzyme Mix (BsaI-HFv2)
  -  4.5  $\mu\text{L}$  Nuclease-free H<sub>2</sub>O

## 2.2 Connection system

INSERT NUMBER	SUGGESTED ASSEMBLY PROTOCOL
For 1 Insert	37°C, 5 min (cloning) or 37°C, 1 hr (library preparation) ® 60°C, 5 min
For 2–10 Inserts	(37°C, 1 min ® 16°C, 1 min) x 30 ® 60°C, 5 min
For 11–20+ Inserts	(37°C, 5 min ® 16°C, 5 min) x 30 ® 60°C, 5 min

## Transformation

### 3 Transformation

- 3.1 Take 50  $\mu\text{L}$  of competent cells in a sterile Eppendorf tube.
- 3.2 Add 2~5  $\mu\text{L}$  of the plasmid to be transformed into each tube, gently mix the plasmid and competent cells with a pipette tip, and ice bath for 25 minutes.
- 3.3 Heat shock the ice-bath mixture in a 42°C water bath for 45s, do not shake the centrifuge tube.
- 3.4 Take out the ice bath and cool for 2 minutes.
- 3.5 Add 800  $\mu\text{L}$  of non-resistant LB medium preheated to 37°C to each tube, shake gently at 37°C, 150r/min for 1h.
- 3.6 Centrifugally concentrate to 100 $\mu\text{L}$ , spread it on an LB plate containing antibiotics, and incubate inverted at 37°C for 12-16h.

### 3.7 Make a negative control at the same time (use distilled water instead of plasmid).

#### Colony PCR screening of positive clones

#### 4 Colony PCR screening of positive clones

4.1 Add 25ul of sterile water to the 96-well plate, pick about 15-20 colonies from the successfully transformed DH5α bacterial plate with a pipette tip, and mix them evenly. (Make a corresponding mark on the picking plate).

4.2 Put the mixed bacteria solution into 10ul to 200ul PCR tube and mark it.

4.3 Take 0.1ul of bacterial solution in a 96-well plate and add it to the antibiotic-containing plate with a pipette tip and mark it.

4.4 After mixing, the bacteria solution is used for PCR, the system is as follows(  20 μL )

 5 μL 2×PCR Mix

 0.5 μL 10μM Forward Primer

 0.5 μL 10 μM Reverse Primer

 2 μL Bacteria solution

 2 μL dd H2O

The procedure is as follows

STEP	TEMP	TIME
Initial Denaturation	94°C	30 seconds
Denaturation	94°C	3min
Annealing	*55-60°C	30seconds
Extension	72°C	1min/kb
25–35 Cycles		
Final Extension	72°C	2 minutes
Hold	4–10°C	

4.5 Verification result of 1% agarose gel electrophoresis.

