



Apr 22, 2020

HubMap UF TMC Tissue Dissociation to Single Cell

Marda Jorgensen¹, Maigan Brusko¹

¹University of Florida

1 Works for me dx.doi.org/10.17504/protocols.io.bd9vi966

Human BioMolecular Atlas Program (HuBMAP) Method Development Community

GUIDELINES

Perform all steps using appropriate aseptic technique.

MATERIALS

NAME	CATALOG #	VENDOR
Forceps	F3001	Gold Biotechnology
Scissors	HCT7.1	Carl Roth
DMEM	11885	Gibco - Thermo Fischer
Sterile deionized H2O		
10X Collagenase/Hyaluronidase in DMEM		Stemcell Technologies
Miltenyi C Tube		
Sterile Petri Dishes		
Disposable Scalpels 10 Blade		
FBS USDA Origin	25-550	
Penicillin/Streptomycin	30-002-CL	Corning
0.5M EDTA	BM-150	Boston Bioproducts
Ammonium Chloride	A661-500	Fisher Scientific
NaHCO	S233-50	Fisher Scientific
Disposable Filter Flask with PES Membrane		
gentleMACS C-Tube	130-093-334	Miltenyi Biotec
gentleMACS Dissociator	130-093-235	Miltenyi Biotec
500um Cell Strainer	43-50500-01	pluriSelect
100um Cell Strainer	22363549	Fisher Scientific
40um Cell Strainer	22363547	Fisher Scientific
15ml Conical Tubes	21103	Olympus
Nexcelom Cellometer		
Nexcelom Cellometer Slides		

Citation: Marda Jorgensen, Maigan Brusko (04/22/2020). HubMap UF TMC Tissue Dissociation to Single Cell. <https://dx.doi.org/10.17504/protocols.io.bd9vi966>

This is an open access protocol distributed under the terms of the **Creative Commons Attribution License** (<https://creativecommons.org/licenses/by/4.0/>), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

NAME ▾	CATALOG # ▾	VENDOR ▾
Nexcelom ViaStain AO/PI		
1mL Cryovials Internally Threaded with O-ring	24-200P	Olympus
Cryostar CS10	C2874	Sigma Aldrich
Cool Cells	BCS405G	biocision
LabXpert Label Maker	060900	Brady Worldwide, Inc.
LabXpert Label Maker Stickers	X-120-492	Brady Worldwide, Inc.

SAFETY WARNINGS

Use universal safety precautions when handling human samples and employ personal protective equipment (e.g., face mask with shield, gloves, lab coat or apron).

Preparation

- 1 To make ACK Lysis Buffer:
 - 1.1 Add 8.02 grams of NH_4Cl , 0.84 grams of NaHCO_3 , and 5uL of 0.5M EDTA into 1000mL of deionized water.
 - 1.2 Filter buffer through 22um PES Filter. Aliquot in 50mL conicals and store at room temperature for long-term storage.
- 2 To make the cDMEM:
 - 2.1 Add 50mL of FBS, 10mL of Penicillin/Streptomycin to 440mL of DMEM.
 - 2.2 Filter through a filter flask and store at 4°C.

Procedure

- 3 Obtain 1cm x 1cm x 1cm tissue cube (thymus and spleen) or whole/partial lymph node from organ dissection protocol. Place in 15mL conical tube with 5mL cDMEM, retain on ice.
- 4 Prepare 10mL digest mix using 1mL of 10x Collagenase/Hyaluronidase in DMEM and 9mL dMEM with 5% FBS. Warm to 37°C in water bath.
- 5 Mince tissue crosswise with two scalpels in sterile petri dish until tissue fragments are approximately 2 mm in all dimensions.
- 6 Transfer tissue fragments to C tube, incubate at 37°C in shaking incubator for 30 minutes.
- 7 Place C tube on gentleMACS dissociator, run protocol mspIn0201.
- 8 Transfer cell suspension to 50 ml conical over 500 uM mesh filter and wash over filter with DMEM until tube is full. Centrifuge 350 x g for 10 minutes beginning at RT, step down to 8°C.
- 9 Discard supernatant. Resuspend cell pellet in cold cDMEM, decant over 100 uM filter into new conical tube. Centrifuge 350 x g for 10 minutes, 8°C.

- 10 Discard supernatant and retain cell pellet. Resuspend in 10 ml cold HBSS, decant over 40 μ M filter into new conical tube. Retain cell suspension on ice while performing cell count.
- 11 Take a 20ul aliquot for counting on a Nexcelom Cellometer.
- 12 Add 20ul of Viastain AO/PI to aliquot and mix thoroughly with pipette. Add 20ul to counting slide and count with "immune cells low RBC" with dilution factor 2.0. Record total yield and viability in case worksheet.
- 13 Centrifuge remaining cell suspension at 350 x g for 10 minutes at 8°C. Decant supernatant and dissociate cell pellet. Resuspend dropwise in Cryostor 10 at 25 million cells/ml.
- 14 Aliquot suspension to labeled cryovial and place in CoolCell. Transfer to -80°C freezer for 18-24 hours prior to transfer to liquid nitrogen cryounit.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited