



FEB 01, 2024

Purification of GST-LC3B

Elias Adriaenssens¹

¹Sascha Martens lab, University of Vienna, Max Perutz Labs - Vienna

ASAP Collaborative Research Network



Elias Adriaenssens

Sascha Martens lab, University of Vienna, Max Perutz Labs - ...

ABSTRACT

This protocol describes the purification of GST-LC3B.

ATTACHMENTS

[ppgyb32rp.pdf](#)

OPEN ACCESS



DOI:

dx.doi.org/10.17504/protocols.io.3byl4qnbjvo5/v1

Protocol Citation: Elias Adriaenssens 2024. Purification of GST-LC3B. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.3byl4qnbjvo5/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

Created: Jan 31, 2024

Last Modified: Feb 01, 2024

PROTOCOL integer ID: 94565

MATERIALS

Keywords: ASAPCRN

Funders Acknowledgement:

Aligning Science Across
Parkinson's (ASAP)
Grant ID: ASAP-000350
Marie Skłodowska-Curie MSCA
Postdoctoral fellowship
Grant ID: 101062916

Materials:

- pGEX-4T1 vector (LC3B-encoding vector is available from Addgene)
- Glutathione Sepharose 4B beads (GE Healthcare)
- Amicon filter (Merck Millipore)
- Superdex 75 16/600 column (Cytiva)

Lysis buffer:

A	B
HEPES pH 7.5	50 mM
NaCl	300 mM
MgCl ₂	2 mM
βmercaptoethanol	2 mM
cOmplete EDTA-free protease inhibitors (Roche)	
DNase (Sigma)	

Wash buffer:

A	B
HEPES pH 7.4	50 mM
NaCl	300 mM
DTT	1 mM

High-salt wash buffer:

A	B
HEPES pH 7.4	50 mM
NaCl	700 mM
DTT	1 mM

SEC buffer:





A	B
HEPES pH 7.4	25 mM
NaCl	150 mM
DTT	1 mM

Equipment

Purification of GST-LC3B

8h 30m 30s

1 To purify GST-LC3B, insert human LC3B cDNA in a pGEX-4T1 vector. Note, the last five amino acids of LC3B are deleted to mimic the cleavage by ATG4.

2 After the transformation of the pGEX-4T1 vector encoding GST-LC3B in E.coli Rosetta (DE3) pLysS cells, grow the cells in LB medium at  37 °C until an OD₆₀₀ of 0.8-1 and induce the protein expression with  1 millimolar (mM) IPTG for  04:00:00 at  37 °C .

4h

3 Collect cells by centrifugation and resuspend in lysis buffer.




Lysis buffer:


A	B
HEPES pH 7.5	50 mM
NaCl	300 mM
MgCl ₂	2 mM
βmercaptoethanol	2 mM
cOmplete EDTA-free protease inhibitors (Roche)	
DNase (Sigma)	

4 Sonicate cell lysates.

1m

4.1 Sonicate cell lysates twice for  00:00:30 (1/2).

30s

4.2 Sonicate cell lysates twice for  00:00:30 (2/2).

30s

5 Clear lysates by centrifugation at 140000 x g, 4°C, 00:30:00 in a Beckman Ti45 rotor.

30m



6 Collect and incubate the supernatant with pre-equilibrated Glutathione Sepharose 4B beads (GE Healthcare) for 02:00:00 at 4 °C with gentle shaking to bind GST-LC3B.

2h



7 Centrifuge the samples to pellet the beads and remove the unbound lysate.



8 Wash beads twice with wash buffer, once with high salt wash buffer and two more times with wash buffer.



Wash buffer

A	B
HEPES pH 7.4	50 mM
NaCl	300 mM
DTT	1 mM

High-salt wash buffer

A	B
HEPES pH 7.4	50 mM
NaCl	700 mM
DTT	1 mM

9 Elute the proteins Overnight with 20 millimolar (mM) reduced L-glutathione in

8h



50 millimolar (mM) HEPES pH 7.4, 300 millimolar (mM) NaCl, 1 millimolar (mM) DTT buffer.

10 Collect the supernatant, filter through a 0.45 µm syringe filter, concentrate using a 10 kDa cut-off Amicon filter (Merck Millipore), and load onto a pre-equilibrated Superdex 75 16/600 column (Cytiva).

11 Elute the proteins with SEC buffer.


SEC buffer

A	B
HEPES pH 7.4	25 mM
NaCl	150 mM
DTT	1 mM

12 Analyze fractions by SDS-PAGE and Coomassie staining.



13 Pool fractions containing purified GST-LC3B protein.

14 After concentrating the purified protein, aliquot the protein and snap-frozen it in liquid nitrogen. Store proteins at  -80 °C .

