

Sep 13, 2024

RNAi induction

 In 1 collection

DOI

dx.doi.org/10.17504/protocols.io.5qpvobb7xl4o/v1

Ben Jenkins¹

¹University of Oxford



Ben Jenkins

University of Oxford

OPEN  ACCESS



DOI: **dx.doi.org/10.17504/protocols.io.5qpvobb7xl4o/v1**

Protocol Citation: Ben Jenkins 2024. RNAi induction. protocols.io **<https://dx.doi.org/10.17504/protocols.io.5qpvobb7xl4o/v1>**

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

Created: February 07, 2022

Last Modified: September 13, 2024

Protocol Integer ID: 57888

Abstract

Back-dilution, induction, and creating frozen *E. coli* feeding stocks for downstream RNAi



Preparing stocks

15m

- 1 In a sterile hood, streak out HT115 *E. coli* frozen stock onto LB Agar plate (containing $[M]$ 0.05 mg/mL Ampicillin and $[M]$ 0.0125 mg/mL Tetracycline)

Ampicillin selects for the L4440 plasmid containing the knock-down construct
In a sterile hood, streak out HT115 *E. coli* strain (deficient in RNase III)

- 2 Incubate overnight at 37°C

Set up

15m

- 3 In a sterile hood, pick a single colony to inoculate 10 mL LB containing $10\text{ }\mu\text{L}$ Ampicillin and $25\text{ }\mu\text{L}$ Tetracycline.

Final concentrations Ampicillin ($[M]$ 0.05 mg/mL), and Tetracycline ($[M]$ 0.0125 mg/mL)

- 3.1 Streak used inoculation loop onto a fresh LB Agar plate (containing $[M]$ 0.05 mg/mL Ampicillin and $[M]$ 0.0125 mg/mL Tetracycline) to create a back-up stock plate

Incubate LB Agar plate overnight at 37°C

- 4 Incubate LB overnight with shaking at 180 rpm , 37°C

Induction

5h 30m

- 5 In a sterile hood, add $600\text{ }\mu\text{L}$ of overnight *E. coli* pre-culture to 15 mL LB containing $15\text{ }\mu\text{L}$ Ampicillin and $37.5\text{ }\mu\text{L}$ Tetracycline.

Final concentrations Ampicillin ($[M]$ 0.05 mg/mL), and Tetracycline ($[M]$ 0.0125 mg/mL)

Conduct remainder of experiment under sterile conditions


- 6 Incubate with shaking at 180 rpm , 37°C , 02:30:00

2h 30m





- 7 Add $60\text{ }\mu\text{L}$ ($[M]$ 0.4 millimolar (mM)) IPTG






This induces template expression within the L4440 plasmid construct

8 Incubate with shaking at  180 rpm, 37°C, 03:00:00


3h

9 In a sterile hood, prepare  150 mL NCL containing  300 μ L Ampicillin,  600 μ L IPTG, and  30 μ L β -sitosterol


Ampicillin ( 0.1 mg/mL), IPTG ( 0.4 millimolar (mM)), and β -sitosterol ( 0.0008 mg/mL)


Match OD and freeze

11m


10 Spin down *E. coli* at  3200 x g, 00:02:00

2m



11 Remove supernatant, and add  20 mL MQ

12 Spin down *E. coli* at  3200 x g, 00:02:00

2m

13 Resuspend with  3 mL NCL


NCL (containing Ampicillin, IPTG, and β -sitosterol) prepared earlier

14 Add  800 μ L of NCL and  200 μ L resuspended *E. coli* to a cuvette, and measure OD₆₀₀

Use OD dilution calculator:

<https://docs.google.com/spreadsheets/d/1hJUsvOjwcuHkSeFQaAbNhn1ZXleTpNXN8RAR4gZzWDE/edit?usp=sharing>


Make sure to pipette *E. coli* up and down 3 times before removing from falcon, and again after adding to cuvette for OD measurement


15 In a sterile hood, add required volumes of resuspended *E. coli* and NCL to give  1.5 mL of *E. coli* at OD 3

Use OD dilution calculator as above




Larger experiments may require greater volumes, scale up your stock accordingly
Conduct remainder of the protocol under sterile conditions

16 Add  37.5 µL of 80% glycerol per 1.5mL *E. coli* (to give 2%), and vortex

17 Add  60 µL of *E. coli* per well to a 96-well plate

Ensure *E. coli* are in a consistent layout for every 96-well plate
This is enough for 20 aliquots per *E. coli* stock
Do not use the outside wells

18 Add lids to each stock plate, and store at  -20 °C until required