

Jul 04, 2024

Protein Aggregation Capture

This protocol is a draft, published without a DOI.

Mark Kinnin¹

¹Queen's University Belfast



Mark Kinnin

Queen's University Belfast

OPEN  ACCESS



Protocol Citation: Mark Kinnin 2024. Protein Aggregation Capture. protocols.io <https://protocols.io/view/protein-aggregation-capture-dgpg3vjw>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: In development

We are still developing and optimizing this protocol

Created: July 04, 2024

Last Modified: July 04, 2024

Protocol Integer ID: 102856

Abstract

PAC Adapted from Batth et al., 2019



Protocol materials



Tris **P212121** Step 10



TFA **Thermo Scientific** In [2 steps](#)



MagReSyn Hydroxyl Beads **Resyn Biosciences** Step 1



Acetonitrile Step 5



70% ACN **Thermo Scientific** In [2 steps](#)



100% ACN **Thermo Scientific** Step 8

















70% Ethanol Step 9



Trypsin **P212121 Catalog #RP-T70010** Step 11



Equilibration

- 1 Add  5 μL of  MagReSyn Hydroxyl Beads **Resyn Biosciences** to 1.5mL LoBind Eppendorf tube
 - 1.1 Place tube on magnet and allow Microparticles to clear  00:00:10 10s
 - 1.2 Remove storage solution and discard
- 2 Equilibrate Microparticles by adding  100 μL  70% ACN **Thermo Scientific** and mix by gentle agitation
 - 2.1 Place tube on magnet and allow Microparticles to clear  00:00:10 10s
 - 2.2 Remove supernatant and discard
- 3 Add  100 μL  70% ACN **Thermo Scientific** and mix by gentle agitation
 - 3.1 Place tube on magnet and allow Microparticles to clear  00:00:10 10s
 - 3.2 Remove supernatant and discard
- 4 Add  2 μL sample in triplicate
- 5 Add  Acetonitrile **Contributed by users** to a concentration of 70% and mix once by pipette to create a uniform suspension



6 Allow to rest ⌚ 00:10:00

10m



Wash - DO NOT REMOVE FROM MAGNET

10s

7 Place tube on magnet and allow Microparticles to clear ⌚ 00:00:10

10s



7.1 Remove supernatant and discard

8 Add 1 mL 100% ACN **Thermo Scientific** and incubate for ⌚ 00:00:10

10s

8.1 Remove supernatant

9 Add 1 mL 70% Ethanol **Contributed by users** and incubate for ⌚ 00:00:10

10s

9.1 Remove supernatant

Digestion

10s

10 Remove tube from magnet and add 8.5 Tris **P212121** and mix thoroughly

11 Add 5 µg Trypsin **P212121 Catalog #RP-T70010** and incubate for ⌚ 04:00:00

4h

12 Quench to a final concentration [M] 1 % volume TFA **Thermo Scientific**

13 Mix and place samples on magnet for ⌚ 00:01:00 and transfer supernatant to tubes

1m

14 Wash beads with 50-100 µL [M] 1 % volume TFA **Thermo Scientific** for ⌚ 00:02:00 with continuous mixing


2m



14.1 Remove supernatant and pool with eluate from Step 13

15 Centrifuge tubes  20000 x g, 00:10:00

10m

16 Transfer supernatant to new tube and freeze down at  -80 °C