



Oct 18, 2021

Background control of Gotcha RCA

Chia-Hsien Shih¹¹Chung Shan Medical University

1


dx.doi.org/10.17504/protocols.io.by5tpy6n

CSMU_Taiwan



Chia-Hsien Shih
Chung Shan Medical University

This protocol aims to verify that GotCha design works as expected.

DOI

dx.doi.org/10.17504/protocols.io.by5tpy6n

Chia-Hsien Shih 2021. Background control of Gotcha RCA. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.by5tpy6n>



_____ protocol ,

Oct 17, 2021









Oct 18, 2021

54163



Preparation

- 1 Add **5 μ L** of functional beads(Gotcha) into eppendorf
- 2 Centrifuge for **15000 rpm, 00:05:00** and remove supernatant. Make sure that eppendorf^{5m} should put on DynaMag when removing supernatant.

Protocol of Control group without miRNA


- 3 Add  **3 µL** of 10X phi29 polymerase reaction buffer into eppendorf with functional beads
- 4 Add  **23.4 µL** RNase-free water
- 5 Add  **3 µL** of 2mM dNTPs
- 6 Add  **0.6 µL** of 10U/ul phi29 polymerase
- 7 Pipetting to mix well
- 8 Incubate for  **02:00:00** at  **Room temperature** 2h
- 9 Add  **1.5 µL** of 20X evagreen dye
- 10 Add  **1.5 µL** of 0.5M EDTA buffer to suspend the reaction

Protocol of Control group without phi29 polymerase

- 11 Add  **3 µL** of 10X phi29 polymerase reaction buffer into eppendorf with functional beads
- 12 Add  **21 µL** RNase-free water

- 13 Add  **3 μ L** of 2mM dNTPs
- 14 Add  **3 μ L** of 100nM miRNA
- 15 Pipetting to mix well
- 16 Add  **1.5 μ L** of 20X evagreen dye
- 17 Add  **1.5 μ L** of 0.5M EDTA buffer to suspend the reaction

Measuring

- 18 Load  **20 μ L** of reaction solution into 384-well plate
- 19 Measure the fluorescence excitation and emission intensity