



Aug 22, 2020

Isolation of Shigella pathogens from oysters

Sade Aisha Folashade John¹, Patrick E. Akpaka¹, Chandrashekhar Unakal¹, Arvind Kurhade¹, Angel Justiz-Vaillant¹

¹University of the West Indies St. Augustine

1 Works for me dx.doi.org/10.17504/protocols.io.bj6pkrdn

University of the West Indies angel.vaillant@sta.uwi.edu

Angel Justiz-Vaillant
University of the West Indies St. Augustine

DOI

dx.doi.org/10.17504/protocols.io.bj6pkrdn

PROTOCOL CITATION

Sade Aisha Folashade John, Patrick E. Akpaka, Chandrashekhar Unakal, Arvind Kurhade, Angel Justiz-Vaillant 2020. Isolation of Shigella pathogens from oysters. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bj6pkrdn>

LICENSE

———— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Aug 22, 2020

LAST MODIFIED

Aug 22, 2020

PROTOCOL INTEGER ID

40879

DISCLAIMER:

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to [protocols.io](#) is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with [protocols.io](#), can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

- 1 The filtered homogenate was streaked using standardized loops on Salmonella-Shigella (SS) agar, the loops were flamed periodically to ensure sterility. This was done in duplicate. The plates were then incubated at 34°C. After an overnight incubation at 34°C, the plate with the best significant and adequate colonies was used to test for the presence of Shigella. Those colonies that were morphologically characteristic for Shigella were gram stained:
- 2 Each slide was labelled according to location and given a number. A circle using a wax pencil was drawn on the underside of each slide.

- 3 Drop of sterile saline solution was placed into the circle and using a sterile loop, a single colony from the plate was taken and mixed into it.
- 4 The smear was allowed to air dry and then heat fixed by holding the slide on one side with a sterile forcep and passing the entire slide through the flame of a Bunsen burner two or three times with the smear side up. Making sure that overheating does not occur.
- 5 The heat fixed slide was placed on a straining tray.
- 6 The slide was gently flood with crystal violet and left to stand for 1 minute.
- 7 Using a wash bottle of distilled water, the slide was tipped slightly and gently rinsed.
- 8 Then using Gram's iodine, the slide was gently flood and let stand for 1 minute.
- 9 The slide was tipped slightly and gently rinsed with distilled was using a wash bottle.
- 10 95% ethyl alcohol was used to decolorize the smear by rinsing the slide for 5-10 seconds.
- 11 The slide was then immediately rinsed with water.
- 12 Using Safranin, the slide was gently flood to counterstain and let stand for 45 seconds.
- 13 The slide was tipped slightly and gently rinsed with distilled was using a wash bottle.
- 14 The slide was then blow dried with bibulous paper.
- 15 Using a light microscope under oil immersion, the smear was observed to the determine whether the colonies were gram negative or gram positive.

- 16 This was done per presumptive *Shigella* colony observed on each plate.
- 17 A single colony of presumptive *Shigella* was then subcultured on SS agar and incubated at 34°C prior to biochemical testing. These pure single colonies were further analysed using biochemical tests (Indole, Citrate, Urease, Motility and Triple Sugar Iron Agar). The tests were done together with a control *Shigella* colony taken from the laboratory in order to make a comparison. A colony of the organism was then saved in a cryovial containing Nutrient Agar -70 °C. This was done per *Shigella* spp. found.