



Mar 17, 2022

Multiplexed snRNA-seq from frozen human brain samples

Marcos Nascimento¹

¹[Department of Neurological Surgery/ Eli and Edythe Broad Center of Regeneration Medicine and St em Cell Research, University of California San Francisco]



dx.doi.org/10.17504/protocols.io.kqdg369keg25/v1

AAB Lab



DOI

dx.doi.org/10.17504/protocols.io.kqdg369keg25/v1

Marcos Nascimento 2022. Multiplexed snRNA-seq from frozen human brain samples. **protocols.io**

https://dx.doi.org/10.17504/protocols.io.kqdg369keg25/v1

____ protocol,

Nov 16, 2020

Mar 17, 2022

44568

Sectioned or finely chopped frozen human tissue (10-30 mg).

Solutions (Detailed recipes here)

- Lysis Buffer (*LB*, 3ml/sample)
- Wash and Resuspension Buffer (WRB, 7ml/sample)
- Sucrose Buffer (10ml/sample)

2000U of RNAse inhibitor/sample

Preparation 10m

1 Check if the Rotor SW32Ti rotor is in the cold room

motocols.io

1

Citation: Marcos Nascimento Multiplexed snRNA-seq from frozen human brain samples https://dx.doi.org/10.17504/protocols.io.kqdq369keq25/v1

- 3 Add RNAse inhibitor (Sigma cat. # 3335402001) to LB and WRB (0.2U/ul).
- 4 Put lysis buffer and sucrose solutions on ice.

Nuclei Isolation

4h

- Use glass dounce homogenizer (Thomas Scientific; Catalog # 3431D76; size A). Put douncer on ice, pipette 1mL of lysis in the douncer. Transfer tissue either using spatula or P1000 pipette with cut tip and additional lysis buffer. Bring the total volume of lysis buffer in the douncer to 3mL. § On ice
- 6

Dounce tissue on ice with 10 strokes or until no chunks of tissue are visible. § On ice

- 7 Transfer homogenized tissue in lysis buffer into a labeled thick wall ultracentrifuge tube on ice (Beckman Coulter; 355631). § On ice
- 8 Carefully pipette 9 mL of Sucrose solution to the bottom of the tube containing Lysis buffer. Be careful not to introduce bubbles. You should see two clearly separated phases: sucrose on the bottom and cloudy homogenate on top. § On ice
- 9 When you are done with all samples weigh them and bring to the same weight by adding Lysis buffer.

10



2h 30m

Load the samples to SW32Ti rotor (needs to be swing bucket). If using less than 6 samples still balance with empty buckets.

3107163 rcf, 4°C, 02:30:00 , 29500 RPM on SW32Ti Rotor

11 /

20m

After the spin, transfer samples on ice and carefully remove the supernatant using a P200 tip cut at an angle and vacuum. Make sure not to touch the bottom (stick to the wall and tilt the tube), but remove all the liquid. Carefully pipette 200uL of WRB on the bottom. Wait 20 min on ice.

- 12 Meanwhile, transfer materials to the tissue culture room. Prepare eppendorf tubes with 10ul of DAPI for each sample.
- Add 800ul of WRB (for a total of 1ml of WRB) and resuspend cells. § On ice
- 14 Filter twice using Miltenyi Pre-separation filters (30um). (130-041-407)
- 15 Add 10ul of each sample to 10ul of DAPI. Count nuclei in each sample using a hemocytometer.

You should have at least 10⁵ nuclei/sample



1mg of human cortex typically yields ~10⁴ cells

CellPlex Barcoding

1h 30m

16

Thaw Cell Multiplexing Oligo at room temperature. Vortex for 5 sec and centrifuge for 5 sec.

17 Centrifuge nuclei using a swing bucket rotor **\$500 rcf, 4°C, 00:10:00**

10m

27

being too weak to demultiplex samples by itself, you can use it to match donors identified on

freemuxlet with each sample.