



Aug 23, 2020

© Copy of Copy of ELISA for quantification of human C6 in serum or plasma.

Angel A Justiz-Vaillant¹, Belkis Ferrer-Cosme²

¹University of the West Indies St. Augustine; ²"Saturnino Lora Torres' Provincial Teaching Clinical Surgical Hospital. Cuba

In Development dx.doi.org/10.17504/protocols.io.bj77krrn

Angel Justiz-Vaillant

DOI

dx.doi.org/10.17504/protocols.io.bj77krrn

University of the West Indies St. Augustine

PROTOCOL CITATION

 $\label{lem:composition} Angel\,A\,Justiz-Vaillant,\,Belkis\,Ferrer-Cosme\,2020.\,Copy\,of\,Copy\,of\,ELISA\,for\,quantification\,of\,human\,C6\,in\,serum\,or\,plasma..\,\textbf{protocols.io}$

https://dx.doi.org/10.17504/protocols.io.bj77krrn

 ι.	ь.	IN	⊢

This is an open access protocol distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Aug 23, 2020

LAST MODIFIED

Aug 23, 2020

PROTOCOL INTEGER ID

40927

- 1 An anti-human C6 coating antibody is adsorbed onto the microwells by incubation overnight at 4°C with carbonatebicarbonate buffer.
- 2 Add 50 µl of human serum or plasma. Human C6 present in the serum or plasma binds to antibodies adsorbed into the microwells.
- 3 The microplate is blocked with 3% non-fat milk-PBS buffer and later wash to remove unbound proteins.
- 4 Fifty (50) μl of biotin-conjugated anti-C6 antibody is added. The optimal dilution must be investigated.
- 5 The microplate is rewashed with PBS-Tween 20 buffer, pH 7.4.

Citation: Angel A Justiz-Vaillant, Belkis Ferrer-Cosme (08/23/2020). Copy of Copy of ELISA for quantification of human C6 in serum or plasma.. https://dx.doi.org/10.17504/protocols.io.bj77krrn

6	One hundred μ I of streptavidin-HRP conjugate is added and it binds to the biotin-conjugated anti-C6 antibody. The optimal dilution of this conjugate must be investigated.
7	The plate is washed following incubation to remove the unbound Streptavidin-HRP.
8	Add 100 μl of 3,3',5,5'- tetramethylbenzidine (TMB; Sigma-Aldrich) into each well.
9	Incubate the microwells in the dark for 15 min.
10	A colored product is formed in proportion to the quantity of C6 present in the sample or standard.
11	The reaction is terminated by addition of 100 $\mu l3MH2SO4$ and the absorbance is measured at 450 nm.
12	A standard curve is made from 7 human C6 standard dilutions and the human C6 sample concentration is determined.