



Mar 01, 2022

Rotarod

Haley Geertsma¹¹University of Ottawa

1

dx.doi.org/10.17504/protocols.io.b5srq6d6**Haley Geertsma**
University of Ottawa

This protocol is used to test mice in the Rotarod behavioural test.

DOI

dx.doi.org/10.17504/protocols.io.b5srq6d6

Haley Geertsma 2022. Rotarod. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.b5srq6d6>



_____ protocol ,

Mar 01, 2022

Mar 01, 2022

58929

IITC Life Sciences Rotarod

- 1 Habituate mice in the testing room for 60 minutes. 1h
- 2 Place each mouse on the rod and turn on the rotarod machine at 4-40rpm over 5 minutes. 5m
Record how long each mouse can stay on the pole.
- 3 Once testing is complete, place mice back into their home cage and wait for 10 minutes. 10m

4 Repeat this for 4 trials per day then repeat for 3 days (12 trials total).

3d