



OCT 09, 2023

OPEN ACCESS



DOI:
dx.doi.org/10.17504/protocols.io.261gedwxwv47/v1

Protocol Citation: NUS iGEM 2023. GFP-sacB Characterization.
protocols.io
<https://dx.doi.org/10.17504/protocols.io.261gedwxwv47/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
 We use this protocol and it's working

Created: Oct 08, 2023

Last Modified: Oct 09, 2023

GFP-sacB Characterization

NUS iGEM¹

¹National University of Singapore



夏雨 陈

ABSTRACT

2023 NUS-Singapore iGEM team followed this protocol to characterise their New Composite Part "GFP-sacB" with IPTG and sucrose solution at various concentrations.

MATERIALS

- LB media
- M9 Media
- Correct Antibiotics
- IPTG Solution
- Sucrose Solution
- DI Water




SAFETY WARNINGS






- Proper lab PPE must be worn at all times.
- Since cells are used in this protocol, a Biosafety Cabinet (BSC) is required to ensure safety.

Keywords: sacB, GFP,
Negative Selection, Negative
Selection Marker,
Characterization

Cell Inoculation and Incubation (Day Before Characterization)

- 1 Inoculate cells with GFP-sacB gene from the cell stock.
- 2 Add  5 mL of LB media and  5 μ L of the appropriate antibiotic to a Falcon tube.
- 3 Incubate the Falcon tube at  37 °C in an incubator.


Sample Preparation (96-well Plate)

- 4 Decide the concentrations of IPTG and sucrose solutions for characterization.
- 5 Prepare new Falcon tube(s).
- 6 Add  5 mL of M9 media,  5 μ L of the appropriate antibiotic, and  100 μ L of cells cultured the previous day.
- 7 Add the required volume of IPTG to reach the desired concentration.

8 Incubate the cells for 02:00:00 at 37 °C .

9 Prepare a sterile 96-well plate.

10 Design an appropriate plate map, each sample (with a particular IPTG and sucrose concentration) must be repeated 3 times. Example of plate map:



One of the actual plate map used by the NUS-Singapore iGEM 2023 team when characterizing the GFP-sacB gene.

11 Adjust sucrose concentration by adding DI water and sucrose solution to each well (final volume of 20 µL) according to the plate map.


12 Add 200 µL of cultured cells to each well already containing DI water and sucrose solution (final well volume of 220 µL).

- 13 Add at least 3 wells of  220 µL of M9 media as the blank (negative control).

Characterization and Plate Reader Reading

- 14 Place the 96-well plate into the plate reader.
- 15 Create a protocol in the plate reader's software according to the following setting:

A	B
Plate Type	96 WELL PLATE
Set Temperature	Setpoint 37°C
	Preheat before moving to next step
Start Kinetic	Runtime 6:10:00 (HH:MM:SS), Interval 0:30:00, 13 Reads
Shake	Orbital: Continuous
	Frequency: 282 cpm (3 mm)
Read 1	Absorbance Endpoint
	Full Plate
	Wavelengths: 600
	Read Speed: Normal, Delay: 100 msec, Measurements/Data Point: 8
Read 2	Fluorescence Endpoint
	Full Plate
	Excitation: 485, Emission: 528
	Optics: Top, Gain: 100
	Light Source: Xenon Flash, Lamp Energy: High
	Read Speed: Normal, Delay: 100 msec, Measurements/Data Point: 10
	Read Height: 7 mm
End Kinetic	

- 
- 16** Start the continuous reading.
 - 17** After the reading is complete, remove the 96-well plate from the plate reader.
 - 18** Save the data as a CSV file for future analysis.
 - 19** Discard the used 96-well plate in a biohazard bin.