



Jun 21, 2020

Cleaning Microscopes (Covid-19 Shutdown Re-Entry, Srivastava Lab)

In 1 collection

Srivastava Lab¹¹Harvard UniversityOther dx.doi.org/10.17504/protocols.io.bhngj5bw

Re-Entering Labs Post COVID-19 Shutdown

Mansi Srivastava

ABSTRACT

This is based on: [Recommendations for operating Core Facilities in a research environment during the SARS-CoV-2 pandemic](#):

Disclaimer: The SOP presented here has been designed by the Srivastava Lab at Harvard University and is being shared outside of Harvard for informational purposes only. If you choose to reuse or repurpose this SOP in another location, please note that you do so at your own risk; you should ensure that any local guidance is also adhered to. None of the authors, contributors, administrators, or anyone else associated with protocols.io, can be held responsible for your use of the information contained in or linked to from these web pages.

DOI

dx.doi.org/10.17504/protocols.io.bhngj5bw

PROTOCOL CITATION

Srivastava Lab 2020. Cleaning Microscopes (Covid-19 Shutdown Re-Entry, Srivastava Lab). protocols.io
dx.doi.org/10.17504/protocols.io.bhngj5bw

COLLECTIONS ⓘ



Collection of SOPs for Covid-19 Shutdown Re-Entry, Srivastava Lab

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Jun 18, 2020

LAST MODIFIED

Jun 21, 2020

OWNERSHIP HISTORY

Jun 18, 2020 Lenny Teytelman protocols.io

Jun 19, 2020 Mansi Srivastava

PROTOCOL INTEGER ID

38312

PARENT PROTOCOLS

GUIDELINES

Disclaimer: The SOP presented here has been designed by the Srivastava Lab at Harvard University and is being shared outside of Harvard for informational purposes only. If you choose to reuse or repurpose this SOP in another location, please note that you do so at your own risk; you should ensure that any local guidance is also adhered to. None of the authors, contributors, administrators, or anyone else associated with protocols.io, can be held responsible for your use of the information contained in or linked to from these web pages.

- 1 When using the microscope and adjacent computer, **always wear gloves**
- 2 Wipe hard surfaces (knobs, body of scope, screens) with 70% ethanol before use
- 3 Place a new piece of plastic wrap over oculars before use



cleaning oculars too frequently can cause the rubber to degrade

- 4 When done, dispose of plastic wrap and wipe hard surfaces with 70% ethanol