



DEC 27, 2022

🌐 Protocol 2: Culturing *Spizellomyces punctatus* (Sp) prior to transformation day

📁 In 1 collection

Sarah M Prostack¹, Edgar M Medina¹, Erik Kalinka¹, Lillian Fritz-Laylin¹

¹University of Massachusetts at Amherst



fritzlaylinlab.umass

ABSTRACT

In order to increase transformation efficiency, a more synchronous culture of Sp is needed. This process involves subculturing active Sp plates in the days leading up to transformation day. At least two generations of Sp should be grown on antibiotic free K1 medium to ensure no trace antibiotics are present in the culture for transformation day. Proper sterile technique should be followed at all times when maintaining cultures of Sp; use either a laminar flow hood or work in the sterile area around an open flame.

ATTACHMENTS

[Spizellomyces transformation steps.pdf](#)

OPEN ACCESS

DOI:
[dx.doi.org/10.17504/protocols.io.36wgqjinxvk5/v1](https://doi.org/10.17504/protocols.io.36wgqjinxvk5/v1)

Protocol Citation: Sarah M Prostack, Edgar M Medina, Erik Kalinka, Lillian Fritz-Laylin 2022. Protocol 2: Culturing *Spizellomyces punctatus* (Sp) prior to transformation day. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.36wgqjinxvk5/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

Created: Nov 02, 2022

Last Modified: Dec 27, 2022

PROTOCOL integer ID:
72190

Keywords: *Spizellomyces punctatus* (Sp), transformation

MATERIALS TEXT

Materials

- Active, healthy culture of wild-type *Spizellomyces punctatus* ( *Spizellomyces punctatus* (Koch) Barr **ATCC Catalog #48900**)
- K1 agar plates (1.5% w/v) without any antimicrobial (see recipe)
- DS solution (see recipe)
- 50 mL conicals, sterile (such as  Centrifuge Tubes-Bag 50mL Centrifuge Tube - Bag Sterile **Cell Treat Scientific Products Catalog #229421**)
- 100-1,000 µL micropipette such as  Eppendorf Research Plus Single Channel pipette 100-1000 uL blue operating button f use with 100 **pipette.com Catalog #3123000063 ES-1000**
- Filter tips for the micropipette, sterile such as  TIPONE® FILTER TIPS **USA Scientific Catalog #1122-1830**
- 20 mL luer-lock syringes (such as  BD General Use Syringes **Fisher Scientific Catalog #22-124-967**)
- 25 mm syringe filter holders ( Cole-Parmer Advantec 43303010 Polypropylene Filter Holder for 25 mm Membranes **Fisher Scientific Catalog #NC9972954**), preloaded with Grade 1 Whatman paper ( Cytiva Whatman™ Qualitative Filter Paper: Grade 1 Circles **Fisher Scientific Catalog #09-927-223**), sterile
- Laminar flow hood, or open flame to maintain sterility
- 70% (v/v) ethanol for maintaining sterility (if using laminar flow hood)
- Incubator at  28 °C


BEFORE START INSTRUCTIONS

In order to prepare enough plates to provide an adequate amount of zoospores for transformation, subculture one plate of Sp per plasmid to be transformed. One active plate is enough to subculture 2-3 new plates. Complete this protocol roughly 36 hours and then again roughly 18 hours prior to the intended transformation time.

ATTACHMENTS

[Spizellomyces transformation steps.pdf](#)

Steps

- 1 Flood enough active Sp plates each with  1 mL of DS to fit your needs.

Note

One active plate is enough to seed 2-3 new plates and you will need one plate of active Sp per plasmid to be transformed on transformation day.

- 2 Holding the plate at an angle, gently wash the DS over the agar.





- 2.1 Wash the  1 mL of DS over the agar. (1/3)

- 2.2 wash the  1 mL of DS over the agar. (2/3)

- 2.3 Wash the  1 mL of DS over the agar. (3/3)


- 3 Collect all zoospores into a 50 mL conical.

- 4 Filter all zoospores into a new 50 mL conical using the sterile 25 mm syringe filter preloaded with Whatman Grade 1 filter paper.

- 5 Inoculate the proper amount of new K1 plates for your needs with  250 µL to  500 µL of filtered zoospores.

Note

The amount of spores transferred to new plates depends on the density of the plate being harvested, be careful not to overcrowd the plates, as zoospore release will decrease.

6 Add  500 μL of fresh DS to the newly inoculated plates.




7 Gently tilt the plates to spread the liquid across the entire surface of the agar.

8 Incubate at  30 $^{\circ}\text{C}$ in a humidity chamber.



Note

Using a humidity chamber is standard practice, as it prevents the agar from drying out. Placing the plates into a plastic bin with a beaker of ~  20 mL of water should be enough.

Complete this protocol in its entirety 36 and 18 hours prior to intended transformation time.