

Oct 01, 2024

Library cloning

 In 1 collection

DOI

dx.doi.org/10.17504/protocols.io.5qpvkrrzl4o/v1

Raining Wang¹, Melinda Wheelock¹

¹University of Washington

IGVF



Sarah Heidl

University of Washington, Starita Lab

OPEN  ACCESS



DOI: **dx.doi.org/10.17504/protocols.io.5qpvkrrzl4o/v1**

Protocol Citation: Raining Wang, Melinda Wheelock 2024. Library cloning. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.5qpvkrrzl4o/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

Created: May 09, 2024

Last Modified: October 01, 2024

Protocol Integer ID: 99788











Abstract

This is the protocol for inserting the library gene of interest into the pre-barcoded vector created in the previous protocol (Barcoded Vector Cloning).







Materials

Materials

- Barcoded library plasmid
- Negative control from Golden Gate reaction
-  SapI - 1,250 units **New England Biolabs Catalog #R0569L**
- rCutsmart Buffer (supplied with SapI restriction endonuclease)
- rSap
-  Molecular Biology Grade Water **Thomas Scientific Catalog #C001X09**
-  Rigid Strip 8-Strip Tubes **Thomas Scientific Catalog #1148A28**
- T4 DNA ligase buffer (supplied with T4 DNA ligase)
-  T4 ligase **New England Biolabs Catalog #M0202M**
-  Zymo DNA Clean & Concentrator Kit **Zymo Research Catalog #D4014**
-  Qubit™ 1X dsDNA BR Assay Kit **Thermo Fisher Scientific Catalog #Q33231**
- Qubit dsDNA 1X HS Assay Kit (Thermo Fisher Scientific, Catalog # Q33266)
-  Qubit™ Assay Tubes **Invitrogen - Thermo Fisher Catalog #Q32856**
- Ice
-  NEB 10-beta Electrocompetent E.coli - 6x0.1 ml **New England Biolabs Catalog #C3020K**
- SOC (supplied with NEB cells)
- pUC19 positive control (supplied with NEB cells)
- Electroporation cuvettes, 1mm (Thomas Scientific, Catalog # 1181C48)
- Glass culture tubes (15ml)
- 500-mL Erlenmeyer flasks
- LB liquid media
-  Carbenicillin **Gold Biotechnology Catalog # C-103-5**
- LB + Ampicillin agarose plates
-  1.7 ml microcentrifuge tube **Thomas Scientific Catalog #1159M35**
- Glass beads

Equipment

-  C1000 Touch™ Thermal Cycler **Bio-Rad Laboratories Catalog #1851148**
-  NanoDrop® 2000 Spectrophotometer **Thermo Fisher Catalog #ND-2000**
-  Qubit Fluorometer **Invitrogen - Thermo Fisher Catalog #Q32866**
-  MicroPulser Electroporator **Bio-Rad Laboratories Catalog #1652100**
- Shaking incubator
- Standing incubator

50- μ L reaction:

A	B	C	D	E
Reagent	Vol for 1 reaction (uL)	Volume for 5 reactions (uL)		
H2O			Error	1.05
DNA (1ug)			# reactions	5
rCutsmart Buffer	5	26.25		
Sapl	1	5.25		
Total volume:	50			

Ligation reaction:

A	B	C	D	E	F	G	H
			Total samples #	4	-ctl	Err	1.05
Component		Final conc.	1Rxn (uL)	Total			
H2O			7.027266667	29.51452	14.5		
T4 DNA ligase buffer		1X	2.5	10.5	2.5		
T4 DNA ligase (2,000U/ul)			0.5	2.1	0.5		
Plasmid (75ng/ul)	5400bp		1	4.2	1		
Sapl (10U/ul)		15	1.5	6.3	1.5		
Sub-total			11.02726667	46.31452			
Insert F1 (75ng/ul)	plasmid 4182bp/1656bp	174ng/ul (1:3mol ratio)	2.32	9.744	0		
Insert F2 (75ng/ul)	plasmid 3144/621bp	131.7ng/ul (1:3mol ratio)	1.756	7.3752	0		
Insert F3 (75ng/ul)	Plasmid 4781/2258bp	199ng/ul (1:3mol ratio)	2.653333333	11.144	0		
Insert 4 (50ng/ul)	Library 892bp	37.17ng/ul (1:3mol ratio)	0.7434	3.12228	0		
Total			20	#REF!			

Additional digestion reaction:

A	B	C	D	E	F
		Total samples #	5	Err	1.05
Component	Final conc.	1Rxn (uL)	Total		
H2O		2.5	13.125		


	A	B	C	D	E	F
	rCutSmart Buffer		1	5.25		
	Sapl (10U/ul)	10U	1.5	7.875		
	Sub-total		#REF!	#REF!		
	Total		5	#REF!		

Before start

- Ensure there are enough maxi-prep kits available to use.



Predigest the barcoded backbone

- 1 Set up a  50 μ L reaction as outlined below:

A	B	C	D	E
Reagent	Vol for 1 reaction (uL)	Volume for 5 reactions (uL)		
H2O			Error	1.05
DNA (1ug)			# reactions	5
rCutsmart Buffer	5	26.25		
SapI	1	5.25		
Total volume:	50			


- 2 Set the thermal cycler for the program below:

A	B	C
Step	Temperature	Time (hh:mm:ss)
1	37°C	03:00:00
2	65°C	00:20:00
3	10°C	hold

Note

- Spike in additional 1ul of SapI enzyme every 60 mins, to each sample
- At the 2:30 mark, add 1ul of rSap- this is to dephosphorylate the 5' end to re-ligate the ends.

- 3 Run all digested samples on a 0.5% TBE gel to ensure complete digestion of the backbone.

- Also run  1 μ g of original starting template as a visual reference for the supercoiled plasmid.

- 4 Gel extract the linearized product.

5 Use the Clean & Concentrator kit on the gel extraction product.

- To avoid issues with freeze-thaws, aliquot the cleaned, linearized product in small volumes.

Golden gate assembly of variant library inserts

4h

6

Note

This adds the variant library inserts to the pre-barcoded vector using SapI.

Set up reaction as outlined below:

6.1

Note

Here we normally use 1:3 nmol backbone to each insert ratio, below is an example with a backbone and 4 insert frag. please note three of the fragments are on the helper plasmid, thus the calculation is done on the whole helper plasmid.

6.2 Here is a [ligation calculator](#) by NEB to use if you don't have a formula set up.

6.3 Here is [NEB's general protocol for GG with SapI](#).

6.4 Neg Ctl: Using the vector with no insert.

A	B	C	D	E	F	G	H
			Total samples #	4	-ctl	Err	1.05
Component		Final conc.	1Rxn (uL)	Total			
H2O			7.027266667	29.51452	14.5		
T4 DNA ligase buffer		1X	2.5	10.5	2.5		
T4 DNA ligase (2,000U/ul)			0.5	2.1	0.5		



A	B	C	D	E	F	G	H
Plasmid (75ng/ul)	5400bp		1	4.2	1		
SapI (10U/ul)		15	1.5	6.3	1.5		
Sub-total			11.02726667	46.31452			
Insert F1 (75ng/ul)	plasmid 4182bp/1656bp	174ng/ul (1:3 mol ratio)	2.32	9.744	0		
Insert F2 (75ng/ul)	plasmid 3144/621bp	131.7ng/ul (1:3 mol ratio)	1.756	7.3752	0		
Insert F3 (75ng/ul)	Plasmid 4781/2258bp	199ng/ul (1:3 mol ratio)	2.653333333	11.144	0		
Insert 4 (50ng/ul)	Library 892bp	37.17ng/ul (1:3 mol ratio)	0.7434	3.12228	0		
Total			20	#REF!			

7 Set a thermocycler for the following program:

- Spike in 1ul of SapI every hour.

A	B	C
Step	Temperature	Time (hh:mm:ss)
1	37°C	00:05:00
2	16°C	00:05:00
3	Go to step 1 for 15 total cycles	
4	60°C	00:15:00
5	4°C	hold

8 Set up additional digestion reaction as follows:

A	B	C	D	E	F
		Total samples #	5	Err	1.05
Component	Final conc.	1Rxn (uL)	Total		
H2O		2.5	13.125		
rCutSmart Buffer		1	5.25		
SapI (10U/ul)	10U	1.5	7.875		
Sub-total		#REF!	#REF!		
Total		5	#REF!		



9 Incubate at 37 °C for at least 02:00:00 , or Overnight .

4h



- Spike in SapI every hour.

10 Purify the product with the Zymo Clean & Concentrator kit.

- Pool 2 sample reactions (40 µL) for 4 columns.

11 Measure the final eluate with Nanodrop to look for a normal purity curve.

- 260/280 ratio should be:

12 Quantify the product concentration with Qubit BR assay.

- The final eluate should be >5ng/uL for the samples and very low for the negative control.

Transformation with electro-competent cells

17h 56m

13 Prepare cells for transformation

13.1 Warm SOC to 37 °C .

13.2 Chill empty cuvettes On ice for 00:15:00 .

15m

13.3 Thaw NEB electrocompetent cells (100 uL/tube) on ice for approximately 00:05:00 .



5m

13.4 To chilled cuvettes, add 26 µL of thawed cells.






- Pipette to first stop only - do not create air bubbles as this will cause the electroporation to fail.
- Can use long gel-loading tips to do this more easily.



13.5 Add  1 μL -  2 μL of cleaned DNA from golden gate reaction and controls to the corresponding cuvette containing cells.

- Pipette to first stop only - do not create air bubbles.
- Do not add more than 3 μL of plasmid.




13.6 Flick the cuvette gently, then lightly tap on the bench. Place the cuvette back  On ice .

13.7 Incubate  On ice for  00:30:00 .

30m





14 Prepare for recovery


14.1 While the cells are incubating  On ice , prepare glass recovery tubes with  900 μL SOC each, and warm to  37 $^{\circ}\text{C}$.

14.2 Set the MicroPulser Electroporator to Manual for 2.0 ms.

- The electroporator is automatically set for 6 kV.

14.3 Prepare six overnight growth Erlenmeyer flasks with  300 mL of LB +  30 μL carbenicillin for each sample.

- Use 500 mL flasks for enough air circulation.

14.4 Label the LB + ampicillin dilution plates for each sample. Place in  37 $^{\circ}\text{C}$ after labeling to pre-warm.

- For Negative Golden Gate control, need 2 dilution plates
- For library and pUC19, need 3 dilution plates each.

15 Electroporation

15.1 Starting with the negative control cuvette, wipe the outside of the cuvette dry to remove condensation and ice.

15.2 Place the cuvette in the electroporator, ensuring the metal on the cuvette makes contact with the electrodes, and press the Pulse button.

- If the electroporator makes a loud noise and displays “Arc” on the digital readout, the sample has failed and is no longer good.


15.3 Remove the cuvette from the electroporator and place back  On ice .

15.4 Repeat with all cuvettes prepared.



16 Sample recovery

16.1 Remove the glass recovery tubes containing the SOC from  37 °C .

16.2 Transfer the cells from each cuvette to the glass recovery tubes.






16.3 From the warmed bottle of SOC, add  100 µL to wash the electroporation cuvette, then transfer from the cuvette to the corresponding glass recovery tube. Repeat this for all samples.

- Long gel-loading tips may be helpful for sample recovery from the cuvettes.

16.4 Place the glass recovery tubes into the shaking incubator for  01:00:00 at  37 °C .

1h

17 Prepare dilution tubes and plates

17.1 As the cells are recovering, prepare microcentrifuge tubes for serial dilution of the library and controls. We recommend to start with a bottleneck scheme by dividing the 1ml recovery culture to 6 portions:  500 µL ,  250 µL ,  125 µL , 2x  50 µL , and  25 µL .

- For each portion, prepare 4 tubes (D0-D3) with the listed volume of LB. Label three plates (D1-D3) for each bottleneck portion (total of 18 plates).

A	B	C	D	E	F
	Starting (S)	D0	D1	D2	D3
Volume	300ml	300ul from S	100ul of D0 + 900ul of LB	100ul of D0 + 900ul of LB	100ul of D0 + 900ul of LB
Plating		No plating	100ul	100ul	100ul

A	B	C	D	E	F
Dilution		0.7361111 11	1:30,000	1:300,000	1:3,000,000

17.2 For Puc, prepare two tubes and plates as the following Negative control and H₂O.

A	B	C	D
	Starting (S)	D0	D1
Volume	1ml	20ul from S +180 ul of LB	100ul of D0 + 900ul of LB
Plating		100ul	100ul
Dilution		1: 100	1:1,000




17.3 For Negative control and H₂O, plate  100 mL from the  1 mL recovery culture to the respective plates.

17.4



3m

Note

Prepare bottleneck for overnight growth and plating dilution plates.

- After the 1 hour recovery, add the appropriate portion(as below) from the  1 mL culture to the respective 500-mL flask based on the bottlenecking scheme.
- Place all 500 flasks in the  37 °C shaking incubator for  00:03:00 to ensure it is well mixed. No need to prepare overnight growth flasks for Puc19, negative control, and H₂O controls.

A	B	C	D	E
	D0 (1:1000)	D1 (1:30,000)	D2 (1:300,000)	D3 (1:3,000,000)
500ul				
250 ul				
125ul				
50ul-1				
50ul -2				
25ul -1				

17.5 Mix each dilution well. Plate each dilution according to the dilution scheme above. Incubate  Overnight at  37 °C .

3m



17.6 Incubate

🕒 Overnight

 growth flasks in

🌡️ 30 °C

 for



🕒 14:00:00

 -

🕒 16:00:00

 .

1d 6h 3m



Estimating barcode coverage for bottlenecking

18 Count colonies on dilution plates to calculate the approximate variant to barcode coverage

For examples:
Library X : 205 x 22 variant changes = 4510 aa variants. We want to have to pick the flask that produces between 45,100 -135,300 colonies. This indicates between 10-30x coverage over the variants.

Note

Best to average all three dilution plates for each bottleneck portion to ensure the most accurate estimation proceed to the maxi preps if the coverage is good.

A	B	C	D	E
	D0 (1:1000)	D1 (1:10,000)	D2 (1:100,000)	D3 (1:1,000,000)
500ul				
250 ul				
125ul				
50ul				
50ul -1				
25ul -2				