

Oct 30, 2020

# BSCI:414--Lab 9: Cloning the SARS-CoV-2 Spike Protein into a E. coli Protein Expression Vector

Harley King<sup>1</sup>

<sup>1</sup>NIST Center for Neutron Research National Institute of Standards and Technology Gaithersburg, MD 20899 Department of Materials Science and Engineering University of Maryland College Park, MD 20742-2115

1 Works for me This protocol is published without a DOI.

USG Fall 2020 BSCI:414



Harley King

NIST Center for Neutron Research National Institute of Stand...

## ABSTRACT

The SARS-CoV-2 spike protein does not express solubly in E. coli cells. This may be due to E. coli's inability to glycosylate or create disulfide bonds. The spike protein can be cloned into a vector containing the periplasmic localization peptide (pelB) like Novagen's pET-26 so the protein is transported to the oxidizing environment of the periplasm. Disulfide bonds may also be favored by expressing in modified E. coli cells like NEB's SHuffle strain. Glycosylation may be improved by expressing in a modified E. coli cell line capable of N-linked glycosylation like CLM24.

## PROTOCOL CITATION

Harley King 2020. BSCI:414--Lab 9: Cloning the SARS-CoV-2 Spike Protein into a E. coli Protein Expression Vector. **protocols.io**  
<https://protocols.io/view/bsci-414-lab-9-cloning-the-sars-cov-2-spike-protei-bn9pmh5n>

## LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Oct 30, 2020

## LAST MODIFIED

Oct 30, 2020

## PROTOCOL INTEGER ID

44047

## ABSTRACT

The SARS-CoV-2 spike protein does not express solubly in E. coli cells. This may be due to E. coli's inability to glycosylate or create disulfide bonds. The spike protein can be cloned into a vector containing the periplasmic localization peptide (pelB) like Novagen's pET-26 so the protein is transported to the oxidizing environment of the periplasm. Disulfide bonds may also be favored by expressing in modified E. coli cells like NEB's SHuffle strain. Glycosylation may be improved by expressing in a modified E. coli cell line capable of N-linked glycosylation like CLM24.


## Explore the Spike Protein

- 1 Explore the structure of the SARS-CoV-2 spike protein using [Cov3D](#).




Explore the spike protein using this link: <https://cov3d.ibbr.umd.edu/viewer7k90>

- 2 Read abstract of [Zuo et al.](#) concerning previous attempts at expressing SARS-CoV spike protein fused with SUMO tag.



C-terminus of SARS-CoV spike protein expressed more solubly with SUMO tag.

### Discuss Gibson Cloning



Gibson cloning uses three enzymes and overlapping, complimentary bases to ligate fragments into vectors.

### Clone spike -sumo into pET28.

- 4 Using a Gibson cloning strategy in Benchling, clone SARS-CoV-2 spike-SUMO into pET28. Open three plasmids in Benchling. Cleave all other tabs.
1. SumoPro
  2. pET-28 Vector
  3. SARS-CoV2-Spike-delta21