



Oct 14, 2020

Plasmid Extraction

Hung Liang Pai¹, Huan Jui Chang¹¹Chung Shan Medical University**1** Works for me This protocol is published without a DOI.Huan Jui Chang
Chung Shan Medical University

ABSTRACT

This protocol is used for extracting plasmid from E. coli DH5-Alpha strains.

PROTOCOL CITATION

Hung Liang Pai, Huan Jui Chang 2020. Plasmid Extraction. **protocols.io**
<https://protocols.io/view/plasmid-extraction-bh5bj82n>

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Jul 01, 2020

LAST MODIFIED

Oct 14, 2020

PROTOCOL INTEGER ID

38787

GUIDELINES

The manual from manufacturer is provided here.

Presto™ Mini Plasmid Kit.pdf

This protocol is usually performed after overnight incubation of bacteria in a centrifuge tube.

MATERIALS

NAME	CATALOG #	VENDOR
Presto™ Mini Plasmid Kit (PDH100 PDH300)	PDH100	Geneaid

ABSTRACT

This protocol is used for extracting plasmid from E. coli DH5-Alpha strains.


BEFORE STARTING

Follow the "preparation" section of protocol.

Preparation

- 1 Set the dry bath incubator at **37 °C**
- 2 Prepare Presto™ Mini Plasmid Kit
- 3 Take out the centrifuge tube containing medium with bacteria that has been cultured overnight.


Protocol




4 Centrifuge the centrifuge tube under this condition -  **15000 x g, 00:06:00** .

5 Discard the supernatant





Recommendation: Keep the bacterial pellet away from supernatant if plausible to insure minimum loss of bacteria.

6 Add  **200 µl** of PD1 reagent to the centrifuge tube. Vortex the tube to completely mix bacterial pellet with PD1.

7 Add  **200 µl** of PD2 reagent to the centrifuge tube. Shake the tube gently and let it stand for  **00:02:00** to  **00:05:00** .






Warning: Do not vortex the mixture after adding PD2!!!!



8 Add  **300 µl** of PD3 reagent to the centrifuge tube. Shake the tube gently and centrifuge the tube under this condition -  **16000 x g, 00:10:00** .






Warning: Do not vortex the mixture after adding PD3!!!!
Tips: Prepare the PDH column when centrifuging.

9 After centrifugation, transfer the supernatant into a PDH column and discard the pellet.
Centrifuge the PDH column under  **15000 x g, 00:00:30** and discard the flowthrough liquid in the collection tube.

10 Add  **400 µl** of Wash 1 Buffer into the PDH column and centrifuge under  **15000 x g, 00:00:30** .
Discard the flowthrough in the collection tube.


11 Add  **600 µl** of Wash 2 Buffer into the PDH column and centrifuge under  **15000 x g, 00:00:30** .
Discard the flowthrough in the collection tube.

12 Centrifuge again under  **15000 x g, 00:03:00** and replace the collection tube with a new eppendorf. (Discard the flowthrough)

- 13 Add  **50 µl** of Elution Buffer into the middle of PDH column and incubate the PDH column in a drybath incubater under  **37 °C** .



Recommendation: Add the Elution Buffer directly onto the silicon membrane to insure optimal eluting efficacy.

- 14 Centrifuge under  **15000 x g, 00:02:00** . The flowthrough in the eppendorf contains plasmids derived from the cultured bacteria.



Recommendation: Proceed to measure the concentration of plasmid.