



AUG 09, 2023

## Mitophagy induction using Oligomycin/Antimycin A

Louise Uoselis<sup>1</sup>

<sup>1</sup>Lazarou Lab, WEHI



Louise Uoselis  
WEHI

### ABSTRACT

Mitophagy induction in HeLa cells using Oligomycin/Antimycin A.

OPEN  ACCESS



#### DOI:

[dx.doi.org/10.17504/protocols.io.14egn32yql5d/v1](https://dx.doi.org/10.17504/protocols.io.14egn32yql5d/v1)

**Protocol Citation:** Louise Uoselis 2023. Mitophagy induction using Oligomycin/Antimycin A. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.14egn32yql5d/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working  
We use this protocol and it's working

**Created:** Aug 09, 2023

**Last Modified:** Aug 09,  
2023

**PROTOCOL integer ID:**  
86203





**Keywords:** ASAPCRN

## Day 1

- 1 Seed cells, aiming for a confluency of 80-90% at the time of treatment the next day.

## Day 2

1h

- 2 Feed cells for  01:00:00 in an appropriate volume of standard growth media. 1h
- 3 To start the treatment, replace the media in each well with standard growth media that contains  10 micromolar ( $\mu\text{M}$ ) Oligomycin,  4 micromolar ( $\mu\text{M}$ ) Antimycin A, and  10 micromolar ( $\mu\text{M}$ ) QVD.
- 4 Harvest the samples after the desired treatment times.