



Mar 25,
2020

Bionano genome mapping from animal tissue

In 1 collection

Graham J Etherington¹, Darren Heavens¹, David Baker¹, Ashleigh Lister¹, Rose McNelly¹, Gonzalo Garcia¹, Bernardo Clavijo¹, Iain Macaulay¹, Wilfried Haerty¹, Federica Di Palma¹

¹The Earlham Institute, Norwich Research Park, Norwich, NR4 7UZ, United Kingdom

1 Works for me dx.doi.org/10.17504/protocols.io.bd7ei9je

Graham Etherington

ABSTRACT

Bionano genome map protocol

- 1 The Bionano genome maps were prepared by taking approximately **10 mg** of DNA (in this case extracted from a polecat) from the sample stored in 100% EtOH.
- 2 The IrysPrep Animal Tissue DNA Isolation from Fibrous Tissue protocol was followed using the IrysPrep Animal Tissue DNA kit (RE-013-10) from Bionano Genomics.
- 3 The animal tissue sample was cut up into <3mm pieces for homogenisation and then fixed in a 2% formaldehyde solution in kit-provided Homogenisation Buffer (HB) for **00:30:00** **On ice** before blending using a Qiagen Tissuereuptor.
- 4 Spin at 1500x g for **00:05:00** in centrifuge, remove supernatant and resuspend in around **50 µl** of HB buffer using wide bore tip to produce a total volume of **66 µl**.
- 5 **40 µl** of LMP agarose was then melted at **70 °C** and cooled to **43 °C** before addition to the cell resuspension and mixing using a wide bore tip.
- 6 One plug of around **90 µl** was cast using the Chef Mammalian Genomic DNA Plug Kit (Bio-Rad 170-3591).
- 7 Once cooled to **4 °C** the plug was added to a lysis solution containing **200 µl** proteinase K (QIAGEN 158920) and **2.5 ml** of Bionano lysis Buffer.

- 8 This was incubated at 50°C for 02:00:00 in a thermomixer, making a fresh proteinase K solution and incubating Overnight .
- 9 The 50 ml tubes were then removed from the thermomixer for 00:05:00 before 50 μl RNase A (Qiagen158924) was added and to the tubes, returned to the thermomixer for a further hour at 37°C .
- 10 The plugs were then washed 7 times in the Wash Buffer supplied with the Chef kit and 7 times in 1xTE.
- 11 The plug was removed and melted for 00:02:00 at 70°C followed by 00:05:00 at 43°C before adding 10 μl of 0.2U/ μl of GELase (Cambio Ltd G31200).
- 12 After 00:45:00 at 43°C the melted plug was dialysed on a 0.1 μM membrane (Millipore VCWP04700) sitting in 15 ml of 1xTE buffer in a small petri dish.
- 13 After 00:45:00 the sample was removed with a wide bore tip and mixed gently 00:45:00 and left Overnight at 4°C . A small amount was removed to QC on an Opgen Argus Q-Card and Qubit HS to calculate the DNA concentration.
- 14 300 ng of DNA was taken into the NLRS (Nick, Label, Repair and Stain) reaction using 1 μl Nt.BspQI (NEB R0644S). The optical maps were then generated using the Bionano Irys platform.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited