



Mar 09, 2021

# Thawing MRC5 Sub-Stock

Morrisey Lab<sup>1</sup><sup>1</sup>University of Pennsylvania

1

Works for me

This protocol is published without a DOI.

Ch

Chelseahortman

SUBMIT TO PLOS ONE

## ABSTRACT

### MRC5 Media

450ul DMFM/F12 (Cell Center cat#11320033)

45ul 10% FBS (Invitrogen cat#10437-028)

5ul 1% pen/strep (Cell Center cat#15240096)

## PROTOCOL CITATION

Morrisey Lab 2021. Thawing MRC5 Sub-Stock. [protocols.io](https://protocols.io/view/thawing-mrc5-sub-stock-bs5kng4w)  
<https://protocols.io/view/thawing-mrc5-sub-stock-bs5kng4w>



## LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Mar 09, 2021

## LAST MODIFIED

Mar 09, 2021

## PROTOCOL INTEGER ID

48012

- 1 Set aside 10mL of cold MRC5 media in 15mL conical tube.
- 2 Take 1 MRC5 sub stock aliquot out of liquid nitrogen.
- 3 Warm aliquot in
- 4 bath
- 5 Note\*\* It is important to thaw until last ice crystal and immediately add to MRC5 media

- 6 From the conical tube with 10mL of cold MRC5, using a P1000, add 200ul of MRC5 media to the sub stock aliquot, slowly mix and then add the entire aliquot to 10mL of MRC5 media.
- 7 Mix gently with 10mL serological pipette on slow settings.
- 8 Spin at 250g(rcf) for 10 min.
- 9 Discard supernatant and wash pellet with 1mL cold MRC5 media.
- 10 Add an additional 10mL cold MRC5 media and spin at 250g for 10 min.
- 11 10. Label T75 flask with your initials, date, and Px+1(passage).
- 12 11. Resuspend pellet in 10mL of cold MRC5 media and add to T75 flask.
- 13 12. Make sure to distribute the cells evenly by gently moving the flask up and down, side to side, and put in appropriate incubator.
- 14 13. Check on media and fibroblasts every 1-2 days, making sure they are growing.
- 15 14. Split fibroblasts in approximately 5 days.
- 16 15. Make sure to not overgrow and keep until passage 10.