



2 ▼

Nov 20, 2021

## 🌐 Genotyping Arabidopsis T-DNA lines V.2

**Steven J Burgess**<sup>1</sup><sup>1</sup>University of Illinois at Urbana-Champaign

1

[dx.doi.org/10.17504/protocols.io.bw5npg5e](https://dx.doi.org/10.17504/protocols.io.bw5npg5e)**Burgess Lab UIUC****Steven Burgess**

University of Illinois at Urbana-Champaign

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to [protocols.io](https://protocols.io) is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with [protocols.io](https://protocols.io), can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

This protocol is used for genotyping Arabidopsis seedlings to test for the presence of a transfer DNA (T-DNA) insertion. By using two primer sets it is possible to determine whether a seedling is homozygous, heterozygous or azygous for an insertion in the predicted genomic location.

To identify lines with T-DNA insertions in a gene of interest, you need the Arabidopsis Genome Identifier (AGI) number corresponding to the genomic locus (e.g. RCS1A = AT1G67090), then visit the [Salk Institute T-DNA Express site](#) to find all the mapped insertions at your locus of interest.

Genotyping primers have been pre-designed for each T-DNA line, these can be retrieved from the [Salk Institute T-DNA primer site](#), and ordered at any supplier of DNA oligonucleotides before starting the protocol.

In the US T-DNA lines can be purchased from the [Arabidopsis Biological Resource Center \(ABRC\)](#) and in the UK and EU from the [European Arabidopsis Stock Center \(NASC\)](#).

#### Recommended reading

- <http://signal.salk.edu/tdnaprimers.2.html>

#### Setting up the PCR reaction

Genotyping is performed with the Phire Direct PCR Mix, this includes the polymerase, nucleotides and salts necessary for amplification. We use the “dilution protocol” which involves taking a small leaf disk and homogenizing it in dilution buffer using a gel tip (see manufacturer’s instructions for more details.)

DOI

[dx.doi.org/10.17504/protocols.io.bw5npg5e](https://doi.org/10.17504/protocols.io.bw5npg5e)

Steven J Burgess 2021. Genotyping Arabidopsis T-DNA lines. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.bw5npg5e>



protocol

DOI: 10.1126/science.1086391

Arabidopsis, Genotyping

\_\_\_\_\_ protocol ,

Aug 04, 2021

Nov 20, 2021

52110

- Phire Direct Master Mix (Thermo Fisher Scientific; F170S)
- Genotyping primers (LBb1.3 ATT TTG CCG ATT TCG GAA C)
- Swiss Line Core Sampling Tool 0.5 mm, 0.50 mm I.D., 0.80 mm O.D. (Fisher Scientific; NC1310089) / (Electron Microscopy Sciences; 69039-05)
- 100bp ladder (New England Biolabs; N0467S)
- 10x Tris-Acetate-EDTA (TAE) Buffer (see [Sigma website for recipe](#))
- 10,000x SYBR<sup>TM</sup> Safe DNA stain (Thermo Fisher Scientific; S33102)
- 6x Gel Loading Dye (New England Biolabs; B7024S)

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to [protocols.io](#) is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with [protocols.io](#), can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

#### Prepare primer working solution

- 1 Re-suspend lyophilized primers in dH<sub>2</sub>O to a stock concentration of 100 mM. Note: primer sequences can be obtained from SALK T-DNA express if you have the T-DNA accession number (<http://signal.salk.edu/tdnaprimers.2.html>)
- 2 Create a 10 mM working solution of the primer by diluting 10 µL of the stock into 90 µL of water.

#### Preparing the PCR template

- 3 Using a 0.5 mm sampling device take a single leaf punch and add to 20 µL of dilution buffer.
- 4 5 µL of DNARelease Additive

- 5 Mash the sample with the end of a pipette tip
- 6 Incubate sample at room temperature for 5 mins
- 7 Incubate sample at 98 °C for 2 mins
- 8 Spin down sample for 1 min max speed in a mini centrifuge

#### Setting up the PCR reaction

- 9 Add the following components to a PCR tube. Set up one reaction with the diluted template from the putative T-DNA insertion line and one with WT arabidopsis as a negative control. The reaction will result in amplification if the T-DNA insert is present.

A	B
Component	Amount
2x Phire Master Mix	10 µL
Diluted template	1 µL
LBb1.3 primer [10 mM]	1 µL
RP [10 mM]	1 µL
dH2O	7 µL

Note: LBb1.3 (ATTTTGCCGATTTTCGGAAC) is the recommended sequencing primer for SALK lines. If you have SAIL, Wisc\_Lox or GABI-KAT lines you will need to use a different border primer

- 10 In a second PCR tube set up the following reaction. Set up one reaction for the putative T-DNA insertion line and one for WT arabidopsis as a positive control. This reaction tests for the absence of a TDNA insertion. If a T-DNA insertion is present there will be no amplification product as the distance between the L and R primers. The WT control should give a product, and heterozygous lines will also give a product as there is one WT allele.

A	B
Component	Amount
2x Phire Master Mix	10 µL
Diluted template	1 µL
LP [10 mM]	1 µL
RP [10 mM]	1 µL
dH2O	7 µL

## 11 Run the following PCR program

A	B	C	D
Step	Temperature	Time	Cycles
Initial denaturation	98	5 min	1
Denaturation	98	5s	40
Annealing	*	5s	
Extension	72	20s	
Final Extension	72	1 min	1

\*indicates annealing temperature to be adjusted to primer

## 12 Analyze the resulting products by gel electrophoresis