



Nov 18, 2020

Isolation and Fixation of Nuclei from the Mouse Brain for Dip-C

Longzhi Tan¹¹Stanford University

1 Works for me dx.doi.org/10.17504/protocols.io.bpsxmfn



Longzhi Tan
Stanford University

DOI

dx.doi.org/10.17504/protocols.io.bpsxmfn

PROTOCOL CITATION

Longzhi Tan 2020. Isolation and Fixation of Nuclei from the Mouse Brain for Dip-C. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bpsxmfn>

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Nov 17, 2020

LAST MODIFIED

Nov 18, 2020

PROTOCOL INTEGER ID

44599

Reagents

1 Prepare 1.5 M sucrose (40 mL):

☒ [sucrose Sigma](#)

- 20.538 g Aldrich Catalog #84097 (final: 1.5 Molarity (M) , 51.3 Mass / % volume)
- 40 mL water
- Heat and vortex to mix.
- Filter to sterilize.
- Store at 4 °C .

2 Prepare Nuclei Isolation Medium 1 (45 mL ; note that Tris is replaced with HEPES):

- 7.5 mL 1.5 M sucrose (final: 250 Milimolar (mM) , 8.56 Mass / % volume)
- 562.5 µl [KCl \(2 M\), RNase-free Thermo](#) (final: 25 Milimolar (mM))
- 450 µl [HEPES \(1 M\) Thermo](#) (final: 10 Milimolar (mM))
- 450 µl Fisher Catalog #15630080

1M MgCl₂ Invitrogen - Thermo

- 225 µl Fisher Catalog #AM9530G
- 36.2625 mL water
- Vortex to mix.
- Store at 4 °C .

(final: 5 Millimolar (mM))

3 Prepare 1% BSA in PBS (1.2 mL per sample; recipe below for 10 mL):

- 0.1 g
- Bovine Serum Albumin (BSA) Fraction V—Molecular Biology Grade Gemini Bio-Products Catalog #700-106P

PBS, pH 7.4 Thermo

- 10 mL Fisher Catalog #10010023
- Vortex to mix.
- Rotate until fully dissolved.
- Chill On ice .
- Store at -20 °C if needed.

Nuclei Isolation

4 KIMBLE 2mL Glass Dounce Tissue Grinder

Chill a Set Sigma Catalog #D8938

(or larger sizes:

KIMBLE Dounce tissue grinder set 7 mL complete Sigma

Aldrich Catalog #D9063

KIMBLE Dounce tissue grinder set 15 mL complete Sigma

Aldrich Catalog #D9938

KIMBLE Dounce tissue grinder set 40 mL complete Sigma

Aldrich Catalog #D9188

KIMBLE Dounce tissue grinder set 100 mL complete Sigma

Aldrich Catalog #D0189

, or

) On ice .

5 Prepare 1 mM DTT:

- 1 mL water
- 1M DL-Dithiothreitol solution (DTT) Sigma
- 1 µl Aldrich Catalog #646563
- -20 °C)
- Vortex to mix.


(aliquoted and stored at

6 Prepare Nuclei Isolation Buffer without Triton (6 mL per sample):

- 6 mL Nuclei Isolation Medium 1
- 6 µl 1 mM DTT (final: 5 Micromolar (µM))
- Vortex to mix.
- Chill On ice .


7 Prepare Nuclei Isolation Buffer with Triton (2 mL per sample):

-  **2 mL** Nuclei Isolation Buffer without Triton

 Triton X-100, 10% solution

-  **20 µl** Sigma Catalog #93443

(final: $[M]0.1\%$ (v/v))

- Vortex to mix.
- Chill  **On ice** .


8 Add  **2 mL** ice-cold Nuclei Isolation Buffer with Triton to the homogenizer.

9 Dounce the tissue with 5 strokes of the loose pestle (A), and 15 strokes of the tight pestle (B).

10 Transfer the homogenate to a tube.


11 Centrifuge at  **100 x g, 4°C, 00:08:00** .

12 Carefully remove supernatant without disrupting the soft pellet.

13 Resuspend in  **2 mL** Nuclei Isolation Buffer without Triton.

14 Centrifuge at  **100 x g, 4°C, 00:08:00** .

15 Carefully remove supernatant without disrupting the soft pellet.

16 Resuspend in  **2 mL** Nuclei Isolation Buffer without Triton.

17  Falcon 40 µm Cell

Filter with Strainer Corning Catalog #352340

or

 Corning™ Falcon™ Test Tube with 35µm Cell Strainer Snap

Cap Corning Catalog #352235

Fixation 10m

18 To every  **1 mL** of cells, add  **66.7 µl**

10m

19 Rotate at  **Room temperature** for  **00:10:00** .

20 Add  **200 µl** 1% BSA in PBS.

21 Invert to mix.

22 Centrifuge at  **1000 x g, 4°C, 00:05:00** .

23 Remove supernatant.

24 Resuspend in  **1 mL** ice-cold 1% BSA in PBS.


25  **C-Chip disposable**

Count with hemacytometer **INCYTO Catalog #DHC-N01**
cells per tube.

, and aliquot to up to 0.5 million

26 Centrifuge at  **1000 x g, 4°C, 00:05:00** .

27 Remove supernatant.

28 Store at  **-80 °C** .