



Oct 14, 2022

Radiolabeled polyamine uptake in cells

Marine Houdou¹, Nathalie Jacobs¹, Peter Vangheluwe¹¹KU Leuven1 *Works for me* Sharedx.doi.org/10.17504/protocols.io.yxmvm2x85g3p/v1 Nathalie Jacobs

ABSTRACT

This protocol provides a technique to determine the radiolabeled polyamine uptake capacity in cells, via the acquisition of disintegrations per minute (DPM) using a Liquid Scintillation Counter.

DOI

dx.doi.org/10.17504/protocols.io.yxmvm2x85g3p/v1

PROTOCOL CITATION

Marine Houdou, Nathalie Jacobs, Peter Vangheluwe 2022. Radiolabeled polyamine uptake in cells. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.yxmvm2x85g3p/v1>



LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Oct 04, 2022

LAST MODIFIED

Oct 14, 2022

PROTOCOL INTEGER ID

70812

GUIDELINES

















Proper guidelines for working with radiolabeled materials should be followed at all times.





SAFETY WARNINGS

Radiation hazards

- 1 Cells are seeded in 12-well plates, such that 70-80% confluency is reached on the day of the assay. Seed out 2 'treatment' wells, 1 'quick wash' well and 1 'untreated' well per cell line and treatment dose.

A	B	C	D
treatment	treatment	quick wash	untreated

- 2 Remove the culturing medium in 'treatment wells' and add  **500 µL** of medium, containing the desired concentration of radiolabeled polyamines. Leave regular culturing medium in the 'quick wash' and 'untreated' wells.
- 3 Incubate  **37 °C**  **00:30:00** 30m
- 4 4 minutes before reaching the 30-minute mark, perform a Quick-wash step as follows: add  **500 µL** of medium, containing the desired concentration of radiolabeled polyamines in the 'quick wash' wells.  **00:00:00** immediately remove the treatment medium, and wash with  **500 µL** of DPBS (-/-) containing  **50 micromolar (µM)** of the respective unlabeled (=cold) polyamine at  **4 °C** . Continue with 2 more washing steps with  **800 µL** of DPBS (-/-)  **4 °C** .
- 5 Proceed with washing the other wells ('treatment' and 'untreated'). Remove the treatment medium, and wash with cold  **500 µL** of DPBS (-/-) containing  **50 micromolar (µM)** of the respective cold polyamine at  **4 °C** . Continue with 2 more washing steps with  **800 µL** of DPBS (-/-) at  **4 °C** .
- 6 After the last wash, add  **200 µL** 0.1% SDS-DPBS (-/-) per well to lyse the cells.

- 7 Incubate  **00:10:00** at  **Room temperature** .
- 8 Scrape the cells and pipette the lysates into scintillation vials, that are filled with  **7 mL** of scintillation fluid (Ecolite (+), MP). Rinse each well with  **200 µL** DPBS (-/-) and pipette into the respective scintillation vial.
- 9 Mix scintillation vials well prior acquisition of disintegrations per minute (DPM) in the Liquid Scintillation Counter.