




SEP 18, 2023

## Whole Organoids Harvesting Procedure (Cultrex-modified)

Gabriela Vallejo Flores<sup>1</sup>,

Annika Fendler<sup>1</sup>

<sup>1</sup>Charite



Gabriela Vallejo Flores

OPEN ACCESS



**Protocol Citation:** Gabriela Vallejo Flores, Annika Fendler 2023. Whole Organoids Harvesting Procedure (Cultrex-modified).

**protocols.io**  
<https://protocols.io/view/whole-organoids-harvesting-procedure-cultrex-modif-cz4cx8sw>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** In development  
 We are still developing and optimizing this protocol

**Created:** Sep 18, 2023

**Last Modified:** Sep 18, 2023

### ABSTRACT

This protocol is use for whole organoids isolation, after isolation organoids may be:

- Resuspended in basement membrane matrix for further organoid culture.
- Resuspended in freezing medium for cryopreservation.
- Processed for biochemical analysis (such as RT-PCR, MS-PCR, sequencing, Western Blot, ELISA, or IHC)

### MATERIALS



Cultrex Organoid Harvesting Solution Bio-Techne Catalog #3700-100-01



Corning® Matrigel® Growth Factor Reduced (GFR) Basement Membrane Matrix Corning Catalog #356231



Bovine Serum Albumin Merck MilliporeSigma (Sigma-Aldrich) Catalog #A9418



1X Dulbecco's Phosphate Buffered Saline (DPBS) Thermo Fisher Scientific Catalog #14190094

- 24 multiwell plate
- 48 multiwell plate
- 15 ml tubes
- Ice

### BEFORE START INSTRUCTIONS

- Pre-warm seeding plate
- work with 1% BSA/PBS coated tubes and pipettes
- use pre-cooled pipette tips for matrigel
- precool centrifuge
- precool PBS and Cultrex™ Organoid Harvesting Solution
- Working on ice

**Keywords:** Whole,  
 Organoids, Harvesting,  
 Cultrex

**Sample description**

1

	Sample-ID	Purpose/Wells/Passage/Medium	Wells/Passage/Medium/Remarks

Table 1. Sample description

**Organoids harvesting**

2

Aspirate cell culture media and gently wash each well with 10 volumes cold (2-8 °C) PBS (Table 2). Be careful not to disrupt the basement membrane matrix containing organoids.




A	B	C
96-well plate	5µl	50µl
48-well plate	25µl	250µl
24-well plate	50µl	500µl

Table 2.


3

Add 10 volumes of cold (2-8 °C) Cultrex™ Organoid Harvesting Solution to each well (Table 2). Dislodge the dome with the pipett and transfer to a 15 ml Falcon.


4

Incubate for  01:00:00 to  01:30:00 , shaking at  40 rpm, 2-8°C on IKA rocker 3D on a cold pack in a Polystyrene box. This incubation is complete when the basement membrane matrix is no longer visible.

2h 30m

5 Centrifuge the tube at  500 x g, 2-8°C, 00:05:00 in a swinging bucket rotor to pellet the organoids. Aspirate the supernatant.

5m

6 Wash organoids with 10 volumes of cold (2-8 °C) PBS (Table 2), centrifuge at  500 x g, 2-8°C, 00:05:00 in a swinging bucket rotor to pellet the organoids. Aspirate the PBS.

5m

## Cell seeding Mullenders

15m

7 Pre-warm a 48, 24 or 6 wells plate and add to the cell pellet the amount of matrigel described in Table 3. Homogenize the matrigel with the cells and dropped in the center of the well.

	A	B	C	D
	48	25µl	5x10 <sup>4</sup>	250µl
	24	50µl	2.5x10 <sup>5</sup>	500µl
	6	5 drops of 50µl	2.5x10 <sup>5</sup> per drop	2500 µl

Table 3.

8 Flip the plate and incubate  00:15:00 at  37 °C 5%CO<sub>2</sub>.

15m

9 Add  250 µL ,  500 µL or  2500 µL of medium per well.

10

11

