

6



Aug 23, 2020

© ELISA for measurement of serum macrophage migration inhibitory factor (MIF).

Angel A Justiz-Vaillant¹, Belkis Ferrer-Cosme²

¹University of the West Indies St. Augustine; ²"Saturnino Lora Torres" Provincial Teaching Clinical Surgical Hospital. Cuba

1	Works for me	dx.doi.org/10.17504/protocols.io.bj7ekrje				
Uni	versity of the We	est Indies	angel.vaillant@sta.uwi.edu			
	Angel Justiz-V	aillant				

DOI

dx.doi.org/10.17504/protocols.io.bj7ekrje

University of the West Indies St. Augustine

PROTOCOL CITATION

Angel A Justiz-Vaillant, Belkis Ferrer-Cosme 2020. ELISA for measurement of serum macrophage migration inhibitory factor (MIF). . **protocols.io**

https://dx.doi.org/10.17504/protocols.io.bj7ekrje

C	Е	N	S	E

This is an open access protocol distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Aug 23, 2020

LAST MODIFIED

Aug 23, 2020

PROTOCOL INTEGER ID

40902

- An anti-human MIF coating antibody is adsorbed onto the microwells by incubation overnight at 4°C with carbonatebicarbonate buffer.
- 2 Add 50 μl of human serum or plasma. Human MIF present in the serum or plasma binds to antibodies adsorbed into the microwells.
- 3 The microplate is blocked with 3% non-fat milk-PBS buffer and later wash to remove unbound proteins.
- 4 Fifty (50) µl of biotin-conjugated anti-MIF antibody is added. The optimal dilution must be investigated.
- The microplate is rewashed with PBS-Tween 20 buffer, pH 7.4.

Citation: Angel A Justiz-Vaillant, Belkis Ferrer-Cosme (08/23/2020). ELISA for measurement of serum macrophage migration inhibitory factor (MIF).ÃÂ . https://dx.doi.org/10.17504/protocols.io.bj7ekrje

6	One hundred μ I of streptavidin-HRP conjugate is added and it binds to the biotin-conjugated anti-MIF antibody. The optimal dilution of this conjugate must be investigated.
7	The plate is washed following incubation to remove the unbound Streptavidin-HRP.
8	Add 100 µl of 3',3',5',5'- tetramethylbenzidine (TMB; Sigma-Aldrich) into each well.
9	Incubate the microwells in the dark for 20 min.
10	A colored product is formed in proportion to the quantity of MIF present in the sample or standard.
11	The reaction is terminated by addition of 100 $\mu l3MH2SO4$ and the absorbance is measured at 450 nm.
12	A standard curve is made from 7 human MIF standard dilutions and the human MIF sample concentration is determined.
13	For better results place the microplate on a microplate shaker in every incubation.