



JAN 28, 2024

OPEN ACCESS



DOI:

dx.doi.org/10.17504/protocols.io.ewov1qbrogr2/v1

Protocol Citation: Chakrakodi N Varun, Dhanya K, Ashwini M Ananda, Reeta Mani 2024. Rabies Virus Sequencing using Illumina- MiSeq. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.ewov1qbrogr2/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Jan 20, 2024

🌐 Rabies Virus Sequencing using Illumina- MiSeq

Chakrakodi N Varun¹, Dhanya K¹, Ashwini M Ananda¹, Reeta Mani¹

¹Department of Neurovirology, NIMHANS

Department of Neurovirology, National Institute of Mental Health and Neuro Sciences



Chakrakodi N Varun

Department of Neurovirology, NIMHANS

ABSTRACT

This Rabies whole genome sequencing protocol has been derived and modified from the Illumina COVIDSeq RUO sequencing pipeline. The protocol has been modified and optimised for sequencing Rabies virus (RABV). The methodology uses RABV-specific primers that have been designed in-house. In brief, the RNA is extracted from samples and converted to cDNA. RABV sequencing library is generated and sequenced using MiSeq.

GUIDELINES

Ensure all Biosafety guidelines are followed as per the working policy.

Please ensure you keep a log of sample processing and use a template to locate sample wells.

The protocol is based on sequencing a total of 96 samples, which includes a known positive and negative control. At least one negative control is mandatory to assess possible contamination. Positive control with a known sequence is mandatory to estimate the quality of the run.

The complete library preparation process takes approximately two days at our laboratory. It is advisable to use "Safe Stops" as needed.

Last Modified: Jan 28, 2024

PROTOCOL integer ID: 93845

Keywords: Sequencing, Rabies Virus, Lyssavirus Rabies, Illumina MiSeq, NGS

MATERIALS

The following reagents available from Illumina have been used in this protocol.

- EPH3 (Elution Prime Fragment 3HC Mix)
- FSM (First Strand Mix)
- RVT (Reverse Transcriptase)
- IPM (Illumina PCR Mix)
- EBLTS (Enrichment BLT)
- TB1 (Tagmentation Buffer 1)
- Nuclease-free water (NFW)
- ST2 (Stop Tagment Buffer 2)
- TWB (Tagmentation Wash Buffer)
- EPM (Enhanced PCR Mix)
- Index adapters (Illumina-PCR Indexes)
- ITB (Illumina Tune Beads)
- Ethanol
- RSB (Resuspension Buffer)

SAFETY WARNINGS

! Use appropriate PPE as needed.

BEFORE START INSTRUCTIONS

All the processes should be performed in Biosafety cabinets. Ensure you have separate Biosafety cabinets for RNA extraction, Reagent preparation, and Template addition.

Samples and Extraction



27m 30s

- 1
- Samples:**

The following samples (human or animal sources) may be used for viral RNA extraction for Rabies lyssavirus.



 - 1. Brain Tissue or Nuchal skin:** Homogenise the tissue by crushing a small piece in a sterile environment. Transfer the contents to a vial and vortex and spin down the sample. Retrieve the supernatant.
 - 2. Saliva**
- 2
- RNA Extraction:**




The RNA is extracted using the QIAamp Viral RNA Mini Kit, as per the procedure outlined in the QIAamp Viral RNA Mini Handbook (<https://www.qiagen.com/us/resources/download.aspx?id=c80685c0-4103-49ea-aa72-8989420e3018&lang=en>).

2.1 Add  140 µL of the supernatant or saliva sample to the Buffer AVL–carrier RNA in the microcentrifuge tube. Mix by pulse-vortexing for  00:00:15 . 15s

2.2 Incubate at  Room temperature for  00:20:00 . 20m


2.3 Briefly centrifuge the tube to remove drops from the inside of the lid.




2.4 Add  560 µL ethanol (96–100%) to the sample, and mix by pulse-vortexing for  00:00:15 . After mixing, briefly centrifuge the tube to remove drops from inside the lid. 15s




2.5 Carefully apply  630 µL of the solution from step 2.4 to the QIAamp Mini column (in a 2 mL collection tube) without wetting the rim. Close the cap, and centrifuge at  8000 rpm for  00:01:00 . Place the QIAamp Mini column into a clean 2 mL collection tube, and discard the tube containing the filtrate. 1m




2.6 Carefully open the QIAamp Mini column, and repeat the step.

Note


If the original sample volume exceeds  140 µL , repeat this step until the lysate has been loaded onto the spin column.




2.7 Carefully open the QIAamp Mini column, and add  500 µL Buffer AW1. Close the cap, and centrifuge at  8000 rpm for  00:01:00. Place the QIAamp Mini column in a clean 2 mL collection tube (provided), and discard the tube containing the filtrate. 1m

2.8 Carefully open the QIAamp Mini column, and add  500 µL Buffer AW2. Close the cap and centrifuge at full speed  14.000 rpm for  00:03:00. 3m

2.9 Place the QIAamp Mini column in a clean 1.5 mL microcentrifuge tube. Discard the old collection tube containing the filtrate. Carefully open the QIAamp Mini column and add  60 µL of Buffer AVE. Close the cap and incubate at  Room temperature for  00:01:00. 1m

Note

If a very low viral load is expected, such as in a salivary sample, it is advised to elute in approximately  20 µL.

2.10 Centrifuge at  8000 rpm for  00:01:00. The eluted RNA can be used immediately or stored at ≤  -80 °C. 1m

Note

It is advised to check the quality of extracted RNA using Nanodrop and run an RABV RT-PCR to determine the Ct value. Sequencing does not work well if extraction quality is low or the Ct is > 30.

cDNA Conversion

25m

3 The extracted RNA is annealed using random hexamers to prepare for cDNA synthesis during this process. 5m

Thaw the EPH3 (Elution Prime Fragment 3HC Mix) at  Room temperature.

Note


Use nuclease-free water as a Negative Control in one or more wells depending on the user's requirement.

Label a new PCR plate as **CDNA**.

Add 8.5 µl EPH3 to each well.

Add 8.5 µl eluted sample to each well.


Seal and shake at 1600 rpm for  00:01:00 .

Centrifuge at 1000 rpm for  00:01:00 .

Set up a PCR (RABV Anneal program) as follows:

 65 °C for  00:03:00

Hold at  4 °C

- 4 This step reverse transcribes the RNA fragments primed with random hexamers into first-strand cDNA using  22m reverse transcriptase.


Thaw the FSM (First Strand Mix) and RVT (Reverse Transcriptase) reagents in  Room temperature

For 96 samples, prepare a master mix in a 1.7 ml tube as follows.

FSM-  720 µL

RVT-  80 µL


Add 8 µl master mix to each well of the CDNA plate.

Seal and shake at 1600 rpm for  00:01:00 .

Centrifuge at 1000 rpm for  00:01:00 .

Set up a PCR (RABV FSM program) as follows:

Choose the preheat lid option

Set the reaction volume to  25 µL

 25 °C  00:05:00

50 °C

00:10:00

80 °C

00:05:00

Hold at 4 °C

Targeted Amplification of cDNA

10m 15s

5 Preparation of RABV Primer Pools

- 5.1** Prepare the primer pool mix as provided in the RABV_PrimerPool sheets. Store the Rabies Primer Pools (Odd and Even Primer mix) at -20 °C

RABV_PrimerPool.xlsx 12KB

Note

1. The RABV_Primer Pools were designed using Primal Scheme (<https://primalscheme.com/>)
2. The RABV_PrimerPool.xlsx contains two sheets, for odd and even primal pool mix. The quantity of 100 micromolar (μM) stock solution to be taken for generating a primer pool and the resulting final concentration in (uM) is provided in the the excel sheet.

- 5.2** Prepare two master mixes as follows.

1. Odd Rabies PrimerPool Master Mix

IPM (Illumina PCR Mix): 1260 μL

Odd PrimerPool Mix: 361.20 μL

Nuclease Free Water: 394.8 μL

2. Even Rabies PrimerPool Master Mix

IPM (Illumina PCR Mix): 1260 µL
 Odd PrimerPool Mix: 361.20 µL
 Nuclease Free Water: 394.8 µL

Note

The Master mix is calculated to accommodate 96 reactions. If using lesser samples calculate accordingly.

5.3 Amplification PCR

10m 15s

Label two PCR plates as

- (i) RABV_Odd plate
- (ii) RABV_Even Plate

Add 20 µL of Odd Rabies PrimerPool Master Mix to each well of the RABV_Odd plate

Add 5 µL of cDNA synthesised in the previous step to the corresponding well of the RABV_Odd plate.

Add 20 µL of Even Rabies PrimerPool Master Mix to each well of the RABV_Even plate

Add 5 µL of cDNA synthesised in the previous step to the corresponding well of the RABV_even plate.

Seal and shake at 1600 rpm for 00:01:00 .

Centrifuge at 1000 rpm for 00:01:00 .

Set up two PCR's (RABV amplification program) as follows:

Choose the preheat lid option

Set the reaction volume to 25 µL

98 °C for 00:03:00

35 cycles of:

98 °C for 00:00:15

🌡️ 63 °C for ⌚ 00:05:00
Hold at 🌡️ 4 °C

Note

If you are stopping, seal the plate and store at 🌡️ -20 °C for up to 3 days

Tagment of PCR Amplicons

7m

- 6** This step uses EBLTS (Enrichment Bead-Linked Transposomes) to tagment PCR amplicons, which is a process that fragments and tags the PCR amplicons with adapter sequences.

- 6.1** Thaw EBLTS and TB1 Buffer at 🌡️ Room temperature

Note

If the RABV_Odd plate and RABV_Even Plate were stored thaw at 🌡️ Room temperature, shake the plates at 1600 rpm for 1 minute and centrifuge at 1000 x g for 1 minute before starting.

Label a new PCR plate as TAG.

Transfer 🧴 10 µL from each well of the RABV_Odd plate to the corresponding well of the TAG plate.

Transfer 🧴 10 µL from each well of the RABV_Even Plate to the corresponding well of the TAG plate

- 6.2** Prepare **Tagmentation Master Mix** in a 15 ml tube as follows.

2m

TB1 🧴 1008 µL

EBLTS 🧴 336 µL

Nuclease Free Water 🧴 1680 µL

Note

Ensure that the beads are uniformly mixed before use. Pipette mix if needed and pulse centrifuge before use.

Add 30 µl master mix to each well in the TAG plate.

Seal and shake at 1600 rpm for 00:01:00

Centrifuge at 1000 rpm for 00:01:00

6.3 Set up PCR (TAG program) as follows:

5m

Choose the preheat lid option

Set the reaction volume to 50 µl

55 °C for 00:05:00

Hold at 10 °C

Post Tagmentation Clean Up

12m

7 This step washes the adapter-tagged amplicons before PCR amplification.

7.1 Vortex ST2 (Stop Tagment Buffer 2) and TWB (Tagmentation Wash Buffer) before use.

10m

Centrifuge the TAG plate at 500 x g for 00:01:00 .

Add 10 µL ST2 to each well of the TAG plate.

Seal and shake at 1600 rpm for 1 minute.

Incubate at Room temperature for 00:05:00 .

Centrifuge at for 00:01:00 .

Place on the magnetic stand and wait until the liquid is clear (00:03:00).

Note


If the liquid is not clear at this point continue to place it on the magnetic stand for another ~2 minutes. Inspect for bubbles on the seal. If present, centrifuge at 500 x g for 1 minute, and then place on the magnetic stand (~3 minutes).

7.2 Wash the beads as follows:

2m

Remove from the magnetic stand.

Add  100 μ L TWB to each well.

Seal and shake at 1600 rpm for  00:01:00 .

Centrifuge  for  00:01:00 .

Place on the magnetic stand and wait until the liquid is clear (~3 minutes).

Note

If the liquid is not clear at this point continue to place it on the magnetic stand for another ~2 minutes. Inspect for bubbles on the seal. If present, centrifuge at 500 x g for 1 minute, and then place on the magnetic stand (~3 minutes).

Remove and discard all supernatant from each well.

Wash beads a second time.

Leave supernatant in the plate for the second wash to prevent beads from overdrying.

Amplify Tagmented Amplicons and Indexing

11m 50s

- 8** This step amplifies the tagmented amplicons using a PCR program. The PCR step adds prepared 10 base pair Index 1 (i7) adapters, Index 2 (i5) adapters, and sequences required for sequencing

8.1 Prepare the Master mix as follows, in a 15 ml tube.

EPM (Enhanced PCR Mix):  2016 μ L

Nuclease-free water:  2016 μ L


Note


Centrifuge the TAG plate, keep the on magnetic stand and remove remaining TWB. Do not leave any residual TWB in the wells.


8.2

Add  40 μL PCR Master Mix to each well.

11m 50s

Add  10 μL index adapters to each well of the PCR plate.


Seal and shake at 1600 rpm for  00:01:00

If the liquid is visible on the seal, centrifuge at  500 x g for 1 minute.

Inspect to make sure beads are resuspended.

Set up a PCR (Amplification and Indexing PCR) as follows:



Choose the preheat lid option and set to  100 $^{\circ}\text{C}$

Set the reaction volume to  50 μL

 72 $^{\circ}\text{C}$ for  00:03:00

 98 $^{\circ}\text{C}$ for  00:03:00

7 cycles of:

 98 $^{\circ}\text{C}$ for  00:00:20

 60 $^{\circ}\text{C}$ for  00:00:30

 72 $^{\circ}\text{C}$ for  00:01:00

 72 $^{\circ}\text{C}$ for  00:03:00

Hold at  10 $^{\circ}\text{C}$

Pool and Clean Up Libraries

8m 30s


- 9 This step combines libraries from each 96-well sample plate into one 1.7 ml tube. Libraries of optimal size are then bound to magnetic beads, and fragments that are too small or large are washed away.

9.1




Centrifuge the TAG plate at  500 x g for  00:01:00 .



6m

Place on the magnetic stand and wait until the liquid is clear (~3 minutes).

Transfer  5 μL library from each well of the TAG plate into a 1.7 ml tube.
Vortex the tubes to mix, and then centrifuge briefly.
Add 0.9x of IPB.

Note




Assuming that  5 μL from 96 wells were pooled, give a total volume of  480 μL add  432 μL of IPB (480x 0.9).

Vortex to mix.
Incubate at  Room temperature for  00:05:00 .
Centrifuge briefly.

Place on the magnetic stand and wait until the liquid is clear (~5 minutes).
Remove and discard all supernatant.

9.2 Wash beads as follows.

30s

Keep on the magnetic stand and add  1000 μL fresh 80% Ethanol.
Incubate at  Room temperature  00:00:30 .
Remove and discard all supernatant.


Wash beads a second time.
Centrifuge briefly.
Use a 20 μL pipette to remove all residual EtOH.


Note

The final

9.3 Add 55 μL Resuspension Buffer (RSB)

2m

Vortex to mix, and then centrifuge briefly.
Incubate at room temperature for  00:02:00 .
Place on the magnetic stand and wait until the liquid is clear (~2 minutes).

Transfer  50 μ L supernatant to a fresh new microcentrifuge tube.

Note

The final library can be at  -20 °C for up to 30 days.

We check the quality of the final library prepared using a Tapestation and the library is quantified using a Qubit Fluorometer.