



Jan 25, 2022

protein extraction from cell pallets

 protein extraction from cell palletsJacob Waldbauer¹, Lichun Zhang¹, Jacob Waldbauer¹, Lichun Zhang¹¹University of Chicago

1

dx.doi.org/10.17504/protocols.io.b4btqsnn

Lichun Zhang

Use LDS buffer to extract protein from cell pallets

DOI

dx.doi.org/10.17504/protocols.io.b4btqsnn

Jacob Waldbauer, Lichun Zhang, Jacob Waldbauer, Lichun Zhang 2022. protein extraction from cell pallets. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.b4btqsnn>



protein extraction from cell pallets, Lichun Zhang

 protocol ,

Jan 25, 2022

Jan 25, 2022

57427

- 1 Make 1X LDS buffer. In 40mL of 1X buffer, add Tris HCl 0.666g, Tris base 0.682g, LDS 0.8g, EDTA 0.006g, Glycerol 4g.
- 2 make 1M DTT stock solution in LC/MS water.
- 3 Add 200~500uL of 10mM DTT in 1X LDS buffer to the cell pallet tube.

- 4 Heat sample up at 95°C for 20 minutes.
- 5 Cool down the sample to room temperature. Meanwhile, make a 0.5M iodoacetamide stock solution.
- 6 Add IAM to 50mM/60mM in the sample, incubate in dark for 30-60minutes.
- 7 Add DTT to 50mM/60mM. Sample is ready to continue with FASP.