



Version 2 ▾

Jun 23, 2022

# nCoV-2019 McGill Artic PCR Protocol, V4.1 at 63C V.2

Version 1 is forked from [nCoV-2019 McGill Artic PCR Protocol, 5 ul RT and V3 only + LA1 at 63C](#)

Sarah J Reiling<sup>1</sup>, Kayleigh Loranger<sup>1</sup>, Anne-Marie Roy<sup>1</sup>, Shu-Huang Chen<sup>1</sup>, Ioannis Ragoussis<sup>1</sup>

<sup>1</sup>McGill University

1 Works for me

Share

[dx.doi.org/10.17504/protocols.io.ewov18e4ygr2/v2](https://dx.doi.org/10.17504/protocols.io.ewov18e4ygr2/v2)

McGill Genome Centre

Kayleigh Loranger

## ABSTRACT

This is the updated SARS-Cov-2 PCR Protocol, with the ARTIC V4.1 primers, that is currently being used at the McGill Genome Center.

## DOI

[dx.doi.org/10.17504/protocols.io.ewov18e4ygr2/v2](https://dx.doi.org/10.17504/protocols.io.ewov18e4ygr2/v2)

## PROTOCOL CITATION

Sarah J Reiling, Kayleigh Loranger, Anne-Marie Roy, Shu-Huang Chen, Ioannis Ragoussis 2022. nCoV-2019 McGill Artic PCR Protocol, V4.1 at 63C.

**protocols.io**

<https://dx.doi.org/10.17504/protocols.io.ewov18e4ygr2/v2>

Version created by Kayleigh Loranger



## FORK NOTE

## FORK FROM

Forked from [nCoV-2019 McGill Artic PCR Protocol, 5 ul RT and V3 only + LA1 at 63C](#), Sarah Reiling

## LICENSE

———— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Jun 23, 2022

## LAST MODIFIED

Jun 23, 2022

## PROTOCOL INTEGER ID

65186

## MATERIALS TEXT

### MATERIALS

[Q5 High-Fidelity 2X Master Mix - 500 rxns](#) **New England**

**Biolabs Catalog #M0492L**

[nuclease-free water](#) **Contributed by users**

[Fresh 80% Ethanol](#) **Contributed by users**

[Quant-iT™ PicoGreen™ dsDNA Assay Kit](#) **Invitrogen - Thermo**

**Fisher Catalog #P11496**

[AmpureXP beads](#) **Beckman**

**Coulter Catalog #A63880**

## Primer pool preparation


### 1 PRIMER POOL PREPARATION

If required resuspend lyophilised primers at a concentration of 100 µM each

V4.1 only primers for this protocol were designed using [Primal Scheme](#) and generate overlapping 400 nt amplicons. V4.1 was added to the V4 primer set for optimization. Primer names and dilutions are listed below.

[primer-schemes/nCoV-2019/V4.1 at master · artic-network/primer-schemes · GitHub](#)

For information on V4 primers visit [Optimization of the SARS-CoV-2 ARTIC Network V4 Primers and Whole Genome Sequencing Protocol - PMC \(nih.gov\)](#)

- 2 The V4.1 pre-pooled primers in  **1.5 mL** Eppendorf labelled tubes are labelled “Pool 1 (100µM)” or “Pool 2 (100µM)”. The primers do not require additional preparation.

If the V4.1 primers are not pre-pooled, follow the below pipetting scheme to make the master mix by adding the following primers to the V4 primer pools.

Added to pool 1:

SARS-CoV-2\_23\_RIGHT\_alt1  
SARS-CoV-2\_27\_RIGHT\_alt1  
SARS-CoV-2\_79\_RIGHT\_alt1  
SARS-CoV-2\_89\_LEFT\_alt1  
SARS-CoV-2\_89\_RIGHT\_alt1

Added to pool 2:

SARS-CoV-2\_10\_LEFT\_alt1  
SARS-CoV-2\_10\_RIGHT\_alt1  
SARS-CoV-2\_76\_LEFT\_alt1  
SARS-CoV-2\_76\_RIGHT\_alt1  
SARS-CoV-2\_88\_LEFT\_alt1  
SARS-CoV-2\_90\_RIGHT\_alt1

The guide to pooling volumes are as follows;

2x volume:

SARS-CoV-2\_1\_LEFT & SARS-CoV-2\_1\_RIGHT  
SARS-CoV-2\_7\_LEFT & SARS-CoV-2\_7\_RIGHT  
SARS-CoV-2\_13\_LEFT & SARS-CoV-2\_13\_RIGHT  
SARS-CoV-2\_17\_LEFT & SARS-CoV-2\_17\_RIGHT  
SARS-CoV-2\_27\_LEFT & SARS-CoV-2\_27\_RIGHT  
SARS-CoV-2\_45\_LEFT & SARS-CoV-2\_45\_RIGHT  
SARS-CoV-2\_59\_LEFT & SARS-CoV-2\_59\_RIGHT  
SARS-CoV-2\_60\_LEFT & SARS-CoV-2\_60\_RIGHT  
SARS-CoV-2\_61\_LEFT & SARS-CoV-2\_61\_RIGHT  
SARS-CoV-2\_64\_LEFT & SARS-CoV-2\_64\_RIGHT  
SARS-CoV-2\_79\_LEFT & SARS-CoV-2\_79\_RIGHT  
SARS-CoV-2\_90\_LEFT & SARS-CoV-2\_90\_RIGHT  
SARS-CoV-2\_91\_LEFT & SARS-CoV-2\_91\_RIGHT


1x volume: All other primers including alts (from V4.1).

Primers should be pooled in the **mastermix cabinet** which should be cleaned with decontamination wipes and UV sterilised before and after use.

## Multiplex PCR









3 

### MULTIPLEX PCR

In the **extraction and sample addition cabinet** add  **5 µL** RT product to each tube and mix well by pipetting.

The **extraction and sample addition cabinet** should be cleaned with decontamination wipes and UV sterilised before and after use.

4 In the mastermix hood set up the multiplex PCR reactions as follows in 0.2mL 8-strip PCR tubes:

Component [10 uM]	Pool 1 [10 uM primer]	Pool 2
Q5 Hot Start High-Fidelity 2X Master Mix	 <b>12.5 µL</b>	 <b>12.5 µL</b>
Primer Pool 1 or 2 (10µM pool 1+2)	 <b>3.7 µL</b>	 <b>3.7 µL</b>
Nuclease-free water	 <b>3.8 µL</b>	 <b>3.8 µL</b>
<b>Total</b>	 <b>20 µL</b>	 <b>20 µL</b>

Add 20ul of PCR mastermix to the **5 ul RT product** of step 10.

A PCR mastermix for each pool should be made up in the **mastermix cabinet** and aliquoted into PCR strip tubes. Tubes should be wiped down when entering and leaving the mastermix cabinet.

5 Pulse centrifuge the tubes to collect the contents at the bottom of the tube.

6 Set-up the following program on the thermal cycler:

Step	Temperature	Time	Cycles
Heat Activation	98 °C	00:00:30	1
Denaturation	98 °C	00:00:15	36
Annealing	63 °C	00:05:00	36
Hold	4 °C	Indefinite	1

Cycle number should be 25 for Ct 18-21 up to a maximum of 36 cycles for Ct 35

#### PCR clean-up

7

#### PCR CLEANUP

Combine the entire contents of "Pool 1" and "Pool 2" PCR reactions for each biological sample into to a single 1.5 mL Eppendorf tube.

8 Clean-up the amplicons using the following protocol:

Add an equal volume (1:1) of AmpureXP beads to the sample tube and mix by pipetting.

Incubate for 5 min at room temperature.

Pellet on magnet for 5 min. Remove supernatant.

Add 200 ul of 80% ethanol to the pellet and wash twice.

Let the beads dry for 3 min.

Add 30 ul elution buffer and resuspend the beads. Incubate for 3 minutes.

Pellet on magnet for 5 min. Remove and keep eluate (30 ul).

Amplicon clean-up should be performed in the **post-PCR cabinet** which should be cleaned with decontamination wipes and UV sterilised before and after use.

#### Amplicon Quantification and normalisation

9

#### AMPLICON QUANTIFICATION AND NORMALIZATION

Quantify the amplicon pools using a fluorimetric dsDNA assay. (e.g: PicoGreen with a standard curve 0-200ng)

We expect following concentrations:

**Pool 1+2 combined:**

100-150 ng/ul for Ct 14-24

30-80 ng/ul for Ct 25-29

10-30 ng/ul for Ct 30-36

10



***Nextera Flex Library Prep:***

After quantification of Pool 1+2, take a new plate and add 150 ng of Pool 1+2 and add up with nuclease-free water to a total volume of **30 ul (= 5 ng/ul)**.

***Nanopore Library Prep:***

After quantification of Pool 1+2, take a new plate and add 200 ng of Pool 1+2 and add up with nuclease-free water to a total volume of **20 ul (= 10 ng/ul)**.