

Aug 29, 2024

# Immunofluorescent staining of dopaminergic neurons in brains of *D. melanogaster*

DOI

[dx.doi.org/10.17504/protocols.io.eq2lyweyevx9/v1](https://dx.doi.org/10.17504/protocols.io.eq2lyweyevx9/v1)

Natalie Kaempf<sup>1,2,3</sup>, Uli Pech<sup>1,2</sup>, Patrik Verstreken<sup>1,2</sup>

<sup>1</sup>VIB-KU Leuven Center for Brain & Disease Research, 3000 Leuven, Belgium;

<sup>2</sup>KU Leuven, Department of Neurosciences, Leuven Brain Institute, 3000 Leuven, Belgium;

<sup>3</sup>Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD 20815

ASAP Collaborative Rese...



Natalie Kaempf

VIB KU Leuven

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.eq2lyweyevx9/v1](https://dx.doi.org/10.17504/protocols.io.eq2lyweyevx9/v1)

**Protocol Citation:** Natalie Kaempf, Uli Pech, Patrik Verstreken 2024. Immunofluorescent staining of dopaminergic neurons in brains of *D. melanogaster*. [protocols.io](https://dx.doi.org/10.17504/protocols.io.eq2lyweyevx9/v1) <https://dx.doi.org/10.17504/protocols.io.eq2lyweyevx9/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** July 17, 2024

**Last Modified:** August 29, 2024

**Protocol Integer ID:** 103595

**Keywords:** ASAPCRN, Immunohistochemistry, dopaminergic neurons, *Drosophila* brains



**Funders Acknowledgement:**

**Aligning Science Across**

**Parkinson's**

**Grant ID: ASAP-000430**

**EMBO long-term postdoctoral**

**fellowship**

**Grant ID: ALTF\_299-2019**

**Research project, FWO**

**Vlaanderen**

**Grant ID: G0A5219N**

**Research project, FWO**

**Vlaanderen**

**Grant ID: G0B8119N**

**Methusalem project**

**Grant ID: METH/21/05**

**(3M210778)**

**Research project, KU Leuven**

**Parkinson Fonds**

**Grant ID: EQZ-PARFON-O2010**

**Opening the Future grant,**

**Leuven Universiteitsfonds**

**(LUF)**

**Grant ID: EQZ-OPTFUP-O2010**

**Research project, FWO**

**Vlaanderen**

**Grant ID: G031324N**

## Abstract

This protocol is used to visualize for dopaminergic neurons in *Drosophila* brains, but can be adapted for other targets as well.

## Materials

Rabbit polyclonal anti-TH, Millipore, Cat# AB152, RRID: AB\_390204

Mouse monoclonal anti-DLG, DSHB, DSHB Cat# 4F3 anti-discs large; RRID: AB\_528203

Goat anti-Rabbit IgG Alexa Fluor<sup>TM</sup> 488, Life Technologies, Cat# A-11034 RRID: AB\_2576217

Goat anti-Mouse IgG Alexa Fluor<sup>TM</sup> 555, Life Technologies, Cat# A-21424 RRID: AB\_141780

RapiClear 1.47, Sunjin Lab, Cat# RC147001
















Triton X-100 Solution, Sigma Aldrich, Cat# 93443-100ML

Formaldehyde solution (37 wt. % in H<sub>2</sub>O) F1635-500ml, Sigma-Aldrich




## Dissection of fly brains, fixation, block and primary antibody incubation

1h 30m

- 1 Prepare fresh 3.7% paraformaldehyde in 1x PBS, 0.2% Triton X-100 (PBX),  500  $\mu$ L per genotype and store  On ice .
- 2 Dissect fly brains in ice-cold PBS under stereomicroscope, collect with a PBX coated pipette and transfer in the 3.7% paraformaldehyde solution.
- 2.1 Store the dissected fly brains in paraformaldehyde solution on ice, but do not keep fly brains longer than  01:00:00  On ice 1h
- 3 Incubate the dissected fly brains in paraformaldehyde solution on rotation wheel for  00:30:00 at  Room temperature 30m
- 4 Take the paraformaldehyde solution off and wash a first time with PBX 0.2% (  500  $\mu$ L ) and take it off right away
- 5 Wash 3 times with PBX 0.2% (  500  $\mu$ L ) for 15 min on a rotator at  Room temperature
- 6 Block with 10% NGS in PBX (  500  $\mu$ L ) for  01:00:00 on rotator at  Room temperature 1h
- 7 Add primary antibody (in  500  $\mu$ L of 10% NGS in PBX) put it on rotation wheel at  4  $^{\circ}$ C  48:00:00 2d
- 7.1 Use rabbit a-TH, Sigma, 1:200, and mouse a-DLG, DSHB, 1:100


## secondary antibody incubation




2d

- 8 Take the primary antibody solution off and wash a first time with PBX 0.2% (  500  $\mu$ L ) and take



it off right away

9 Wash 3 times with PBX 0.2% (  500 µL ) for 15 min on a rotator at  Room temperature


10 Add secondary antibody (in  500 µL of 10% NGS in PBX), cover with tinfoil and put it on rotation wheel at  4 °C  Overnight



2d

10.1 Use Goat anti-Rabbit IgG Alexa Fluor<sup>TM</sup> 488, Life Technologies, and Goat anti-Mouse IgG Alexa Fluor<sup>TM</sup> 555, Life Technologies, at 1:500 dilution.

## Mounting

15m


11 Take the secondary antibody solution off and wash a first time with PBX 0.2% (  500 µL ) and take it off right away

12 Wash 3 times with PBX 0.2% (  500 µL ) for 15 min on a rotator at  Room temperature

13 prepare slides:

13.1 clean slides with EtOH

13.2 place 2 book binder rings on top of each other

13.3 fill book binder rings with 0.075% PLL solution, leave minimum  00:15:00 , take the solution off and leave a thin layer to dry

15m

14 Mount fly brains by pipetting the brains next to the book binders, fill the book binder chamber with PBX and transfer the brains individually with forceps to the book binder chamber, orient the brains correctly with anterior facing up (antenna lobe up) and softly push down until they stick

15 remove PBX



- 16 add 8  $\mu$ l of mounting medium (RapiClear 1.47) by gently pipetting on top of brains
- 17 add glass coverslip, seal with nail polish and let cure overnight at 🌡️ 4 °C in dark