



2 ▼

Feb 07, 2022

A-Tailing with Taq Polymerase V.2

New England Biolabs¹¹New England Biolabs

1

dx.doi.org/10.17504/protocols.io.be6ajhae**New England Biolabs (NEB)**Tech. support phone: **+1(800)632-7799** email: **info@neb.com****New England Biolabs**
New England Biolabs

This protocol can be used to add As to the blunt-ends of DNA fragments that have been amplified using a high-fidelity polymerase (such as Q5® High Fidelity DNA Polymerase).

DOI

dx.doi.org/10.17504/protocols.io.be6ajhae<https://www.neb.com/protocols/2013/11/01/a-tailing-with-taq-polymerase>New England Biolabs 2022. A-Tailing with Taq Polymerase. **protocols.io**<https://dx.doi.org/10.17504/protocols.io.be6ajhae>

New England Biolabs



A-tail, Taq, Taq pol, Taq polymerase

_____ protocol ,

Apr 15, 2020

Feb 07, 2022

35746

MATERIALS

☒ [ThermoPol Reaction Buffer Pack - 6.0 ml](#) **New England****Biolabs Catalog #B9004S**☒ [Taq DNA Polymerase with ThermoPol Buffer - 400 units](#) **New England****Biolabs Catalog #M0267S**

Please refer to the Safety Data Sheets (SDS) for health and environmental hazards.

1 Clean-up the amplified DNA from the PCR components

This can be done by using a PCR-column purification protocol. This step is essential because the robust exonuclease activity associated with the high-fidelity enzyme will remove any untemplated nucleotides that are added by Taq DNA Polymerase.

2

Set-up the reaction by adding the following components:

A	B
REAGENT	AMOUNT
PCR-amplified DNA	X µl
10X ThermoPol® Buffer (NEB#B9004)	5 µl
1mM dATP	10 µl
Taq DNA Polymerase (NEB#M0267)	0.2 µl
H2O	X µl
<i>Total Reaction Volume</i>	<i>50 µl*</i>

*This volume can be adjusted based on the volume of PCR-amplified DNA that needs to be added).

3

Incubate the reaction at 🔥 **72 °C** for ⌚ **00:20:00** .