



Jul 20, 2022

🌐 Preparation of ATP13A2 microsomes from Sf9 cells

Sue Sim¹, eunyong_park¹¹University of California, Berkeley

1 Works for me

 Sharedx.doi.org/10.17504/protocols.io.kxygxzde4v8j/v1

Sue Sim

ABSTRACT

Isolate microsomes from Sf9 cells expressing ATP13A2

DOI

dx.doi.org/10.17504/protocols.io.kxygxzde4v8j/v1

PROTOCOL CITATION

Sue Sim, eunyong_park 2022. Preparation of ATP13A2 microsomes from Sf9 cells. **protocols.io**<https://dx.doi.org/10.17504/protocols.io.kxygxzde4v8j/v1>

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Jul 17, 2022

LAST MODIFIED

Jul 20, 2022

PROTOCOL INTEGER ID

66880

MATERIALS TEXT

Phosphate-buffered Saline (ph 7.4)

137 mM NaCl

2.7 mM KCl

10 mM Na₂HPO₄

1.8 mM KH₂PO₄

Lysis Buffer

10 mM Tris pH 7.5

0.5 mM MgCl₂

2 mM DTT

0.1 mM PMSF

Plus protease inhibitors (5 µg/mL aprotinin, 5 µg/mL leupeptin, 1 µg/mL pepstatin A)

Resuspension Buffer

10 mM Tris pH 7.5

0.5 M sucrose

2 mM DTT

0.5 mM PMSF

Storage Buffer

0.25 M sucrose

Plus protease inhibitors (5 µg/mL aprotinin, 5 µg/mL leupeptin, 1 µg/mL pepstatin A, 2 mM PMSF)

1 Thaw Sf9 cell pellets at room temperature (typical size around 5g for 0.4 L of culture)

2 All subsequent steps should be carried out at **4 °C**

3 Wash pellet twice with 15 mL of Phosphate-buffered Saline, centrifuge at **1000 x g, 4°C, 00:07:00** in between washes 7m

3.1 Gently resuspend by inverting tube and pipetting

4 Collect cells after final wash by centrifugation at **1500 x g, 4°C, 00:07:00** 7m

- 5 Resuspend cells in 10 mL Lysis Buffer
- 6 Swell cells in Lysis Buffer by incubating on ice for ⌚00:10:00 to ⌚00:15:00 25m
- 7 Lyse with Dounce homogenizer, 40 strokes tight
- 8 Dilute homogenate in equal volume Resuspension Buffer and mix
- 9 Further lyse with Dounce homogenizer, 20 strokes tight
- 10 Spin down at ⌚1000 x g, 4°C, 00:10:00 to remove nuclear fraction and unbroken cells and save supernatant 10m
- 11 Spin down at ⌚10000 x g, 4°C, 00:20:00 (Sorvall SS-34 rotor) to remove mitochondrial-lysosomal fraction and save supernatant 20m
- 12 Transfer supernatant to ultracentrifuge tubes and spin down at ⌚200000 x g, 4°C, 00:35:00 (Beckman Type 45 Ti rotor) to collect microsomes 35m
- 13 Resuspend microsomal pellet in Storage Buffer
- 14 Measure microsome concentrations based on total protein concentration using the Bradford Assay and bovine serum albumin as a standard
- 15 Flash-freeze in aliquots of 25-100 µL at concentrations between 2-5 mg/mL using liquid

nitrogen and store at -80C until use