

Apr 15, 2021

## 4AA Mix

## Erçağ Pinçe<sup>1</sup>

<sup>1</sup>AMOLF

1 Works for me

dx.doi.org/10.17504/protocols.io.mkpc4vn



ABSTRACT

4 amino acid mixture for MMH1

DOI

dx.doi.org/10.17504/protocols.io.mkpc4vn

PROTOCOL CITATION

Erçağ Pinçe 2021. 4AA Mix. **protocols.io** https://dx.doi.org/10.17504/protocols.io.mkpc4vn

LICENSE

This is an open access protocol distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Jan 12, 2018

LAST MODIFIED

Apr 15, 2021

PROTOCOL INTEGER ID

9583

DISCLAIMER:

## DISCLAIMER - FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to <a href="protocols.io">protocols.io</a> is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with <a href="protocols.io">protocols.io</a>, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

## Mix the ingredients

- 1 40 mL ddH<sub>2</sub>O(milliQ)
  - 620 mg L-Histidine [100mM]
  - 596 mg L-Methionine [100mM]
  - 524 mg L-Leucine [100mM]
  - 476 mg L-Threonine [100mM]
  - Stir for overnight at 37.5C.
  - Adjust pH to 7.0.
  - Filter sterilize and store in the cold room(+4°C)

