



Mar 09, 2022

Construction of a Moore Swab

Dilip Abraham¹, Nicola Elviss², Nicholas Feasey³, Nick Grassly⁴,
Jacob John¹, Gagandeep Kang¹, John Scott Meschke⁵,
Venkata Raghava Mohan¹, Satheesh Nair², Jonathan Rigby³,
Rajan Srinivasan¹, Catherine Troman⁴, Christopher Uzzell⁴, Nicolette Zhou⁵

¹Christian Medical College, Vellore, India; ²UK Health Security Agency;

³Malawi Liverpool Wellcome Trust Clinical Research Program, Blantyre, Malawi;

⁴Imperial College London; ⁵University of Washington

3



protocol .

Typhoid Environmental Surveillance

Catherine Troman

Please note, the author list is in alphabetical order and does not reflect contribution.

Moore swabs are gauze swabs used to capture bacteria in water such as sewage and wastewater. The gauze is held together by a string which is then tied to a stationary object in the water and left for 24-72 hours. During this time, as water moves past the swab bacteria are caught and attach to the swab, allowing isolation and identification for monitoring and surveillance of the bacteria of interest. This protocol based on the Moore swab construction described in Sikorski and Levine (2020).

Dilip Abraham, Nicola Elviss, Nicholas Feasey, Nick Grassly, Jacob John, Gagandeep Kang, John Scott Meschke, Venkata Raghava Mohan, Satheesh Nair, Jonathan Rigby, Rajan Srinivasan, Catherine Troman, Christopher Uzzell, Nicolette Zhou 2022. Construction of a Moore Swab. **protocols.io**
<https://protocols.io/view/construction-of-a-moore-swab-bq84mzyw>



protocol

Reviving the "Moore Swab": a Classic Environmental Surveillance Tool Involving Filtration of Flowing Surface Water and Sewage Water To Recover Typhoidal Salmonella Bacteria Michael J. Sikorski, Myron M. Levine Applied and Environmental Microbiology Jun 2020, 86 (13) e00060-20; DOI: 10.1128/AEM.00060-20

_____ protocol ,

Jan 07, 2021

Mar 09, 2022

46076

The material used to tie the swab together may depend on the nature of the water it is being placed in. For example, if it is faster flowing or there are many solids in the water then you may require a stronger fishing line rather than twine to reduce loss of swabs.

Cotton gauze (6 x 48 inches)

A length of twine or fishing line (nylon thread)

:

Please note, the author list is in alphabetical order and does not reflect contribution.

- 1 Fold the gauze eight times length-wise to form an 8-ply square
- 2 Tie the gauze around the middle with twine/fishing line leaving a long tail for attaching the swab once it has been placed.