



Oct 21, 2021

# Luminol Calibration

**Michael Burgis**<sup>1</sup><sup>1</sup>Philipps-University Marburg[dx.doi.org/10.17504/protocols.io.bzdfp23n](https://dx.doi.org/10.17504/protocols.io.bzdfp23n)**Michael Burgis**

Calibration curve using luminol in order to standardize luminescence data

DOI

[dx.doi.org/10.17504/protocols.io.bzdfp23n](https://dx.doi.org/10.17504/protocols.io.bzdfp23n)Michael Burgis 2021. Luminol Calibration. **protocols.io**<https://dx.doi.org/10.17504/protocols.io.bzdfp23n>

\_\_\_\_\_ protocol ,

Oct 21, 2021

Oct 21, 2021

54407

## Preparation of the Buffers

- 1 Preparation of 100ml luminol stock solution (1 mM luminol sodium salt; 50 mM sodium carbonate; 300 mM sodium bicarbonate; 5 mM ammonium carbonate) :  
dissolve the following components in ddH<sub>2</sub>O  
19,94mg of Luminol sodium salt  
529,94mg of sodium carbonate  
2,52g of sodium bicarbonate  
48mg of ammonium carbonate

Preparation of the coppersulfate solution (1,5 mM CuSO<sub>4</sub>):  
dissolve 23,94mg of CuSO<sub>4</sub> in ddH<sub>2</sub>O

## Preparation of dilution series

- 2 Perform a dilution series in a 384 well plate similar to the protocol provided by the iGEM

foundation on Fluorescein calibration:

- 3 Add 100µl of the CuSO<sub>4</sub> solution in the first well
- 4 Fill every second well in that row with 50µl of water (8 rows in total)
- 5 Transfer 50µl of the first well into the next and carefully pipette up and down to ensure proper mixing of the solutions
- 6 Repeat this until you reach the second last well and then dispose the last 50µl

#### Perform the measurement

- 7 Premix the luminol buffer with hydrogen peroxide in a ratio of 4:1
- 8 Transfer 10µl of this solution into empty wells on the plate
- 9 Use a multichannel pipette to add the inducing copper sulfate solution to the luminol solution
- 10 Transfer the well plate into the plate reader and shake the plate for 30 seconds until you start the measurement
- 11 Measure the luminescence with an exposure time of 250ms
- 12 Plot the data in a linear correlation and use this calibration curve in order to normalize your luminescence data

