



Version 3

Nov 30, 2020

Quick Protocol for Monarch® Plasmid Miniprep Kit (NEB #T1010) V.3

New England Biolabs¹¹New England Biolabs

1

Works for me

dx.doi.org/10.17504/protocols.io.bp9fmr3n**New England Biolabs (NEB)**Tech. support phone: +1(800)632-7799 email: info@neb.comDanielle Freedman
NEB

ABSTRACT

This is the quick version of the Monarch® Plasmid DNA Miniprep Kit Protocol (NEB #T1010). For the full protocol, please click [here](#).

EXTERNAL LINK

<https://www.neb.com/protocols/2015/12/08/quick-protocol-for-monarch-plasmid-miniprep-kit-t1010>

DOI

dx.doi.org/10.17504/protocols.io.bp9fmr3n

EXTERNAL LINK

<https://www.neb.com/protocols/2015/12/08/quick-protocol-for-monarch-plasmid-miniprep-kit-t1010>

PROTOCOL CITATION

New England Biolabs 2020. Quick Protocol for Monarch® Plasmid Miniprep Kit (NEB #T1010).

protocols.io<https://dx.doi.org/10.17504/protocols.io.bp9fmr3n>

Version created by Lenny Teytelman



EXTERNAL LINK

<https://www.neb.com/protocols/2015/12/08/quick-protocol-for-monarch-plasmid-miniprep-kit-t1010>

KEYWORDS

Monarch Kit

LICENSE

———— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Nov 30, 2020

LAST MODIFIED

Nov 30, 2020

OWNERSHIP HISTORY

Nov 30, 2020



Danielle Freedman

NEB

PROTOCOL INTEGER ID

45063

GUIDELINES

For detailed protocol and more information, visit www.neb.com/T1010

The full protocol is available [here](#).

The video protocol is available [here](#).

MATERIALS TEXT

MATERIALS

 **Monarch® Plasmid Miniprep Kit** **New England**

Biolabs Catalog #T1010


ABSTRACT

This is the quick version of the Monarch® Plasmid DNA Miniprep Kit Protocol (NEB #T1010). For the full protocol, please click [here](#).

BEFORE STARTING

- Please review the important information under the “Guidelines” tab before beginning.
- All centrifugation steps should be carried out at 16,000 x g (~13,000 RPM).
- Add 4 volumes of ethanol (≥ 95%) to one volume of Plasmid Wash Buffer 2.
- If precipitate has formed in Lysis Buffer (B2), incubate at 30–37°C, inverting periodically to dissolve.
- Store Plasmid Neutralization Buffer (B3) at 4°C after opening, as it contains RNase A.

- 1 Pellet 1–5 ml (not to exceed 15 OD units) bacterial culture by centrifugation for 30 seconds. Discard supernatant.

 **00:00:30**

1.5 ml of culture is sufficient for most applications. Ensure cultures are not overgrown (12-16 hours is ideal).

- 2 Resuspend pellet in 200 µl Plasmid Resuspension Buffer (B1).



Vortex or pipet to ensure cells are completely resuspended. There should be no visible clumps.

- 3 Add 200 µl Plasmid Lysis Buffer (B2), gently invert tube 5–6 times, and incubate at room temperature for 1 minute. Do not vortex.

 **00:01:00 incubation**

Color should change to dark pink, and solution will become transparent and viscous. Do not vortex.

- 4 Add 400 µl of Plasmid Neutralization Buffer (B3), gently invert tube until neutralized, and incubate at room temperature

for 1 minute. Do not vortex.

🕒 **00:02:00 incubation**

Sample is neutralized when color is uniformly yellow and precipitate forms. Do not vortex.

- 5 Centrifuge lysate for 2–5 minutes.

🕒 **00:02:00 centrifugation**

For culture volumes >1 ml, we recommend a 5 minute spin to ensure efficient RNA removal by RNase A. Pellet should be compact; spin longer if needed.

- 6 Carefully transfer supernatant to the spin column and centrifuge for 1 minute. Discard flow-through.

🕒 **00:01:00 centrifugation**

- 7 Re-insert column in the collection tube and add 200 µl of Plasmid Wash Buffer 1. Centrifuge for 1 minute.

🕒 **00:01:00**

Discarding the flow-through is optional.

- 8 Add 400 µl of Plasmid Wash Buffer 2 and centrifuge for 1 minute.

🕒 **00:01:00 centrifugation**

- 9 Transfer column to a clean 1.5 ml microfuge tube.

Use care to ensure that the tip of the column does not come into contact with the flow-through. If there is any doubt, re-spin the column for 1 minute.

- 10 Add ≥ 30 µl DNA Elution Buffer to the center of the matrix. Wait for 1 minute, then spin for 1 minute to elute DNA.

🕒 **00:01:00 wait**

🕒 **00:01:00 spin**

Nuclease-free water (pH 7–8.5) can also be used to elute the DNA. Yield may slightly increase if a larger volume of DNA Elution Buffer is used, but the DNA will be less concentrated. For larger size DNA, (≥ 10 kb), heating the elution buffer to 50°C prior to use can improve yield.