


FEB 04, 2023

🌐 General subculture protocol for HEK-Blue cells

Andreas Sagen¹

¹The National Institute of Occupational Health in Norway



Andreas Sagen

University of Oslo, The National Institute of Occupational H...

ABSTRACT

HEK-Blue is a product from Invivogen, which provide reporter cells for endotoxin-testing among others. Here is a generalized protocol for culturing cells with reporter characteristics.

GUIDELINES

The protocol has been proven to work with different HEK-Blue cells, but not all. Take time to optimize.

MATERIALS

Centrifuge
Laminar flow cabinet
CO2 incubator

SAFETY WARNINGS



Use appropriate PPE when working with selective antibiotics.

BEFORE START INSTRUCTIONS

Prepare complete growth medium

OPEN ACCESS

Protocol Citation: Andreas Sagen 2023.
General subculture protocol for HEK-Blue cells.
protocols.io
<https://protocols.io/view/general-subculture-protocol-for-hek-blue-cells-cna3vagn>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: In development
We are still developing and optimizing this protocol

Created: Jan 25, 2023

Last Modified: Feb 04, 2023



PROTOCOL integer ID:
75835

Keywords: HEK, Endotoxin, TLR, TOLL-Like receptor

Complete growth medium (CCM)

3h 15m

- 1 For the HEK-Blue cells to work correctly, it has to contain a selective agent.

1.1 In  450 mL DMEM, add  50 mL FBS, and selective antibiotics (cell line dependent; Remove equal volume DMEM/FBS):

Zeocin:  0.5 mL

HEK-Blue selection:  2.0 mL

Blasticidin:  1.5 mL

Materials:

 DMEM high glucose GlutaMAX **Gibco - Thermo Fischer Catalog #31966021**

 ULTRA-LOW ENDOTOXIN FETAL BOVINE SERUM (FBS) **BioWest Catalog #S1860**

 Zeocin **InvivoGen Catalog #ant-zn**

 HEK-Blue Selection **InvivoGen Catalog #hb-sel**

 Blasticidin **InvivoGen Catalog #ant-bl**

1.2 Aliquot complete growth medium in  50 mL or  100 mL units and store refrigerated at  4 °C or frozen at  -20 °C



Note

Create descriptive labels containing the following information: Cell line identifier, DMEM, 10% FBS, selective-antibiotic, Date

Subculture process

40m

2 HEK-Blue cells are easy to subculture. They do not require trypsin, only PBS to dissociate from growth surface. Different cell lines of HEK-Blue grow with different doubling times, and the initial seeding density has to be optimized.



2.1 Add  11 mL warm CIM in a T-75 flask and incubate for  00:30:00 in a 5% CO₂ incubator


30m


2.2 Remove supernatant from T-75 and wash with PBS

Materials:

 PBS pH 7.2 **Gibco - Thermo Fischer Catalog #20012019**

2.3 Add  6 mL PBS and incubate for  00:05:00 5m

2.4 Transfer content to centrifugation tube and centrifuge with  200 rcf, 00:05:00 5m

2.5 Discard supernatant and resuspend pellet in  1 mL warm CIM

2.6 Split between 1:2 to 1:6 of resuspended suspension from centrifugation tube to T-75 and incubate in a 5% CO₂ incubator

Note

Incubate for 2-3 days, then change medium. Passage when cells reach 80-90% confluence