

Apr 26, 2024

# 🌐 Crystallization of Enterovirus coxsackievirus A16 2A protease

DOI

[dx.doi.org/10.17504/protocols.io.3byl49kdzgo5/v1](https://dx.doi.org/10.17504/protocols.io.3byl49kdzgo5/v1)



ryan Lithgo<sup>1,2</sup>, Peter Marples<sup>1,2</sup>, Lizbé Koekemoer<sup>3</sup>, Daren Fearon<sup>1,2</sup>

<sup>1</sup>Diamond Light Source; <sup>2</sup>Research Complex at Harwell; <sup>3</sup>Centre of Medicines Discovery, University of Oxford

ryan Lithgo: The principle crystallographer on the Coxsackievirus A16 project.;

ASAP Discovery



Lizbé Koekemoer

University of Oxford

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.3byl49kdzgo5/v1](https://dx.doi.org/10.17504/protocols.io.3byl49kdzgo5/v1)

**Protocol Citation:** ryan Lithgo, Peter Marples, Lizbé Koekemoer, Daren Fearon 2024. Crystallization of Enterovirus coxsackievirus A16 2A protease. [protocols.io https://dx.doi.org/10.17504/protocols.io.3byl49kdzgo5/v1](https://dx.doi.org/10.17504/protocols.io.3byl49kdzgo5/v1)

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** April 26, 2024

**Last Modified:** April 26, 2024

**Protocol Integer ID:** 98841

**Keywords:** crystallisation, XChem, ASAP, AViDD, CMD, Diamond Light Source, i04-1, Coxsackievirus, A16

**Funders Acknowledgement:**

**National Institutes of  
Health/National Institute Of  
Allergy and Infectious  
Diseases (NIH/NIAID)**

**Grant ID: Grant ID:  
U19AI171399**

**Disclaimer**

The content is solely the responsibility of the authors and does not necessarily represent the official views of the National Institutes of Health.

**Acknowledgements:**

Diamond Light Source Ltd, Harwell Science and Innovation Campus, Didcot OX11 0QX, UK  
Research Complex at Harwell, Harwell Science and Innovation Campus, Didcot OX11 0FA, UK  
Oxford Lab Technologies crystal shifter <https://doi.org/10.1107/S2059798320014114>


**Abstract**

Picornaviridae coxsackievirus A16 is the causative agent of paediatric hand-foot-and-mouth disease, and a target for pandemic preparedness due to the risk of higher order complications in a large-scale outbreak. The 2A protease of the virus is responsible for self-cleavage from the poly protein, allowing for correct folding and assembly of capsid proteins in the final stages of viral replication. Inhibition deranges capsid folding and assembly, preventing formation of mature virions in host cells and making the protease a valuable target for antiviral activity. This protocol was used to grow coxsackievirus A16 crystals that were applied high-throughput crystallographic fragment screening on the target.

**Materials**

<https://swissci.com/product/3-lens-crystallisation-plate/> **Codes:**

*Midi:* UVXPO-3LENS 3W96T-PS 3W96T-UVP

[M] 1 Molarity (M) MES  6.7 , Molecular Dimensions, Catalog # MD2-013-PH 6.7


50% w/v PEG 20000, Molecular Dimensions, Catalog # MD2-250-16

Purified SARS CoV-2 Coxsackievirus A16 protein ( [M] 20 mg/mL ) in [M] 10 millimolar (mM) HEPES,  7.5 ,

[M] 0.5 Molarity (M) NaCl, 5% glycerol, [M] 0.5 millimolar (mM) TCEP

Protein construct <https://www.addgene.org/204809/>

## Safety warnings

 Follow all handling warning for the chemicals used in the crystallisation screen composition.



## Equipment needed

- 1 **Formulatrix Rock Imager** (or incubator of choice)  
**SPT mosquito**

### Equipment

**Mosquito HV**

NAME

High Volume 16-Channel Robotic Liquid Handler

TYPE

SPT LabTech

BRAND

3097-01057

SKU

<https://www.sptlabtech.com/products/liquid-handling/mosquito-hv/><sup>LINK</sup>

P100 8 multi-channel pipette

**SwissCI 3 lens plate**

## Crystallization experiment

1d

- 2 **Prepare seed stock:**

### Protocol



NAME

**Diamond XChem Seeding Protocol**

CREATED BY

Peter Marples

**PREVIEW**

1: 1000 dilution  Sample seeds

- 3 **Protein and buffer requirements:**



43.2  $\mu$ L



20 mg/mL



Sample



3.36 mL Crystallization screen



14.4 µL Sample seeds, dilution 1:1000

#### 4 **Crystallisation screen composition:**

13.5 % PEG 20000

[M] 0.1 Molarity (M) MES 6.7

#### **Stock solutions used:**

[M] 1 Molarity (M) MES 6.7

50% w/v PEG 20000

#### **Note**

The crystallisation screen can be stored in a duran bottle or aliquoted into 96 deep well block for easy dispensing into SwissCI 3 lens plates.

For long term storage keep the Crystallisation screen in the fridge at 4°C.

5 Dispense 35 µL Crystallisation screen into SwissCI 3 lens plate reservoir wells using a 100 µl multi-channel pipette.

Dispense 150 nL [M] 20 mg/mL Sample to each lens using the SPT mosquito.

Dispense 150 nL Crystallisation screen to each lens using the SPT mosquito.

Dispense 50 nL Seeds to each lens using the SPT mosquito.

**Drop ratio:** 3:3:1 ratio (150 nl Sample : 150 nl reservoir solution: 50 nl seeds)

**Final drop volume:** 350 nl

6 Incubate at 20 °C for 24:00:00 h in Formulatrix Rock Imager.

1d

**Imaging Schedule:** The first images are taken after 12hrs and the imaging schedule follows a Fibonacci sequence of days for further collections.

7 Crystal typically form after ~24hrs

### Expected result

Crystals typically reach their maximum size after ~24 h.


**Morphology:** typically rectangles.

**Size:** ~75  $\mu\text{m}$  in length and ~10  $\mu\text{m}$  in width, depth of the crystals is ~10  $\mu\text{m}$ , giving a rectangular appearance

**Average resolution:** 1.6  $\text{\AA}$

**Space group:** C2

**Unit cell:** 86  $\text{\AA}$ , 57  $\text{\AA}$ , 32  $\text{\AA}$   
90°, 95°, 90°



An example of a drop containing Cocksackievirus A16 crystals.



## Data collection at Synchrotron

8 Diamond Light Source  
**Unattended Data Collection (UDC)**  
**Data Collection Temperature:** 100K  
**Detector:** DECTRIS EIGER2 X 9M  
**Beamline:** I04-1  
**Wavelength:** 0.9212 Å  
**Resolution (Å):** 1.21  
**Beam Size (µm):** 60 X 50  
**Number of images:** 3600  
**Oscillation:** 0.10°  
**Exposure (s):** 0.0020  
**Transmission (%):** 100  
**Flux (ph/s):** 9.50e+11

## Protocol references

N/A