



Oct 16, 2020

# 1M KPO<sub>4</sub> pH6.0

Ida Barlow<sup>1</sup><sup>1</sup>Imperial College London

1

Works for me

dx.doi.org/10.17504/protocols.io.bnh2mb8e

## Behavioural Genomics



Ida Barlow

### ABSTRACT

Protocol for making phosphate salt solution for adding to NGM

### DOI

dx.doi.org/10.17504/protocols.io.bnh2mb8e

### PROTOCOL CITATION

Ida Barlow 2020. 1M KPO<sub>4</sub> pH6.0. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.bnh2mb8e>



### LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

### CREATED

Oct 16, 2020

### LAST MODIFIED

Oct 16, 2020

### PROTOCOL INTEGER ID

43290

### MATERIALS

NAME	CATALOG #	VENDOR
Potassium phosphate (dibasic)		P212121
Potassium phosphate (monobasic)		P212121
Potassium hydroxide		P212121

### ABSTRACT

Protocol for making phosphate salt solution for adding to NGM

### Make 1M solutions of monobasic and dibasic salts

#### 1 Prepare 1M solution dibasic salt:

**174.18 g** K<sub>2</sub>HPO<sub>4</sub>
**1 L** water

#### 2 Prepare 1M solutions of monobasic salt:

🪴 **136.09 g**  $\text{KH}_2\text{PO}_4$

🪴 **1 L** water

- 3 Autoclave both solutions and allow to cool

Mix salts to achieve correct pH

- 4 Mix prepared salt solutions to achieve 1M  $\text{KPO}_4$  solution close to pH6:

🪴 **100 mL** 1M  $\text{K}_2\text{HPO}_4$  (dibasic)

🪴 **400 mL** 1M  $\text{KH}_2\text{PO}_4$  (monobasic)

Adjust pH to pH6.0 using pH meter and addition of drops of KOH